
ENGLISH	
CONTOUR DENTAL ALLOY	
CAUTION <div>FEDERAL LAW RESTRICTS THIS DEVICE TO SALE BY OR ON THE ORDER OF A DENTIST.</div>	
INDICATIONS FOR USE <div>Contour dental alloy is a dispersed-phase (admix) alloy for use in the preparation of dental amalgam.</div>	
COMPOSITION <div>Ag = 41% Sn = 31% Cu = 28% (Contour dental alloy is Zinc Free) 70% Spherical 30% Lathe Cut</div>	
INSTRUCTIONS FOR USING CONTOUR WITH THE DENTOMAT AND DUOMAT BULK POWDER AMALGAMATORS <div>Although the composition of Contour dental alloy has been formulated specifically for use in the Dentomat or Duomat Bulk Powder Amalgamators, slight deviations from the standard operating procedures for those amalgamators are required.</div>	
Mixing Ratio <div>Contour requires less mercury than many other alloys and therefore, the amalgamator’s regulating screw should be adjusted for a leaner mixture. It is recommended that the arrow be positioned as shown below.</div> 	
NOTE: The unique character of Contour requires only a very slight adjustment of the regulating screw to change the mix consistency. Avoid excessively wet mixes, which can cause subsequence inconsistent mixes because pre-mixed amalgam has tendency to hang up in dispensing tube.	

WARNING <ul style="list-style-type: none">This product, when used to create a dental amalgam, will be mixed with mercury. Mercury vapors may be harmful if inhaled. See Mercury Health Hazard Data Section below.
PRECAUTIONS <p>Proper care should be taken when handling mercury and dental amalgams. Protective measures such as the following should be taken:</p> <ul style="list-style-type: none">The placement of a dental amalgam in a patient will increase the level of mercury in the body of the patient. The use of a rubber dam may decrease the amount of mercury absorbed by a patient during the removal or placement of an amalgam. Adequate ventilation should be provided when mixing the powder with mercury and using the mixed product (vapor exposure should be kept under 0.05 mg/m³ or any lower limit set by any governmental or regulatory agency governing your use of the product). Store in a cool, well-ventilated place. Do not place dental amalgams in direct contact with other types of metals. Dental amalgams are for single use only. Wear gloves when handling Contour dental alloy, mercury, and dental amalgams. Use an enclosed amalgamator when mixing Contour dental alloy with mercury. Properly dispose of spent capsules and any excess unused amalgam. Use HGX or similar type mercury-absorbing compounds in the event of spillage of mercury or dental amalgams. If a patient experiences a localized hypersensitivity reaction to dental amalgam, the amalgam should be removed.

PRECAUTIONS <p>Proper care should be taken when handling mercury and dental amalgams. Protective measures such as the following should be taken:</p> <ul style="list-style-type: none">The placement of a dental amalgam in a patient will increase the level of mercury in the body of the patient. The use of a rubber dam may decrease the amount of mercury absorbed by a patient during the removal or placement of an amalgam. Adequate ventilation should be provided when mixing the powder with mercury and using the mixed product (vapor exposure should be kept under 0.05 mg/m³ or any lower limit set by any governmental or regulatory agency governing your use of the product). Store in a cool, well-ventilated place. Do not place dental amalgams in direct contact with other types of metals. Dental amalgams are for single use only. Wear gloves when handling Contour dental alloy, mercury, and dental amalgams. Use an enclosed amalgamator when mixing Contour dental alloy with mercury. Properly dispose of spent capsules and any excess unused amalgam. Use HGX or similar type mercury-absorbing compounds in the event of spillage of mercury or dental amalgams. If a patient experiences a localized hypersensitivity reaction to dental amalgam, the amalgam should be removed.	9006
European Union Representative: Kerr Corporation 28200 Wick Road Romulus, MI 48174-2600 Kerr Italia S. r.l. Via Passanti, 332 I-84018 Scafalli (SA) Italy TEL.: + 39-081-851-8311 kerithawe.com 27969A 9109	9006

STATEMENTS BY GOVERNMENTAL AGENCIES REGARDING RESTRICTIONS AND RECOMMENDATIONS ON THE USE OF DENTAL AMALGAMS <div>The U.S. Food and Drug Administration states: "Dental amalgam has been demonstrated to be an effective restorative material that has benefits in terms of strength, marginal integrity, suitability for large occlusal surfaces, and durability. Dental amalgam also releases low levels of mercury vapor, a chemical that at high exposure levels is well-documented to cause neurological and renal adverse health effects. Mercury vapor concentrations are highest immediately after placement and removal of dental amalgam but decline thereafter.</div> <div>Clinical studies have not established a causal link between dental amalgam and adverse health effects in adults and children age six and older. In addition, two clinical trials in children aged six and older did not find neurological or renal injury associated with amalgam use.</div>	
---	--

The developing neurological systems in fetuses and young children may be more sensitive to the neurotoxic effects of mercury vapor. Very limited to no clinical information is available regarding long-term health outcomes in pregnant women and their developing fetuses, and children under the age of six, including infants who are breastfed.

The Agency for Toxic Substances and Disease Registry’s (ATSDR) and the Environmental Protection Agency (EPA) have established levels of exposure for mercury vapor that are intended to be highly protective against adverse health effects, including for sensitive subpopulations such as pregnant women and their developing fetuses, breastfed infants, and children under age six. Exceeding these levels does not necessarily mean that any adverse effects will occur.

FDA has found that scientific studies using the most reliable methods have shown that dental amalgam exposes adults to amounts of elemental mercury vapor below or approximately equivalent to the protective levels of exposure identified by ATSDR and EPA. Based on these findings and the clinical data, FDA has concluded that exposures to mercury vapor from dental amalgam do not put individuals age six and older at risk for mercury-associated adverse health effects.

Taking into account factors such as the number and size of teeth and respiratory volumes and rates, FDA estimates that the estimated daily dose of mercury in children under age six with dental amalgams is lower than the estimated daily adult dose. The exposures to children would therefore be lower than the protective levels of exposure identified by ATSDR and EPA.

In addition, the estimated concentration of mercury in breast milk attributable to dental amalgam is an order of magnitude below the EPA protective reference dose for oral exposure to inorganic mercury. FDA has concluded that the existing data support a finding that infants are not at risk for adverse health effects from the breast milk of women exposed to mercury vapors from dental amalgam.”

- The Health Authorities of Canada indicate:
- Non-mercury filling materials should be considered for restoring the primary teeth of children where the mechanical properties of the material are suitable.
 - Whenever possible, amalgam fillings should not be placed in or removed from the teeth of pregnant women.

- Amalgam should not be placed in patients with impaired kidney function.
- In placing and removing amalgam fillings, dentists should use techniques and equipment to minimize the exposure of the patient and the dentist to mercury vapor and to prevent amalgam waste from being flushed into municipal sewage systems.
- Dentist should advise individuals who may have allergic hypersensitivity to mercury to avoid the use of amalgam. In patients who have developed hypersensitivity to amalgam, existing amalgam restorations should be replaced with another material where this is recommended by a physician.

The health authorities of various countries including Germany, France, the United Kingdom, Norway and Austria have also made statements regarding the use and safety of dental amalgams. These statements include recommendations against the placement or removal of an amalgam in certain individuals such as pregnant and nursing women and persons with impaired kidney function. You should check with the authorities in your country that govern the practice of dentistry and dental materials to determine what recommendations or restrictions apply to the use of dental amalgams.

MERCURY HEALTH HAZARD DATA	
Effects and Hazards of Eye Contact: Irritant	
Acute Exposure: Contact may cause irritation. Mercury is corrosive and may cause corneal injury or burns.	
Chronic Exposure: Mercury may be deposited in the lens of the eye, causing visual disturbances.	
Effects and Hazards of Skin Contact: Irritant/Sensitizer/Neurotoxin/Nephrotoxin	
Acute Exposure: May cause redness and irritation	
Chronic Exposure: Possible sensitization, dermatitis, and swelling. Mercury may be absorbed through the skin causing urinary problems.	
Effects and Hazards of Inhalation: Irritant/Sensitizer/Neurotoxin	
Acute Exposure: Inhalation of mercury vapor can cause cough, fever, nausea, and vomiting.	
Chronic Exposure: Inhalation of high concentrations of mercury vapor over a long period may cause teething, blue lines on gums, pain, and numbness in extremities.	
Effects and Hazards of Ingestion: Neurotoxic/Nephrotoxic	
Acute Exposure: May cause nausea, vomiting, kidney damage and nerve effects.	
Chronic Exposure: Symptoms include central nervous system (CNS) disorders.	

LIMITED WARRANTY AND LIMITATION ON LIABILITY	
Kerr's technical advice, whether verbal or in writing, is designed to assist dentists in using Kerr's product. The dentist assumes all risk and liability for damages arising out of the improper use of Kerr's product. In the event of a defect in material or workmanship, Kerr's liability is limited, at the Kerr's option, to replacement of the defective product or part thereof, or reimbursement of the actual cost of the defective product. In order to take advantage of this limited warranty, the defective product must be returned to Kerr. KERR IS NOT RESPONSIBLE FOR DIRECT, SPECIAL, INCIDENTAL OR CONSEQUENTIAL DAMAGES RESULTING FROM ANY BREACH OF WARRANTY OR CONDITION, OR UNDER ANY OTHER LEGAL THEORY, INCLUDING BUT NOT LIMITED TO LOST PROFITS, DOWNTIME, GOODWILL, DAMAGE TO OR REPLACEMENT OF EQUIPMENT AND PROPERTY. THIS WARRANTY AND REMEDIES SET FORTH ABOVE ARE EXCLUSIVE AND IN LIEU OF ALL OTHER WARRANTIES, REMEDIES AND CONDITIONS, WHETHER EXPRESS OR IMPLIED. KERR SPECIFICALLY DISCLAIMS ANY AND ALL IMPLIED WARRANTIES, INCLUDING, WITHOUT LIMITATION, WARRANTIES OF MERCHANTABILITY AND FITNESS FOR A PARTICULAR PURPOSE.	
CONTRAINDICATIONS	
<ul style="list-style-type: none">Compressive strength <ul style="list-style-type: none">Regular Set: 533MPa @ 24 hrs Dimensional change during hardening <ul style="list-style-type: none">Regular Set: -0.09%	


WORKING TIMES OF DENTAL AMALGAMS MADE USING CONTOUR DENTAL ALLOY *

Pick-up and Condensing Time:	Regular Set - 2.25 mins.
Carving Time:	Regular Set - 9.75 mins.
*All times are minimum	

PHYSICAL PROPERTIES OF DENTAL AMALGAMS MADE USING CONTOUR DENTAL ALLOY	
<ul style="list-style-type: none">Compressive strength <ul style="list-style-type: none">Regular Set: 533MPa @ 24 hrs Dimensional change during hardening <ul style="list-style-type: none">Regular Set: -0.09%	
WARNING	
<ul style="list-style-type: none">This product, when used to create a dental amalgam, will be mixed with mercury. Mercury vapors may be harmful if inhaled. See Mercury Health Hazard Data Section below.	

PRECAUTIONS <p>Proper care should be taken when handling mercury and dental amalgams. Protective measures such as the following should be taken:</p> <ul style="list-style-type: none">The placement of a dental amalgam in a patient will increase the level of mercury in the body of the patient. The use of a rubber dam may decrease the amount of mercury absorbed by a patient during the removal or placement of an amalgam. Adequate ventilation should be provided when mixing the powder with mercury and using the mixed product (vapor exposure should be kept under 0.05 mg/m³ or any lower limit set by any governmental or regulatory agency governing your use of the product). Store in a cool, well-ventilated place. Do not place dental amalgams in direct contact with other types of metals. Dental amalgams are for single use only. Wear gloves when handling Contour dental alloy, mercury, and dental amalgams. Use an enclosed amalgamator when mixing Contour dental alloy with mercury. Properly dispose of spent capsules and any excess unused amalgam. Use HGX or similar type mercury-absorbing compounds in the event of spillage of mercury or dental amalgams. If a patient experiences a localized hypersensitivity reaction to dental amalgam, the amalgam should be removed.	
---	--

PRECAUTIONS <p>Proper care should be taken when handling mercury and dental amalgams. Protective measures such as the following should be taken:</p> <ul style="list-style-type: none">The placement of a dental amalgam in a patient will increase the level of mercury in the body of the patient. The use of a rubber dam may decrease the amount of mercury absorbed by a patient during the removal or placement of an amalgam. Adequate ventilation should be provided when mixing the powder with mercury and using the mixed product (vapor exposure should be kept under 0.05 mg/m³ or any lower limit set by any governmental or regulatory agency governing your use of the product). Store in a cool, well-ventilated place. Do not place dental amalgams in direct contact with other types of metals. Dental amalgams are for single use only. Wear gloves when handling Contour dental alloy, mercury, and dental amalgams. Use an enclosed amalgamator when mixing Contour dental alloy with mercury. Properly dispose of spent capsules and any excess unused amalgam. Use HGX or similar type mercury-absorbing compounds in the event of spillage of mercury or dental amalgams. If a patient experiences a localized hypersensitivity reaction to dental amalgam, the amalgam should be removed.	9006
European Union Representative: Kerr Corporation 28200 Wick Road Romulus, MI 48174-2600 Kerr Italia S. r.l. Via Passanti, 332 I-84018 Scafalli (SA) Italy TEL.: + 39-081-851-8311 kerithawe.com 27969A 9109	9006

STATEMENTS BY GOVERNMENTAL AGENCIES REGARDING RESTRICTIONS AND RECOMMENDATIONS ON THE USE OF DENTAL AMALGAMS <div>The U.S. Food and Drug Administration states: "Dental amalgam has been demonstrated to be an effective restorative material that has benefits in terms of strength, marginal integrity, suitability for large occlusal surfaces, and durability. Dental amalgam also releases low levels of mercury vapor, a chemical that at high exposure levels is well-documented to cause neurological and renal adverse health effects. Mercury vapor concentrations are highest immediately after placement and removal of dental amalgam but decline thereafter.</div> <div>Clinical studies have not established a causal link between dental amalgam and adverse health effects in adults and children age six and older. In addition, two clinical trials in children aged six and older did not find neurological or renal injury associated with amalgam use.</div>	
---	--

The developing neurological systems in fetuses and young children may be more sensitive to the neurotoxic effects of mercury vapor. Very limited to no clinical information is available regarding long-term health outcomes in pregnant women and their developing fetuses, and children under the age of six, including infants who are breastfed.

The Agency for Toxic Substances and Disease Registry’s (ATSDR) and the Environmental Protection Agency (EPA) have established levels of exposure for mercury vapor that are intended to be highly protective against adverse health effects, including for sensitive subpopulations such as pregnant women and their developing fetuses, breastfed infants, and children under age six. Exceeding these levels does not necessarily mean that any adverse effects will occur.

FDA has found that scientific studies using the most reliable methods have shown that dental amalgam exposes adults to amounts of elemental mercury vapor below or approximately equivalent to the protective levels of exposure identified by ATSDR and EPA. Based on these findings and the clinical data, FDA has concluded that exposures to mercury vapor from dental amalgam do not put individuals age six and older at risk for mercury-associated adverse health effects.

Taking into account factors such as the number and size of teeth and respiratory volumes and rates, FDA estimates that the estimated daily dose of mercury in children under age six with dental amalgams is lower than the estimated daily adult dose. The exposures to children would therefore be lower than the protective levels of exposure identified by ATSDR and EPA.

In addition, the estimated concentration of mercury in breast milk attributable to dental amalgam is an order of magnitude below the EPA protective reference dose for oral exposure to inorganic mercury. FDA has concluded that the existing data support a finding that infants are not at risk for adverse health effects from the breast milk of women exposed to mercury vapors from dental amalgam.”

- The Health Authorities of Canada indicate:
- Non-mercury filling materials should be considered for restoring the primary teeth of children where the mechanical properties of the material are suitable.
 - Whenever possible, amalgam fillings should not be placed in or removed from the teeth of pregnant women.

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