PARKELL PRODUCT INFORMATION STOCK NO. S020

OTHER SEE

(A) MSDS INFOR SIDE

(A) Ingredients Instead on MSDS)

CAUTION: US Federal Law restricts this material to sale by or on the order of a licensed professional.

Instructions for Liquid AccuFilm® IV

DESCRIPTION:

AccuFilm® IV is a non-toxic, easy-to-use blue marking liquid that assures clear, crisp markings of prematurities, high-spots and pressure areas on all dental materials and tooth enamel.

DIRECTIONS FOR USE:

- 1. Shake well before using.
- 2. Apply with brush for accurate placement and economical use.
- 3. When used intraorally, dispense a small amount of liquid AccuFilm IV into a dappen dish and apply with suitable disposable brush or paper point to avoid possibility of cross-contamination of stock bottle.
- 4. Allow to dry completely (about 5 seconds) before registration of markings.
- 5. AccuFilm IV is easy to remove with a cotton pellet and clean up solvent provided in package. Or alcohol may be used.

SUGGESTIONS FOR USE:

- 1. Occlusal adjustments in all excursions.
- 2. Balancing side prematurities.
- 3. Gothic arch tracings.
- 4. Contacts of restorations in any dental material.
- 5. Internal binding areas of crowns, inlays or bridges.
- 6. Occlusal contour adjustments of any dental restorative.
- 7. Binding areas on cast partial denture frames.
- 8. Pressure spot indicator on periphery of denture bases.
- 9. Used on dies, internal casting adjustments are more accurate.
- 10. Anywhere clean, crisp, easy-to-read markings are useful.

SUPPLIED AS:

Stock No. S020. Package includes 10ml. AccuFilm® IV with applicator brush, 20ml. clean-up solvent and simple instructions.

Contraindications:

This product is contraindicated for use in or by individuals who have shown hypersensitivity to any of its components.



ORDER TOLL-FREE 1.800.243.7446 USA & CANADA

PARKELL, Inc., 300 Executive Drive, Edgewood, NY 11717 USA Toll-Free 1-800-243-7446 Fax: (631) 249-1242 • E-Mail: sales@parkell.com • www.parkell.com

European Authorized Representative:

(Not a dealer/distributor)

MDCI Ltd. Arundel House, 1 Liverpool Gardens Worthing, West Sussex BN11 1SL, UK A00176revD-1007



MATERIAL	PARKELL INC.	24-HOUR
SAFETY	300 EXECUTIVE DRIVE	EMERGENCY
DATA SEEET	EDGEWOOD, NY 11717	TELEPHONE
השווה שווהם	631-249-1134	1-800-535-5053

SECTION I - PRODUCT IDENTIFICATION

PRODUCT NAME:	ACCU-FILM IV		STOCK NO:	8020
DOT HAZARD LABEL:	Flammable Liquid		UN NUMBER:	UN1170
PROPER SHIPPING NAME:	Ethyl Alcohol Solution		DATE PREPARED: 05/01/09	05/01/09
NFPA CODES:	HEALTH 0	FLAMMABILITY -	3 REACTIVITY -	ITY . 0

SECTION II - HAZARDOUS INGREDIENTS/IDENTITY INFORMATION

HAZARDOUS COMPONENTS	×1	CAS NUMBER	DEL	TLV	প্ল
Ethyl alcohol solution		64-17-5	1000 ppm	1000 ppm	
SECTIC	N III - PHYSICAL	SECTION III - PHYSICAL & CHEMICAL CHARACTERISTICS	CTERISTIC	δί	
BOILING POINT:	172° F	SPECIFIC GRAVITY (H ₂ O = 1):	;;; = 1;	1.38	
VAPOR PRESSURE:	40 mm HG @ 19° F	40 mm HG @ 19° F PERCENT VOLATILES:		%05	_
VAPOR DENSITY (Air = 1):	1.59	EVAPORATION RATE (Butyl Acetate = 1):	ityl Acetate = 1)	: 1.5	
APPEARANCE AND ODOR: Blue liquid, alcoholic odor.	Blue liquid, alcoholic	odor.			

SECTION IV - FIRE AND EXPLOSION HAZARD DATA

FLASH POINT (TCC): 56° F FLAMIN	56° F FLAMMABLE LIMIT (air, % by vol.)	UPPER: 39	3% LOWER:	15%
EXTINGUISHER MEDIA: Carbon dioxide or dry chemical.	r dry chemical.			
SPECIAL FIRE FIGHTING PROCEDURES: During emergency conditions, overexposure to thermal decomposition products may cause health hazard. Self contained breathing apparatus should be worn.	During emergency conditions, ove ined breathing apparatus should be	rexposure to therms worn.	al decomposition	
UNUSUAL FIRE & EXPLOSION HAZARDS: When exposed to flame, emits toxic fumes and gases.	S: When exposed to flame, emits to	xic fumes and gase	Š	
SECTION V - RE	SECTION V - REACTIVITY DATA (PHYSICAL HAZARDS)	AL HAZARDS	7	
STABLLITY: X STABLE CONDITIONS TO AVOID: Heat and flame.	UNSTABLE			
INCOMPATIBILITY (Materials to avoid):	Oxidizers, peroxides, persulfates, light, heat, nitric acid, ammonia and amines, halogens, and halogen compounds.	ight, heat, nitric aci npounds.	d, ammonia and	
HAZARDOUS DECOMPOSITION None PRODUCTS:				
HAZARDOUS POLYMERIZATION: CONDITIONS TO AVOID: NA	MAY OCCUR	×	X WILL NOT OCCUR	8

	MATERIAL SAFETY DATA SHEET (con	ATA SHEET (continued)	
PRODUCT NAME:	ACCU-FILM IV	STOCK NO:	8020

SECTION VI - HEALTH HAZARD DATA

PRIMARY ROUTE(S) OF ENTRY: X EYES X SKIN X INHALATION X INGESTIGON REALTH HAZANGSK (ACRUE & Chronic). Middly sover by inhalation and skin contact. Human system certifies by ingestion and subcutaneous tours step disorders, hallucinations, disorder perceptions, convolsions, moor activity changes, ataxia, coma, antipsycotic, headsethe, pulmonary changes, alteration in gastric secretion, nauses or vomiting, other agreements and any antipsycotic cheadsethe, pulmonary changes, alteration in gastric secretion, nauses or vomiting, other agreements and any acrue of the changes and obly emperature decreases. Any exp and skin intains. Ingestion of large doss sen cause a lockologism. It is a central inervous system	
NHALATION X tact. Human systemic ef ons, convulsions, motor e tric secretion, nausea or An eye and skin irritant. An eye and skin irritant. It is a central nervous	
NHALATION tact. Human sy ons, convulsion oritic secretion, n An eye and skir An it is a centra	
NHAI tact. ons, c stric se An ey	
L Con Con Eptin	
d perc d perc tion in fecreas	
SKIN lation and s, distorte s, altera perature of	
ic by inha ucination ary chang body tem	
EYES fildly tox ders, hall e, pulmon nges and peated in	
nic): N p disor sadacho cle cha	
OF ENTRY: Acute & Chro bus route: slee ntipsycotic, ho menstrual cy,	
ROUTE(S) - AZARDS (, I subcutance cia, coma, a nal changes, an cause ale	
PRIMARY I HEALTH H ingestion and changes, atas gastrointestif large doses of	

SIGNS & SYMPTOMS OF EXPOSURE: Exposure to concentrations over 1,000 ppm may cause headache, irritation of the eyes, nose and throat, and, if continued for an hour, drowsiness and lassitude, loss of appetite and inability to concentrate. Voc JARC MONOGRAPHS? No. NTP CARCINOGENECITY

Tes IMAC MUNOUMAFINS: NO USER		INHALATION: Remove victim to fresh air. If cough or other respiratory symptoms develop, consult medical personna	If not breathing, give artificial respiration, preferably mouth-to-mouth. If breathing is difficult, give
NO NIC	ID PROCEDURES:	im to fresh air. If cough or other respirate	ing, give artificial respiration, preferably n
CARCINOGENECII I.	EMERGENCY AND FIRST AID PROCEDURES:	INHALATION: Remove vici	If not breath

Immediately flush with copious amounts of water, including under eyelids, for at least 15 minutes. Get oxygen. Consult medical personnel. EYES:

Wash skin with copious amounts of soap and water. If irritation exists, get medical attention. Wash contaminated clothing and decontaminate footwear before reuse SKIN

Give several glasses of water. Induce vomiting. Get medical attention immediately. INGESTION:

SECTION VII - SPECIAL PRECAUTIONS AND SPILL OR LEAK PROCEDURES

PRECAUTIONS TO BE TAKEN IN HANDLING & STORAGE: Keep away from heat and sources of ignition OTHER PRECAUTIONS: NE

STEPS TO BE TAKEN IN CASE MATERIAL IS RELEASED OR SPILLED: Absorb spill thoroughly as residue is slippery. Non-corrosive to waste containers. WASTE DISPOSAL METHODS (Consult federal, state, and local regulations): Dispose of material in accordance with all Federal, State and Local regulations.

SECTION VII. SPECIAL PROTECTION INFORMATION/CONTROL MEASURES

RESPIRATORY PROTECTION: NIOSH-approved respiratory protection for organic gases if needed. VENTILATION: Use local exhaust to keep exposures to a minimum.

EYE PROTECTION: Safety glasses or goggle if needed. PROTECTIVE GLOVES: Rubber or PVC if needed.

OTHER PROTECTIVE CLOTHING OR EQUIPMENT: Eyewash station.

WORK/HYGIENIC PRACTICES: Wash hands before eating, drinking or smoking.

The information accumulated herein is believed to be accurate but is not warranted to be whether originating with the company or not. Recipients are advised to confirm in advance of need that the information is current, applicable, and suitable to their circumstances.