<b>3M ESPE</b> <sup>®</sup>
<b>Durelon</b> <sup>™</sup>
<b>Carboxylate Cement</b>
<b>Carboxylatzement</b>
<b>Ciment au carboxylate</b>
<b>Cemento carbossilico</b>
<b>Cemento de carboxilato</b>
<b>Cimento carboxilado</b>
<b>Carboxylatzement</b>
<b>Πολυκαρβοξυλική κονία</b>
<b>Καρχοξυλατήρημα</b>
<b>Karboxylaattipöytäaineen semenntti</b>
<b>Carboxylatzement</b>
<b>Karboxylatzement</b>

Instructions for Use
Gebrauchsinformation
Mode d'emploi
Informazioni per l'uso
Información de uso
Instruções de uso

Gebrauchsinformatie
Πληροφορίες χρήσης
Bruksanvisning
Käyttöinfo
Bruksanvisning
Bruksinformasjon

<b>3M ESPE AG</b>	<b>CE</b>	<b>0123</b>
Dental Products		
D-82229 Seefeld - Germany		

**3M ESPE** Dental Products
St. Paul, MN U.S.A. 55144-1000
3M ESPE Corporate Care/MSDS Information in U.S.A. and Canada 1-900-634-2249.

3M, ESPE and Durelon are trademarks of 3M or 3M ESPE AG. © 3M 2008. All rights reserved.

<b>ENGLISH</b>
----------------

**Product Description**
Durelon<sup>™</sup>, manufactured by 3M ESPE, is a carboxylate cement for fixation and lining. Two variants of the powder, "normal curing" and "fast curing", and two consistencies of the liquid, "universal" and "low viscosity", are available.

Two delivery variants of the "universal" liquid, dropper bottle and dosing syringe, are available. The "low viscosity liquid" is available with a dropper bottle only. Either of the powder variants and the "universal" liquid may be used for lining. Any of the combinations can be employed for fixation.

⚠️ These Instructions for Use should not be discarded for the duration of product use.

**Fields of Application**

- Cementation of inlays, onlays, crowns, and bridges made from metal or metal-ceramics or veneered with composite
- Cementation of inlays, onlays, crowns, and bridges made from composite or ceramics provided these are suitable for conventional cementing
- Cementation of pins and screws provided these are suitable for conventional cementing
- Cementation of orthodontic bands
- Linings

Durelon is not suitable for cementing titanium frameworks.

**Precautionary Measures**
3M ESPE MSDSs can be obtained from www.mmm.com or contact your local subsidiary.

**Preparation**

► To provide for optimal adhesion, carefully clean enamel, dentin, and the restoration to be cemented.

► Blow the cavity dry in only 2–3 intervals with water and oil-free air or dab dry with cotton pellets. Do not overdry!

► **The cavity should be so dry that the surface has a semi-gloss appearance.** As is the case with every type of cement, overdrying can result in post-operative sensitivity.

► Dry the surfaces of the restoration with air.

► Avoid further contamination by saliva.

**Pulp protection**

To prevent pulp irritation, cover areas in vicinity to the pulp by applying small amounts of hardsetting calcium hydroxide material prior to taking an impression for the final restoration.

**Dosing:**

**Powder:**

► Assemble the dosing rod prior to use.

► Press the dosing rod into the powder, and then remove excess powder with some instrument.

► Dose the powder onto a mixing block or glass plate by lightly tapping the dosing rod.
– Avoid underdosing the powder (especially in free dosing), since this may impair the quality of the cement.

**Liquid Dropper Bottle:**

► Holding the bottle vertically, dispense single droplets next to the powder by lightly squeezing the bottle.

► Recap the bottle immediately after use.

**Liquid Dosing Syringe:**

► Insert the enclosed plunger and dose the desired number of scale units.

– Keep the plunger for subsequent use with other syringes.

► Reseal the dosing syringe immediately after use.

**Mixing ratio**

1 dose powder = 1 dosing rod filling.

1 dose liquid = 1 drop from the bottle or 1 small scale unit from the dosing syringe.

Indication	Powder		Liquid	
	normal setting	fast setting	universal	low viscosity
Number of doses for fixations	1	–	2	–
	–	1	2	–
	1	–	–	1
Number of doses for linings	1	–	1	–
	–	1	1	–

The following powder/liquid mixing ratios apply:

for fixations: 0.15 g powder : 0.10 g liquid universal or 0.20 g powder : 0.10 g liquid low viscosity,

for linings: 0.25 g powder : 0.10 g liquid universal.

These quantities are not identical to the quantities required for one application!
**Mixing**
► Mix Durelon at an ambient temperature of 20–25°C/68–77°F with a metal spatula.
► Place all powder into the liquid in one move and mix until a homogeneous consistency is attained.

**Application**

► Protect the working area from water and saliva during the application and setting phase.

► Coat neighboring teeth or crowns with Vaseline to make removal of excess material easier.

► Apply a thin layer of cement to the internal surface of the crown and the prepared tooth.
– Do not overfill the crown as hydrostatic pressure in tightly fitting crowns may lead to discomform of the pulp.

► Then seal the restoration.

► Wash metal instruments with cold water before Durelon has set.

**Times**

The following times are valid at an ambient temperature of 23°C/73°F and 50% relative humidity.

	Mixing	Processing	Setting from
	min <span> </span> :sec	incl. mixing min <span> </span> :sec	start of mixing min <span> </span> :sec
<b>Fixation</b>			
Powder, normal setting Liquid, universal	00:30	02:30	10:00
Powder, normal setting Liquid, low viscosity	00:30	02:00	07:00
Powder, fast setting Liquid, universal	00:30	01:00	05:00

	Mixing	Processing	Setting from
	min <span> </span> :sec	incl. mixing min <span> </span> :sec	start of mixing min <span> </span> :sec
<b>Lining</b>			
Powder, normal setting Liquid, universal	00:30	01:30	06:00
Powder, fast setting Liquid, universal	00:30	00:45	03:30

Processing times are shortened at higher temperatures and prolonged at lower temperatures than those stated (e.g., mixing on a chilled glass plate). Using a larger quantity of powder also shortens the processing time. Exceeding these processing times (apparent by the matte surface of the paste) leads to a loss of adhesion to enamel and dentin.

**Removing Excess Material**

► To simply clean, remove any excess material shortly before definitive setting of the material using an instrument. Be careful to avoid pulling any material from the seam.
► After the setting time has elapsed remove any remaining excess material, e.g., with a probe.

**After Removal of Excess and Setting**

► Examine carefully the sulcus of the treated teeth and the surrounding area; remove any remaining cement residues.

**Notes**

- Dose the liquid directly prior to mixing to avoid the evaporation of ingredients.
- Powder and liquid bottles and the dosing syringe should be carefully resealed after use.

**Storage and Shelf-life**

Store the product at 15–25°C/59–77°F.

Do not use after the expiration date.

**Customer Information**

No person is authorized to provide any information which deviates from the information provided in this instruction sheet.

**Warranty**

3M ESPE warrants this product will be free from defects in material and manufacture. 3M ESPE MAKES NO OTHER WARRANTIES INCLUDING ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. User is responsible for determining the suitability of the product for user's application. If this product is defective within the warranty period, your exclusive remedy and 3M ESPE's sole obligation shall be repair or replacement of the 3M ESPE product.

**Limitation of Liability**
Except where prohibited by law, 3M ESPE will not be liable for any loss or damage arising from this product, whether direct, indirect, special, incidental or consequential, regardless of the theory asserted, including warranty, contract, negligence or strict liability.

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