



Structur 2 SC / Structur Premium

EXCELLENT TEMPORARIES WITH STRUCTUR

Structur

SIMPLE TO USE

Reliable products are required to provide your patients with prosthetic treatments, especially with regard to temporaries. Your patients' masticatory function, phonetics and aesthetics must be guaranteed until the definitive restoration is fitted.

For the patient, tooth-like temporaries are an important step on the way to the definitive treatment. Patients increasingly judge these no longer according to function only (protection of the prepared tooth and guaranteed masticatory function), but also with regard to aesthetics (matching tooth shade and shape).

Long-term temporaries frequently substitute the dental prosthesis for weeks or months when extensive restoration work is carried out. The properties of the material for temporary restorations are of particular importance in these cases.

Tooth shade and natural appearance

Temporary restorations should be "invisible"; they should be indistinguishable from natural teeth. A wide range of eight tooth shades, plus a special incisal shade, together with the natural fluorescence of Structur make it possible to fulfil these requirements.

Tooth-like hardness

Of course, the patient would like to eat with his/her temporaries. Unrestricted food intake is essential for the patient and crucially influences his/her well-being and compliance with the ongoing treatment.

High material stability

The utilized material must be able to withstand chewing pressure for weeks or months, without giving way or breaking when loaded. This applies in particular to crown margins thinning towards the edge and to other delicate structures of the restoration.

Quick and simple application

The fabrication of temporary restorations must be easy to integrate into the process of fabrication of prostheses in the surgery: Simple, bubble-free mixing even of small amounts, i.e. for a single crown, as well as larger amounts for multi-unit treatments, are not a problem even for your assistant, thanks to the presentation in the 1:1 cartridge and the tried and tested mixing tips.

Mixing errors, shade variations and mixed-in air bubbles belong to the past, thanks to the choice of presentation of either the cartridge or QuickMix syringe with the coordinated mixing tips.

Short period of retention in the mouth

Both Structur 2 SC and Structur Premium can be removed from the mouth after only 1.5 minutes, and both will have reached the solidity for finishing and polishing after just another 3 minutes.

Finishing and polishing

Temporary restorations must be finished, whether they are fabricated directly or in the laboratory. The reason is the overfilling of the impression or deep-drawing tray that is essential in the fabrication of the temporary. The resulting excess must be removed and the prepared surfaces smoothed and polished.

For the surgery, the hardness of the material used for temporary restorations is important in two aspects: On the one hand, the initial hardness must be as high as possible after setting of the material to permit finishing and incorporation of the restoration; on the other hand, the patient is unlikely to be gentle with the temporary after leaving the surgery. High stability values of materials for temporary restorations that are reached only hours after setting thus do not reflect the reality in the surgery.

All Structur products can be finished and polished with fine, cross-cut carbide burs, fine metal or sandpaper separating disks, goat hair brushes and textile buffers. Carefully polished, smooth temporaries impede adhesion of plaque and simultaneously facilitate cleaning for the patient.

Structur

QUICK FABRICATION

Tooth-like brilliance

The surfaces of temporary restorations exhibit a texture that is different from that of a natural tooth. Modern developments now also permit tooth-like brilliance and “high gloss” surfaces on temporaries. Easy Glaze, for instance, is very easy to use for this, and it achieves an immediately impressive result.

Corrections and repairs

Margins of crowns in particular, also those on temporary restorations, determine the success or failure of the prosthetic treatment to a large degree.

With subgingivally prepared teeth, there is no support of the gingiva provided by the hard substance of the tooth. If this condition is not taken into account during fabrication of the temporary, then potential consequences are gingival recession, exposed preparation boundaries and thus aesthetic complications from the now visible margin of the restoration. In the anterior region, this is not tolerated by many patients and it thus frequently leads to fabrication of a new temporary, with the associated time and financial expenditure caused to your surgery. In contrast, carefully examining the crown margins and possibly carrying out additional finishing are less costly and time-consuming.

This includes lining the crown margins with a thin layer of material that supports and shapes the gingiva in the area of the preparation boundary. The natural emergence profile on the gingival line thus remains intact. Cutting the area free from the impression before filling is too imprecise and frequently leads to overextended margins on the temporaries and thus also to the most undesirable effect of compression of the gingiva!

Structur 2 SC and Structur Premium are composites and can be seamlessly adapted to the existing conditions in a quick and efficient manner, without additional application of a bond.

The situation is the same with regard to repairs and small modifications that can arise from treatment of more complex prosthetic restorations: Crowns and bridges made from Structur 2 SC or Structur Premium can be quickly and efficiently corrected and repaired with the same material or a light-curing composite without any problems.

2-layer technique

Structur Premium QM in the “incisal” shade is available in the QuickMix syringe and is, above all, designed for restorations in the anterior region.

Structur Premium QM “Incisal” is used together with Structur 2 SC or Structur Premium to fabricate temporaries with a distinct incisal edge and tooth-like light refraction. This technique gives the temporary high transparency and natural appearance.

The effort required for the application of the 2-layer technique is minimal, and this technique can also be quickly learned and carried out by an assistant. Particularly when used in combination with Structur Premium, this creates temporary restorations that exceed your expectations and those of your patients for natural aesthetics, transparency and stability. This applies in particular to the anterior region.

Cost-effective

Structur Premium and Structur 2 SC enable you to control the cost of fabricating temporary restorations: Structur 2 SC is the right choice of product for individual crowns, short bridges in the posterior region with a short to moderate retention time in the mouth, and for cost-conscious patients. For patients with high aesthetic expectations, Structur Premium is the appropriate product for long-term temporaries and for long bridges and restorations in the anterior region. Structur Premium QM in the incisal shade can be combined with both Structur Premium and Structur 2 SC to create a natural appearance for the incisal part of the temporary restoration.

Structur 2 in the QuickMix syringe is perfect where only a small amount is required for smaller temporaries such as individual crowns or partial crowns and for corrections.

Structur

CLINICAL CASES



Source: Dr. Wentaschek, University Medical Center, Department of Prosthodontics, Johannes Gutenberg University Mainz

Severely abraded maxillary front showing spalling and defects on teeth 12-22 (Fig. 1). The situation is aesthetically and functionally unsatisfactory. The maxillary front was rebuilt from wax on the model and transferred to a silicone key (Fig. 2). Restoration of the incisors with the assistance of the composite key (Fig. 3). Preparation of teeth 12-22 for crown treatment after a new silicone impression has been taken (Fig. 4). Fabrication of the provisional restorations with Structur 2 SC. Refining of the Structur crowns with Easy Glaze (Fig. 7) and incorporation. Distinct improvement of the clinical situation through the temporaries, even before further treatment is carried out (Fig. 9).



Source: Dr. Wentaschek, University Medical Center, Department of Prosthodontics, Johannes Gutenberg University Mainz

Teeth 13, 11, 21 and 22 have been prepared for a bridge. Tooth 12 is missing; 11 has been furnished with a post build-up. Temporary restoration made from Structur Premium in situ (Fig. 1); situation without the provisional bridge (Fig. 2). Modelling of the gingiva through accurately designed bridge unit (Fig. 3) with the emergence profile of a natural tooth. Structur Premium supports such measures through its tissue-friendly properties and because it facilitates precise, appropriate finishing.

Structur

EVEN BETTER WITH EASY GLAZE



Source: Dr. Wentaschek, University Medical Center, Department of Prosthodontics, Johannes Gutenberg University Mainz

Unsatisfactory aesthetic situation due to missing incisors and mesial inclination. The black triangle between the central incisors is particularly conspicuous and unpleasing (Fig. 1). The clinical case shows the prepared and thereby axially corrected teeth (Fig. 2), as well as the treatment with temporaries made from Structur 2 SC (Fig. 3). A considerable improvement on the initial situation, with visible aesthetic benefit provided by temporaries.

Easy Glaze®

HIGH-GLOSS SURFACES

Easy Glaze – Visible brilliance

Easy Glaze is a light-curing, high-gloss, nano-filled, protective varnish for surface finishing of composite and glass ionomer materials.

Structur Premium and Structur 2 SC are composites and thus ideally suited for surface finishing with Easy Glaze: after application, the protective varnish with its perfectly coordinated viscosity penetrates the smallest recesses and anomalies, smoothes them out and leaves behind a high-gloss surface after polymerisation with a conventional polymerisation light, without polishing or modification of the occlusion.

Application couldn't be simpler, and it permanently and recognisably improves the aesthetic result: the sealed surface impedes plaque adhesion, facilitates oral hygiene and increases the quality of the temporary to an extent that is visible to your patient.

Product advantages:

- High-gloss, tooth-like surfaces
- Quick and simple application
- Immediately visible results



REF 1016 Bottle 5 ml, accessories

Provicol

UNCOMPLICATED LUTING

Provisional luting of temporaries

A material for luting temporaries must meet several criteria:

- Reliable luting
- Simple mixing and uncomplicated application
- Quick application
- Biocompatibility
- No inhibition of composite cements

The demands on temporary luting materials are contradictory. On the one hand, secure luting is expected; on the other hand, uncomplicated removal of the restoration thus incorporated must be guaranteed.

The preparation wound must be closed without irritation and the material must not impact on subsequent definitive luting carried out with frequently used luting composites. Eugenol therefore has to be ruled out as an ingredient, because it acts as an inhibitor in the polymerisation of composites. The calcium hydroxide (Ca(OH)₂) contained in Provicol has a

bacteriostatic effect and does not impact negatively on the subsequent definitive cementation of the restoration, not even if composite luting cements are used.

Provicol is available in three different types of presentation, enabling you to select the one that can best be integrated into your chairside work routine.

Provicol QM:

- Application from the QuickMix syringe
- Quick and always bubble-free mixing

Provicol C:

- Presentation in 65 g cartridges
- Ideal for many units in surgeries with a high rate of prosthodontics

Provicol:

- Manual-mix version
- The proven Provicol in tubes

Indications

Provicol / Provicol QM / Provicol C

Eugenol-free luting cement with calcium hydroxide for temporary mountings

Product advantages

- Quick application
- Contains calcium hydroxide
- Eugenol-free, thus ideal before adhesive luting with composite-cements



REF 1075 Base 25 g, catalyst 25 g
REF 1078 Bulk pack
5 x 2 tubes base and catalyst, 25 g each

REF 1074 QuickMix syringe 5 ml,
mixing tips type 10

REF 1076 Cartridges 2 x 65 g,
mixing tips type 6

Structur Premium

PROVISIONAL CROWN AND BRIDGE MATERIAL FOR THE HIGHEST STANDARDS

Indications

Structur Premium

All-purpose crown and bridge material for the highest aesthetic standards for temporary crowns, bridges, veneers, inlays and onlays

Structur Premium QM Incisal

Special incisal shade for the fabrication of temporaries in the 2-layer technique

Product advantages

- Quick and easy handling
- Surgery-oriented and hygienic application
- Short setting time in the mouth
- High initial *and* end hardness
- Excellent shape and shade stability
- Durable temporaries
- Outstanding aesthetics
- Wide range of shades
- Natural fluorescence
- Simple to finish and polish
- 2-layer technique through special incisal shade



Structur Premium

REF 1710	Intro Set cartridge 75 g A3, dispenser type 2, mixing tips type 6
REF 1712	Cartridge 75 g A1, mixing tips type 6
REF 1713	Cartridge 75 g A2, mixing tips type 6
REF 1714	Cartridge 75 g A3, mixing tips type 6
REF 1717	Cartridge 75 g A3.5, mixing tips type 6
REF 1718	Cartridge 75 g B1, mixing tips type 6
REF 1715	Cartridge 75 g B3, mixing tips type 6
REF 1716	Cartridge 75 g BL, mixing tips type 6
REF 1719	Bulk pack cartridges 5 x 75 g A3, mixing tips type 6



Structur Premium QM

REF 1726	QuickMix syringe 8 g I (incisal), mixing tips type 10
----------	---

Structur 2 SC

PROVISIONAL CROWN AND BRIDGE MATERIAL THAT MEETS ALL DEMANDS

Indications

Structur 2 SC

All-purpose, economical crown and bridge material for temporary crowns, inlays and onlays

Structur 2 QM

All-purpose, economical crown and bridge material for temporary crowns, inlays and onlays, presented in the QuickMix syringe

Product advantages

- Quick and easy handling
- Economical and hygienic application
- Short setting time in the mouth
- Good shape and shade stability
- Excellent aesthetics
- Wide range of shades
- Natural fluorescence
- Simple to finish and polish



Structur 2 SC

- REF 1478 Intro Set cartridge 75 g A3, dispenser type 2, mixing tips type 6
- REF 1479 Cartridge 75 g A1, mixing tip type 6
- REF 1480 Cartridge 75 g A2, mixing tip type 6
- REF 1481 Cartridge 75 g A3, mixing tip type 6
- REF 1490 Cartridge 75 g A3.5, mixing tip type 6
- REF 1482 Cartridge 75 g B1, mixing tip type 6
- REF 1483 Cartridge 75 g B3, mixing tip type 6
- REF 1486 Cartridge 75 g C2, mixing tip type 6
- REF 1485 Cartridge 75 g BL, mixing tip type 6
- REF 1492 Bulk pack cartridges 5 x 75 g A2, mixing tips type 6
- REF 1491 Bulk pack cartridges 5 x 75 g A3, mixing tips type 6



Structur 2 QM

- QuickMix syringe 8 g, mixing tips type 10
- REF 1309 A1
- REF 1310 A2
- REF 1311 A3
- REF 1312 B1

VOCO GmbH
P.O. Box 767
27457 Cuxhaven
Germany

Tel.: +49 (0)4721-719-0
Fax: +49 (0)4721-719-140

info@voco.com
www.voco.com

Available from: