

Instructions /Applications for Use

Night Guard Fabrication Application:



1 Heats at 135° F and turns clear.



2 Form as anatomy requires with a wet gloved hand.



3 Air for 10 seconds.



4 Custom fit chairside in seconds.



5 Trim with slow speed if required.

Instructions /Applications for Use

Quick and Flexible Bite Registration Application:




photos credit: Dr. Greg Gallant

1 Heat the **Temp Tabs - True Blue** or **Temp Tabs - Natural** (as in night guard). Place the Tab or Tabs over the units selected by the doctor.

2 Have the patient close. This will give the doctor centric occlusion.*

3 To get centric relation the doctor may modify the technique by using an anterior dis-occlusion device (Lucia jig - also using **Temp Tabs**) and then guiding the patient to centric relation.



Two Temp Tabs were used in this demonstration, but one would have been sufficient.

4 **Temp Tabs** may be cooled with air/water.

Other Uses for the Temp Tabs

- Implant Placement Stent
- Matrix for Temporary Material to Create Provisionals
- Implant Abutment Jig
- Implant Surgical Site Splint to Prevent Extrusion of Adjacent Units



* Draw down a small amount of softened Tab on the buccal side interproximal for an instant placement guide!