Biological Monitoring Service

Mail-In

Contents of the Biological Monitoring Service

This Biological Monitoring Service contains test envelopes. Within each envelope are test strip(s) and one (1) control strip.

Instructions For Use

1. Remove test strip(s) from the pocket marked "test strip." (DO NOT remove test strip(s) from the glassine pouch, as this will invalidate the test. (DO NOT REMOVE THE CONTROL STRIP.)
2. Place one (1) test strip in the center of the packs to be sterilized, and the other test strip (if applicable) in the most difficult to sterilize location. If the sterilizer manufacturer identifies another location in the chamber for placement of the test strip, follow those directions. Process as part of a normal load and run a normal sterilization cycle. DO NOT PROCESS THE CONTROL STRIP.
3. Upon completion of the sterilization cycle, remove the test strip(s) and return to the pocket marked "test strip." Seal the flap closed. **Note: Let the strips dry completely before returning to the envelope.**
4. Complete all information requested on the envelope and mail.
NOTE: NOT RESPONSIBLE FOR ENVELOPES NOT PROCESSED DUE TO INCOMPLETE INFORMATION.
5. Test results will be e-mailed or faxed to your office quarterly. Test failures will be reported immediately by e-mail or fax, with a follow-up phone call to verify receipt of the report. Confidential test results are available online at crosstexbms.com.

Test Strips

Bacterial spores test strips are recommended for monitoring steam, dry heat, chemical vapor, and ethylene oxide gas sterilization processes. These spore strips are certified to meet the following standards and guidelines: • Military Specification MIL-S-36586A • Association for the Advancement of Medical Instrumentation, Biological Indicators for Saturated Steam Sterilization Processes in Healthcare Facilities, ANSI Approval Date: March 17, 1986; Designated: ANSI/AAMI ST-19-1994.

Performance Certification

This product can be used for the following sterilization methods. Performance data is available upon request.

- Steam (Gravity) • Steam (Pre-Vac) • Chemical Vapor • Dry Heat • Ethylene Oxide

Patterson Mail-In biological monitoring service provides increased testing accuracy, with third-party documentation. Patterson Mail-In biological monitoring service test results are available online at crosstexbms.com.

Processing Tips

- Biological indicators should be placed into the densest sterilization package located in the coolest spot in the sterilizer chamber.
- Paper/poly pouches should be placed on edge and with sufficient space between the packages to allow the sterilizing agent to contact every surface of articles to be sterilized.
- Use the size of sterilization package that is appropriate to the dimensions of the device or number of instruments to be processed.
- Open all hinged instruments to permit the sterilization agent to contact all surfaces.
- Make sure that your sterilizer is being operated using the correct times, temperatures, and pressures.
- Do not overload the sterilizer chamber.
- Do not interrupt the sterilization cycle while in process.
- Allow packages and biological indicators to dry completely before handling.

Tests are processed by the **Biological Monitoring Service**.

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