

DIRECTIONS FOR USE:

Use fresh pdCARE™ Dual Enzyme Detergent for each set of instruments.
Discard diluted pdCare Dual Enzyme Detergent after each use.

MANUAL CLEANING: Add 1oz. (1 pump yields 1oz.) of concentrate to one gallon of warm water (68°F-104°F). Soak instruments and equipment immediately after use, until soil is dissolved and removed. Soak for a minimum of 1 minute. Soak longer if matter is dried on, until visibly clean. Immerse instruments completely and aspirate solution into any channels and lumens. Use an instrument brush for light mechanical cleaning if necessary. Following soak period, visually inspect for cleanliness and rinse thoroughly with warm water. Rough dry. After cleaning with pdCare Dual Enzyme Detergent, equipment and instruments must be disinfected/sterilized according to appropriate procedures.

ULTRASONIC UNITS: Add 1 oz. of concentrate per gallon of water. Recommended water temperature is 68°F-104°F (20°C-40°C). A soak time of 2 minutes or more is recommended for ultrasonic cleaners. Process instruments according to manufacturer's directions.

WASHER STERILIZERS/DECONTAMINATORS: Attach concentrate container to the unit's automatic dispensing unit.

Note: pdCare Dual Enzyme Detergent is not a sterilant/disinfectant.

PRECAUTIONS:

Protective eye gear, gloves and clothing MUST BE WORN, especially where exposure to blood and body fluid is expected. If accidental eye contact with pdCare Dual Enzyme Detergent concentrate occurs, flush immediately with water for 15 minutes and seek medical attention. Temporary redness of the eye may occur.


PATTERSON
DENTAL

pd CARE™
Dual Enzyme Detergent

Pre-Cleaning Detergent Concentrate

CLEANS FAST: In manual, ultrasonic and automatic cleaning systems, pdCare Dual Enzyme Detergent rapidly dissolves and carries away organic soils.

IMPORTANT: Store in a cool place.

**Net Contents: One Gallon
(128 fluid oz.) (3.785 L)**

Manufactured for:
Fabriqué pour :


PATTERSON
DENTAL

Patterson Dental Supply, Inc.
1031 Mendota Heights Road
Saint Paul, MN 55120

Low foaming, dual enzymatic detergent designed for easy one-step pre-cleaning by disintegrating blood and protein from delicate instruments and equipment.



+D6620708565001

Reorder Number
N° de produit

085-5650

MADE IN U.S.A. - Fabriqué aux É.-U.

FEATURES:

pdCare Dual Enzyme Detergent is a low-foaming, dual enzymatic detergent that effectively cleans away blood, tissue, mucous and other protein-rich body fluids. pdCare Dual Enzyme Detergent is designed for use in Dental, Surgery, Central Processing, or other areas where instruments and equipment may be soiled with organic and inorganic debris. pdCare Dual Enzyme Detergent contains two different proteolytic enzymes that work in a broad range of temperatures and pH conditions. pdCare Dual Enzyme Detergent contains low foaming surfactants that efficiently remove lipids, starches, inorganic debris and prevent debris from redepositing onto instruments. pdCare Dual Enzyme Detergent also contains corrosion inhibitors to protect a broad range of metals from corrosion.

SAFE FOR USERS: pdCare Dual Enzyme Detergent reduces the need for manual handling and scrubbing, especially in hard-to-reach places like channels, lumens and instrument box locks.

SAFE FOR INSTRUMENTS AND EQUIPMENT: pdCare Dual Enzyme Detergent is a neutral pH formula that provides advanced corrosion protection that is safe for rubber, plastic, stainless steel, carbon steel and soft metals.

STORAGE AND DISPOSAL:

Storage: Store in a cool area.

Solution Disposal: For unused solution, flush thoroughly with large quantities of water into sewage disposal system in accordance with federal, state, and local regulations.

For used solution, the waste solution must be characterized by the generator and disposed of in accordance with federal, state and local regulations.

Container Disposal: Do not reuse empty container. Wrap container and put in trash.