Mise en garde :

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Patterson Dental Supply, Inc.
1031 Mendota Heights Road
Saint Paul, MN 55120

PATTERSON® TEMPORARY CROWN AND BRIDGE MATERIAL
Carefully read instructions prior to use

PATTERSON® TEMPORARY CROWN AND BRIDGE MATERIAL is a fluorescent cold polymerizing paste-paste system for the production of temporary crowns, bridges, inlays and onlays. PATTERSON® TEMPORARY CROWN AND BRIDGE MATERIAL consists of base paste and catalyst paste. PATTERSON® TEMPORARY CROWN AND BRIDGE MATERIAL is supplied in 3 shades A1, A2 and A3.

Indications:

PATTERSON® TEMPORARY CROWN AND BRIDGE MATERIAL is indicated for the fabrication of temporary crowns, bridges, inlays and onlays.

Contraindications:

PATTERSON® TEMPORARY CROWN AND BRIDGE MATERIAL contains methacrylates, amines, terpenes, benzoyl peroxide, BHT. PATTERSON® TEMPORARY CROWN AND BRIDGE MATERIAL should not be used in cases of known hypersensitivities (allergies) to any of these ingredients.

Instructions:

Insert the safety cartridge into a type 2 dispenser or a compatible dispenser with a suitable 4:1 style plunger. Before using a new cartridge, remove the cartridge cap and dispense a small amount of material, just until the material flows evenly from both chambers. Discard this material. Insert a type 6 mixing tip and lock with a 90 degree clockwise turn. Gently pressing the dispenser trigger will extrude the material and automatically mix it in the correct ratio. Store the cartridge with the used mix tip - do not remove the used mixing tip from the cartridge after each use. Replace the used mixing tip with a new tip just before each application. Before putting the new mix tip on, visually check the cartridge opening to ensure that the material can flow.

Application:

1. In cases where no prefabricated crowns or similar restorations are used, prepare an impression by leveling undercuts and cutting drainage channels. If the wall of the temporary restoration is too thin, enlarge the impression, especially removing septums in the interdental embrasures. Then clean and air dry the impression.
2. Apply PATTERSON® TEMPORARY CROWN AND BRIDGE MATERIAL directly into the impression. Do not overfill. PATTERSON® TEMPORARY CROWN AND BRIDGE MATERIAL attains an elastic phase after 1:00 - 1:30 minutes from beginning of mixing (approx. 0:30 - 1 minute after inserting into mouth).
3. Remove the temporary from the mouth at least 1:30 minutes (elastic phase) after mixing, to ensure a trouble-free removal. The setting process can be monitored by checking the excess material in the mouth. The setting time of non-used material (in mixing tip) cannot be compared with the setting time of the material in the mouth. At a room temperature of 72°F (23°C) and 50% humidity working times are as follows:

minutes:s	minutes:s	minutes:s	minutes:s	minutes:s
0:00	0:30	1:00	1:30	4:00
application in the mouth	setting in the mouth	removal	complete setting	finishing

Slight differences may result in different temperatures or humidity levels.

4. After approx. 2-3 minutes remove the temporary from the impression and wipe dry. Consider using a solvent, i.e., ethyl alcohol to remove the thin oxygen inhibition layer and allow for easier finishing. The surface of the temporary restoration should appear silky-matte finish.
5. Carefully trim any excess material and/or rough areas near the gingiva.

6. Finishing should only be started after the material sets completely (approx. 5 minutes after mixing).
7. Cement the temporary restoration using an eugenol-free temporary cement.

Caution:

- Before making the temporary, the prepared cores should be insulated for sensitivity using a fluoride varnish, (PATTERSON® 5% Sodium Fluoride Varnish) or other suitable products.
- In order to avoid fractures or injuries from debris, the temporary should not be exposed to excessive loads especially in the molar areas.
- Linings and cores built up with composite (resin-based) materials should be sufficiently insulated with glycerine for easy separation of the temporary from composite resin material.
- Inadequate hygiene or intensive use of some mouthwashes may lead to discoloration.
- Finishing should only be started after complete setting of the material (approx. 4 min after mixing).

Repairs:

If the temporary work breaks despite PATTERSON® TEMPORARY CROWN AND BRIDGE MATERIAL's extreme stability, the following procedures are recommended:

- Air bubbles: fill with PATTERSON® TEMPORARY CROWN AND BRIDGE MATERIAL
- In case of fracture immediately after preparing the temporary: re-connect the fracture with freshly mixed PATTERSON® TEMPORARY CROWN AND BRIDGE MATERIAL
- Restoration of already used temporaries: Roughen the fracture with a bur and place undercuts. Apply a bonding agent to the prepared fractures and cure according to the respective instructions for use. Apply a flowable composite to the prepared areas. Bond the fractured pieces together and cure according to the instructions for use.

Storage:

Store at temperatures between 39°F - 73°F (4°C - 23°C). Do not use after expiration date.

Do not expose PATTERSON® TEMPORARY CROWN AND BRIDGE MATERIAL to direct sunlight or the chairlight.

Product ordering information

Pack Safety Cartridge, 6 mixing tips
 Shade A1 Art. no.: 090-2213
 Shade A2 Art. no.: 090-2221
 Shade A3 Art. no.: 090-2239

This material has been developed solely for use in dentistry. Processing should be done strictly according to the instructions for use. Patterson Dental Supply, Inc. recognizes its responsibility to replace products if proven to be defective. Patterson Dental Supply, Inc. does not accept liability for any damage or loss, directly or indirectly, stemming from the use of or inability to use the products described. Before using, it is the responsibility of the user to determine the suitability of the product for its intended use. The user assumes all risk and liability in connection therewith. Descriptions and data constitute no warranty of attributes and are not binding.

CAUTION: U.S. Federal Laws restrict this device to sale by or on the order of a dentist.

No person is authorized to provide any information that deviates from the information provided in the instructions for use. For questions or comments, please call 1-800-873-7683.

Keep this material out of reach of children.
 For dental use only.

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