# Kerr™

**Composite Material**

**Herculite™ *Ultra Flow***  
Nanohybrid Flowable Composite

**Instructions For Use**

Herculite™ Ultra Flow is a flowable medium viscosity light cure nanohybrid resin restorative suitable for use as a base/ liner for Class I and Class II restorations, for Class III restorations, Class V restorations, and small non-stress bearing restorations. Additional indications include: pit and fissure sealant, repair of enamel defects, repair of porcelain restorations, minor occlusal build-ups in non-stress bearing areas, and incisal abrasions.

## **Instructions for Use**

### **Preparation**

1. Isolate area using a rubber dam. Dry the preparation.
2. Following the manufacturer's instructions, apply a bonding system for either total-etch or self-etch using your preferred OptiBond™ product.
3. Select the appropriate shade of flowable composite.
4. Remove the syringe cap and save. Twist the disposable dispensing tip securely onto the syringe.

### **Placement**

5. Dispense a layer of Herculite Ultra Flow directly onto the preparation along the axial and gingival margins. After extruding material onto preparation, pull back syringe plunger to relieve pressure on the material.
6. Use a brushing action to allow the material to flow into the marginal interfaces and recesses of the preparation.
7. Incremental layers of Herculite Ultra Flow should be no more than 2 mm. Light cure each increment\*.
8. If using Herculite Ultra Flow in a base/liner application, place a restorative composite, such as Herculite Ultra,

directly over the cured Herculite Ultra Flow according to the manufacturer's instructions.

9. When using Herculite Ultra Flow as the final restorative product, finish and polish using your preferred polishing system. Kerr recommends products from Axis Dental.

### **Shading System**

Herculite Ultra Flow is available in shades that match the Vita Shade Guide.

For a list of available shades, please visit

<http://www.kerrdental.com>

### **Notes:**

1. Replace syringe cap after each use to prevent the resin from hardening in the syringe. Syringe dispensing tips are intended for single patient use only, to prevent cross contamination between patients.
2. Refer to outer package for expiration date. Do not use after expiration date.
3. Protective eyewear for patients and staff is recommended.

### **Storage Conditions:**

Store at ambient temperature.

### \*Recommended Cure Times:

Use a curing light with output equal to or greater than 800 mW/cm<sup>2</sup> such as Demi™ Ultra, Demi™ Plus, Demi™, or Optilux™ 501.

	<b>Shades A3 and lighter</b>	<b>Shades A3.5 and darker (except Universal Opaque)</b>	<b>Universal Opaque</b>
<b>Demi Ultra / Demi Plus / Demi</b>	5 seconds	10 seconds	20 seconds
<b>Optilux 501 (Boost Mode)</b>	10 seconds	10 seconds	40 seconds
<b>(Ramp Mode)</b>	20 seconds	20 seconds	40 seconds
<b>(Regular Mode)</b>	20 seconds	20 seconds	40 seconds

For all other lights, see manufacturer's recommendation.

**CAUTION:** Uncured methacrylate resin may cause contact dermatitis and damage the pulp. Avoid contact with skin, eyes and soft tissues. Wash thoroughly with water after contact.

### Limited Warranty - Limitation of Kerr's Liability

Kerr's technical advice, whether verbal or in writing, is designed to assist dentists in using Kerr's product. The dentist assumes all risk and liability for damages arising out of the improper use of Kerr's product. In the event of a defect in material or workmanship, Kerr's liability is limited, at Kerr's option, to replacement of the defective product or part thereof, or reimbursement of the actual cost of the defective product. In order to take advantage of this limited warranty, the defective product must be returned to Kerr. In no event shall Kerr be liable for any indirect, incidental, or consequential damages.

EXCEPT AS EXPRESSLY PROVIDED ABOVE, THERE ARE NO WARRANTIES, BY KERR, EXPRESS OR IMPLIED, INCLUDING WARRANTIES WITH RESPECT TO DESCRIPTION, QUALITY, OR FITNESS FOR A PARTICULAR PURPOSE.

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are registered trademarks in the
## **Kerr Corporation**
Orange, CA 92867; U.S.A.
## **Kerr Distribution Facility: Kerr Australia Pty. Ltd.**
owned by Kerr Corporation.