VINYL POLYSILOXANE IMPRESSION MATERIAL

380 ml CARTRIDGE SYSTEM

INSTRUCTIONS FOR USE:

The Monophase and Heavy Body addition cured silicone elastomeric impression materials in the 380 ml cartridge meet or exceed all requirements for dental impression materials as specified by the American Dental Association Specification No. 19. The enhanced hydrophilic properties improve the exact reproduction of the subgingival areas and other moist surfaces. Both are formulated to be thixotropic, which allows them to stay where placed without running or sagging, yet flow under pressure for easy application in the oral cavity around the preparation and gingival area. The Monophase and Heavy Body are suitable for crowns, bridges, inlays, partial and complete dentures and denture repairs. Monophase and Heavy Body can also be used for relining functional impressions, copper bands and prosthetic restorations and are compatible with all other VPS systems.

PREPARATION FOR USE:

1. Position the cartridge, nozzle side up, with the nose cap tab facing you. While firmly holding the cartridge with both hands, using your thumb, gently push on the tab until the nose cap unlocks. Using your thumb and forefinger slowly lift and remove the nose cap. Once free, discard the nose cap. Do not reuse the nose cap.
2. Load the 380 ml cartridge into the dispenser following the manufacturer's direction.
3. Following the operating instructions supplied by the machine manufacturer, bleed the cartridge by extruding a small amount of material until the base and catalyst are evenly dispensing. Wipe the extruded materials from the end of the cartridge using an upward movement to avoid cross-contamination. This is recommended for first-time use and when the cartridge has been stored upright for a prolonged period of time.
4. Attach the dynamic mixing tip to the extrusion nozzle by aligning the holes of the cartridge with the corresponding mixer ports. Secure the tip to the cartridge using a bayonet locking ring (to lock, rotate ring clockwise approximately a ¼ turn). Note: only use Sulzer/Mixpac dynamic mixers as other manufacturers' tips will not be compatible with this cartridge system.
5. Following the operating instructions supplied by the machine manufacturer, start dispensing material discarding the first half inch of product. Immediately after dispensing the required amount of material, if the dispensing machine does not automatically release pressure from the plungers, manually move the plungers back, relieving pressure on the pistons, thus avoiding excess discharge of material.
6. During the workday, the cartridge may be stored in the machine with the used mixer attached. To start a new impression, remove the locking ring (to unlock, rotate ring counterclockwise approximately ¼ turn), then remove and discard the used mixer. Attach the new dynamic mixer starting at step #4.
7. At the end of the workday or when the cartridges will not be used for an extended period of time, it is recommend that the cartridge be removed from the machine and stored in an upright position.

DIRECTIONS FOR USE:

1. Prepare the teeth in the usual manner for impression taking.
2. In any technique where an impression tray is required, the first step is to coat the tissue side of the preselected tray with vinyl adhesive. The optimum bond occurs when the adhesive film is applied as thinly and evenly as possible. You can use either a perforated or nonperforated metal or plastic impression tray.
3. To load the impression tray, follow the operating steps in the "preparation for use" section above. Fill the tray with the desired amount of material, holding the tray in a tilted position while keeping the tip immersed in the material to avoid air bubbles. Seat the filled impression tray into the mouth passively, maintaining the position until set.
4. When fully set, remove the impression in an axial direction without any rocking motion. The impression can now be gently rinsed and disinfected.
5. Impressions can be poured after 60 minutes or up to 60 days later.

ADDITIONAL INFORMATION:

- As with all VPS systems, avoid contamination with latex gloves, certain retraction agents, or any material that could inhibit the set of the VPS. Do not use with condensation silicones.
- Impressions can be disinfected with any standard disinfecting solution.
- Impressions can be silver or copper plated.
- Multiple pours are acceptable if care is taken.
- Shelf life is two years from date of manufacture when stored at room temperature (74°F or 23°C).

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