VINYL POLYSILOXANE IMPRESSION MATERIAL

DIRECTIONS FOR USE

INTRODUCTION

Addition cured silicone elastomeric impression material (VPS) meets or exceeds requirements for Type 1 impression materials specified by the American Dental Association Specification No. 19. Products are available in multiple viscosities (and as a putty) in both tubes and cartridges for use in various impression techniques.

VPS materials are manufactured to perform according to their specifications for work and set time at 73°F (23°C). Higher temperatures will decrease the work time and set time. Conversely, lower temperatures will extend work/set times. The shelf life of these impression materials is two years. During hot weather VPS materials can be refrigerated to prevent degradation – but should be warmed to room temperature before use.

In techniques where an impression tray is required, the first step is to paint the tissue side of the preselected tray with vinyl adhesive. The optimum bond occurs when the adhesive film is as thin and even as possible. The adhesive should be allowed to dry for 5 minutes.

INSTRUCTIONS FOR USE OF CARTRIDGES

Important: Before attaching a mixing tip to the cartridge follow these steps:

Note: Double barrel cartridges require no hand mixing. Working time begins as the base and catalyst are extruded from the mixing spout.

1. Load the cartridge into the dispensing gun.
2. Remove the plastic cap from the cartridge.
3. After a cartridge is properly positioned in the dispensing gun and the cap is removed, a small amount should be extruded to make certain that the product is exiting the barrels evenly and that there is no plugging.
4. Attach a mixing tip to the cartridge with a 90° rotation. When the dispensing gun trigger is squeezed, the cartridge material will be mixed as it flows through the mixer. The mixing tip can be used directly in the mouth or in a tray or to backfill a syringe. In addition, an intraoral tip can be snapped into the mixer to create a very fine bead of material.
5. If either an intraoral syringe tip attached to the mixer or a single barrel back loaded syringe is being used, the silicone materials should be flowing in front of the syringe tip to minimize the entrapment of air.
6. Once fully set, remove the impression in an axial direction without any rocking motion. The impression can be gently rinsed and then disinfected.
7. Do not remove the used mixer; leave it in place on the cartridge where it acts as a seal for the remaining material. The next time the cartridge is used, the old mixer can be removed and discarded just prior to attaching the new mixer.

NOTES

Vinyl impression materials will not work properly and may never set when sulfur, eugenol or amines are present. Latex examination gloves should not be used as they can leave a residue on the teeth or soft tissue sufficient to inhibit the set of the VPS and should never be used to mix VPS putty. Certain gingival retraction cords containing sulfides can also cause severe inhibition particularly in the sulcus areas.

VPS adhesives are flammable and may contain toxic, volatile materials. DO NOT inhale and DO NOT use around open flames. Always close the adhesive bottle immediately after use. VPS impression material is nontoxic, but will stain clothing on contact.

Federal U.S. law restricts this device to sale by or on the order of a licensed dentist.

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