# VPS PUTTY

## DIRECTIONS FOR USE

### INTRODUCTION

Addition cured VPS materials are manufactured to perform according to their specifications for work and set time at 73°F (32°C). Higher temperatures will decrease the work times and conversely, materials that are lower than 73°F will have extended work times. The shelf life of VPS impression materials is two years, from date of manufacturing. During hot summer months, VPS impression materials can be stored in the refrigerator to prevent degradation — but should be allowed to warm to room temperature before using.

In any technique where an impression tray is required, the first step is to paint the tissue side of the preselected tray with vinyl adhesive. The optimum bond occurs when the adhesive film is as thin and even as possible. The adhesive should be allowed to dry for 5 minutes or per manufacturer's instructions.

### DIRECTIONS FOR USE

Preliminary impression material for use in a two-step (Putty-Wash) technique.

### INSTRUCTIONS FOR TWO-STEP TECHNIQUE

1. Dispense equal amounts of base and catalyst putty with the color-coded scoops provided. Be careful not to interchange scoops or jar caps.
2. Incorporate the two components by hand, kneading until uniform color is achieved. This should be accomplished in 30 seconds or less.
3. Place mixed putty into the prescheduled and adhesive-coated tray. Insert into the mouth and hold passively in position until set. A polyethylene space can be used to create space if desired.
4. When set, the impression should be removed without any rocking motion. The ADA recommends that putty impressions not be used as the final impression.
5. Space must be created for the final wash and this can be accomplished by removing all interproximal material as well as at least 1 mm of cured putty completely around all preparation and adjacent teeth. A vulcanite bur is ideal for cutting putty and creating the necessary space.
6. The primary impression that has been relieved is now ready for the defining wash that will be the final impression.

Type	Work Time*	Total Set Time*
Fast Set	1:00 minute	3:00 minutes
Regular Set	2:00 minutes	4:00 minutes

\*From start of mix

### NOTES

VPS impression material will not work properly and may never set when sulfur, eugenol or amines are present. For example, latex examination gloves — if placed in contact with the teeth or soft tissue or if used to mix putty will probably inhibit the set of the vinyl. Certain gingival retraction cords which contain sulfides can cause sulcus areas to remain unset.

All vinyl adhesives are flammable and contain toxic, volatile materials. DO NOT inhale and DO NOT have any open flames near the adhesive during application to the tray. Always close the adhesive bottle immediately after use.

VPS impression materials are completely nontoxic. However, silicone can stain clothing and when set is chemically resistant.

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