CEREC Cercon 4D™ Abutment Block

Multidimensional Zirconia

 **INSTRUCTIONS FOR USE**

 **MODE D'EMPLOI**

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CEREC Cercon 4D™ Abutment Block

Multidimensional Zirconia

Instructions for Use – ENGLISH

Caution: This is a medical device. For dental use only.
 USA: Rx Only

EN

1. PRODUCT DESCRIPTION

CEREC Cercon 4D™ Abutment Block, Multidimensional Zirconia, for CEREC® are part of the CEREC Cercon 4D™ Abutment System suitable for the fabrication of 2-piece screw retained implant supported hybrid restorations when cemented to prefabricated TiBase structures. CEREC Cercon 4D™ Abutment Blocks can be used to design a zirconia ceramic structure such as a Mesostructure for a hybrid abutment to support a cemented restoration or to fabricate a full contour crown for a fully anatomical hybrid abutment crown final restoration in the posterior and anterior region using a CAD/CAM procedure. CEREC Cercon 4D™ Abutment Blocks are provided with a pre-drilled screw access channel and anti-rotation feature. Indirect restorations are fabricated by grinding CEREC Cercon 4D™ Abutment Blocks using a Dentsply Sirona CAD/CAM system. Note that CEREC® Version 5.2.7 or higher is required.

The CEREC Cercon 4D™ Abutment Blocks are provided in a partially sintered state, then milled enlarged by the CEREC® CAD/CAM system. Restorations are individually processed to specification, and finally, densely sintered in the CEREC SpeedFire® Sintering Furnace. The software determines the sintering program according to the indication.

Restorations made from CEREC Cercon 4D™ Abutment Blocks have to be sintered in dry conditions. The use of stains and/or glaze is optional; however, polishing is required. The CEREC SpeedFire® furnace offers a program with a pre-drying function.

Polished or glazed, fired restorations are cemented with adhesive cement (available separately) to the appropriate Dentsply Sirona TiBase (available separately) for the implant.

CEREC Cercon 4D™ Abutment Blocks are a Type II dental ceramic material, Class 5 pursuant to ISO standard 6872:2015 + Amd.1:2018 (CTE: $10.3 \times 10^{-6}/K$ (25–500 °C)).

1.1 Intended Purpose

Ceramics for fixed dental prosthetic restorations.

1.2 Indications

CEREC Cercon 4D™ Abutment System is intended for use in partially or fully edentulous mandibles and maxillae in support of single cement-retained restorations.

The system comprises three parts:

- CEREC Cercon 4D™ Abutment Block
- TiBase
- CAD/CAM system

The CEREC Cercon 4D™ ceramic structure cemented to the TiBase is recommended for two-piece hybrid abutments for single tooth restorations and hybrid abutment crowns, used in conjunction with endosseous dental implants.

Implant Systems:

- Dentsply Sirona: PrimeTaper EV, OmniTaper EV, AstraTech OsseoSpeed TX, Frialit/XiVE, AstraTech Implant EV, Ankylos
- BioHorizons: Internal connection
- Nobel Biocare: Replace, Replace Select, Nobel Active, NobelReplace Conical Connection, Bränemark, NobelSpeedy Groovy
- Straumann: Tissue Level, Bone Level
- Thommen Medical: Element, Contact Osstem/Hiossen: Osstem TS, USA: Hiossen ET
- Zimmer/Biomet: External hex, Certain, Tapered Screw-Vent
- MIS: C1 Conical Connection, V3 Conical Connection, SEVEN internal hex, M4 internal hex
- Altatec – Camlog

CAD/CAM Systems:

- Sirona Dental CAD/CAM System

Compatible TiBase and Implant Systems**

Implant System		Titanium Base		
Manufacturer/line	Platform	Name	REF	Size
Dentsply Sirona Implants				
Primetaper EV	S	TiBase AT EV 3.6 GH 1 S	6586312	S
	M	TiBase AT EV 4.2 GH 1 L	6586320	L
	L	TiBase AT EV 4.8 GH 1 L	6586338	L
Omnitaper EV	XS	TiBase AT EV 3.0 GH 1 S*	6586304	S
	S	TiBase AT EV 3.6 GH 1 S*	6586312	S
	M	TiBase AT EV 4.2 GH 1 L*	6586320	L
	L	TiBase AT EV 4.8 GH 1 L*	6586338	L
	XL	TiBase AT EV 5.4 GH 1 L*	6586346	L

Implant System		Titanium Base		
Manufacturer/line	Platform	Name	REF	Size
AstraTech OsseoSpeed TX	3.0	TiBase AT TX 3.0 GH 1 S*	6598085	S
	3.5/4.0	TiBase AT TX 3.5/4.0 GH 1 L	6598093	L
	4.5/5.0	TiBase AT TX 4.5/5.0 GH 1 L	6598101	L
	3.5/4.0	TiBase AT OS 3.5/4.0 GH 1 L	6282532	L
	4.5/5.0	TiBase AT OS 4.5/5.0 GH 1 L	6282540	L
Xive	3.4	TiBase XiVE 3.4 GH 1 S*	6597327	S
	3.8	TiBase XiVE 3.8 GH 1 S*	6597335	S
	4.5	TiBase XiVE 4.5 GH 1 L*	6597343	L
	5.5	TiBase XiVE 5.5 GH 1 L*	6597350	L
Frialit/Xive	3.4	TiBase FX 3.4 GH 1 S	6282433	S
	3.8	TiBase FX 3.8 GH 1 S	6282441	S
	4.5	TiBase FX 4.5 GH 1 L	6282458	L
	5.5	TiBase FX 5.5 GH 1 L	6282466	L
AstraTech Implant EV	XL	TiBase AT EV 5.4 GH 1 L	6586346	L
Ankylos	C/X	TiBase ANK C/GH 1 S	6586528	S
		TiBase ANK C/GH 2 S	6586536	S
		TiBase ANK/X GH 1 S	6586544	S
		TiBase ANK/X GH 2 S	6586551	S
BioHorizons				
Internal connection	3.0	TiBase BH 3.0 GH 1 S	6532779	S
	3.5	TiBase BH 3.5 GH 1 L	6532894	L
	4.5	TiBase BH 4.5 GH 1 L	6532951	L
	5.7	TiBase BH 5.7 GH 1 L	6536242	L
Nobel Biocare				
Replace, Replace Select	NP	TiBase NB RS 3.5 GH 1 L	6282474	L
	RP	TiBase NB RS 4.3 GH 1 L	6282482	L
	WP	TiBase NB RS 5.0 GH 1 L	6282490	L
	6.0	TiBase NB RS 6.0 GH 1 L	6282508	L
Nobel Active	NP	TiBase NB A 4.5 GH 1 L	6308188	L
NobelReplace Conical Connection	RP	TiBase NB A 5.0 GH 1 L	6308253	L
Bränemark	NP	TiBase NB B 3.4 GH 1 L	6282516	L
NobelSpeedy Groovy	RP	TiBase NB B 4.1 GH 1 L	6282524	L
Straumann				
Tissue Level	RN (4.8 mm)	TiBase SSO 4.8 GH 1 L	6284249	L
	WN (6.5 mm)	TiBase SSO 6.5 GH 1 L	6284256	L
Bone Level	NC (3.3 mm)	TiBase S BL 3.3 GH 1 L	6308154	L
	RC (4.1 mm/ 4.8 mm)	TiBase S BL 4.1 GH 1 L	6308337	L
Thommen Medical				
Element, Contact	3.5	TiBase TM 3.5 GH 1 S	6531854	S
	4	TiBase TM 4 GH 1 S	6532829	S
	4.5	TiBase TM 4.5 GH 1 S	6532837	S
	5	TiBase TM 5 GH 1 S	6544360	S
	6	TiBase TM 6 GH 1 S	6544378	S
Osstem/Hiossen				
Osstem TS	Mini	TiBase O TS 3.5 GH 1 L	6527035	L
(USA: Hiossen ET)	Regular	TiBase O TS 4.0 GH 1 L	6527043	L
Zimmer/Biomet				
External hex	3.4	TiBase B O 3.4 GH 1 L	6282557	L
	4.1	TiBase B O 4.1 GH 1 L	6282565	L
	5.0	TiBase B O 5.0 GH 1 L	6282573	L

Implant System		Titanium Base			
Manufacturer/line	Platform	Name		REF	Size
Certain	3.4	TiBase B C 3.4 GH 1 S		6308048	S
	4.1	TiBase B C 4.1 GH 1 L		6308097	L
	5.0	TiBase B C 5.0 GH 1 L		6308121	L
Tapered Screw-Vent	3.5	TiBase Z TSV 3.5 GH 1 L		6282581	L
	4.5	TiBase Z TSV 4.5 GH 1 L		6282599	L
	5.7	TiBase Z TSV 5.7 GH 1 L		6282607	L
M.I.S					
C1 Conical Connection	NP	CN-TB001 C1 NP GH 0.5		CN-TB001	L
		CN-TB015 C1 NP GH 1.5		CN-TB015	L
		CN-TB030 C1 NP GH 3*		CN-TB030	L
V3 Conical Connection	NP	VN-TB001 V3 NP GH 0.5		VN-TB001	L
		VN-TB015 V3 NP GH 1.5		VN-TB015	L
		VN-TB030 V3 NP GH 3*		VN-TB030	L
V3 Conical Connection, C1 Conical Connection	SP	CS-TB001 SP GH 0.5		CS-TB001	L
		CS-TB015 SP GH 1.5		CS-TB015	L
		CS-TB030 SP GH 3		CS-TB030	L
C1 Conical Connection	WP	CW-TB001 C1 WP GH 0.5		CW-TB001	L
		CW-TB015 C1 WP GH 1.5		CW-TB015	L
		CW-TB030 C1 WP GH 3		CW-TB030	L
SEVEN internal hex, M4 internal hex	NP	MN-TB001 INT HEX NP GH 0.5		MN-TB001	L
		MN-TBC15 INT HEX NP GH 1.5		MN-TBC15	L
		MN-TBC30 INT HEX NP GH 3*		MN-TBC30	L
	SP	MD-TB001 INT HEX SP GH 0.5		MD-TB001	L
		MD-TBC15 INT HEX SP GH 1.5		MD-TBC15	L
		MD-TBC30 INT HEX SP GH 3		MD-TBC30	L
	WP	MW-TB001 INT HEX WP GH 0.5		MW-TB001	L
		MW-TBC15 INT HEX WP GH 1.5		MW-TBC15	L
		MW-TBC30 INT HEX WP GH 3		MW-TBC30	L
Straumann					
Tissue Level	NNC (3.5 mm)	NNC Variobase C 3.5 GH 1		220.018	S
	RN (4.8 mm)	RN Variobase C GH 1		220.019	L
	WN (6.5 mm)	WN Variobase C GH 1		220.020	L
Bone Level	NC (3.3 mm)	NC Variobase C 3.3 GH 1		220.043	S
	RC (4.1 mm/ 4.8 mm)	RC Variobase C 4.1 GH 1		220.044	L
TLX	NT	NT Variobase C, TAN GH 0*		370.205	S
	RT	RT Variobase C, TAN GH 0*		371.205	L
	WT	WT Variobase C, TAN GH 0*		372.205	L
BLX	RB/WB	RB/WB Variobase C 3.8 GH 1.5*		624.981	S
	RB/WB	RB/WB Variobase C 4.5 GH 1.5*		624.982	L
	WB	WB Variobase C 5.5 GH 1.5*		624.983	L
Camlog					
Camlog	3.3	CAMLOG® Titanium base CAD/CAM, for Ø 3.3 mm GH 0.4*		K2244.3348	S
	3.8	CAMLOG® Titanium base CAD/CAM ,for Ø 3.8 mm GH 0.3		K2244.3848	S
		CAMLOG® Titanium base CAD/CAM PS, for Ø 3.8 mm, GH 0.8*		K2210.3808	S
	4.3	CAMLOG® Titanium base CAD/CAM, for Ø 4.3 mm GH 0.3		K2244.4348	S
		CAMLOG® Titanium base CAD/CAM PS, for Ø 4.3 mm, GH 0.8*		K2210.4308	S
	5.0	CAMLOG® Titanium base CAD/CAM, for Ø 5.0 mm GH 0.3		K2244.5048	L
		CAMLOG® Titanium base CAD/CAM PS, for Ø 5.0 mm, GH 0.8*		K2210.5008	L
	6.0	CAMLOG® Titanium base CAD/CAM, for Ø 6.0 mm GH 0.3		K2244.6048	L

Implant System		Titanium Base	REF	Size
Manufacturer/line	Platform	Name		
ConeLog	3.3	CONELOG® Titanium base CAD/CAM, for Ø 3.3 mm, GH 0.8 mm*	C2242.3308	S
		CONELOG® Titanium base CAD/CAM, for Ø 3.3 mm, GH 2.0 mm*	C2242.3320	S
	3.8	CONELOG® Titanium base CAD/CAM, for Ø 3.8 mm, GH 0.8 mm*	C2242.3808	S
		CONELOG® Titanium base CAD/CAM, for Ø 3.8 mm, GH 2.0 mm*	C2242.3820	S
	4.3	CONELOG® Titanium base CAD/CAM, for Ø 4.3 mm, GH 0.8 mm*	C2242.4308	S
		CONELOG® Titanium base CAD/CAM, for Ø 4.3 mm, GH 2.0 mm*	C2242.4320	S
iSy	5.0	CONELOG® Titanium base CAD/CAM, for Ø 5.0 mm, GH 0.8 mm*	C2242.5008	L
		CONELOG® Titanium base CAD/CAM, for Ø 5.0 mm, GH 2.0 mm*	C2242.5020	L
	4.5	iSy® Titanium base CAD/CAM, Ø 4.5 mm, GH 0.8 mm*	P2244.4408	S
		iSy® Titanium base CAD/CAM, Ø 4.5 mm, GH 2.0 mm*	P2244.4420	S
	5.2	iSy® Titanium base CAD/CAM, Ø 5.2 mm, GH 0.8 mm*	P2244.5008	L
		iSy® Titanium base CAD/CAM, Ø 5.2 mm, GH 2.0 mm*	P2244.5020	L

* Not valid for USA.

** Not all combinations of compatible implant systems and ti-bases might be available in all countries.

1.3 Intended user

Dental professionals (dental technicians and dentists)

1.4 Intended patient population and medical conditions

CEREC Cercon 4D™ Abutment Blocks are intended for patients in need of long-term restorative or prosthodontic dental therapy, or esthetic dental corrections. The use of ceramics is not limited to a specific patient population.

1.5 Contraindications

- Insufficient space available
- For restorations with angulation correction of more than 20° to the implant axis
- For individual tooth restorations with free end saddle
- For restorations whose length exceeds a ratio of 1:1.25 in comparison to the length of the implant.
- Contraindicated for use with TiBases not listed within Compatible TiBases Table
- Temporary Cementation to TiBase
- Parafunction (bruxism)
- Cantilever bridges
- Patients with a substantially reduced residual dentition

1.6 Composition

CEREC Cercon 4D™ Abutment Blocks are a high strength* multidimensional zirconia material that is tooth colored.

*High strength claim is based on material properties and in comparison, with competitive materials based on compliance to biaxial strength related to ISO 6872 standard. The strength claim is not indicative of the clinical use of that of an abutment.

1.7 Compatible stains and glazes

The use of glaze or stains is optional. CEREC Cercon 4D™ Abutment Blocks are compatible with Dentsply Sirona's Universal Stains and Glaze, Dentsply Sirona Universal Spray Glaze, or Dentsply Sirona Universal Spray Glaze Fluo (available separately, see complete Instructions for Use).

1.8 Compatible luting cements

CEREC Cercon 4D™ Abutment Blocks restorations are compatible with Universal/Self-Adhesive and Adhesive Resin Cements. Use of other manufacturer's cements or cement systems with CEREC Cercon 4D™ Abutment Blocks are at the discretion and sole responsibility of the dental practitioner. CEREC Cercon 4D™ Abutment Blocks should be cemented with Calibra Abutment Resin Cement (K240888) for use for extraoral cementation to a TiBase.

1.9 Technical Data

The following specifications apply to CEREC Cercon 4D™ restorations that have been sintered in a CEREC SpeedFire®, an inFire® HTC speed, or an inLab® Profire sintering furnace. Coefficient of thermal expansion (20 – 500 °C): $10.35 \times 10^{-6} \text{ K}^{-1}$

Flexural strength (3-point flexural strength): > 1100 MPa

2. GENERAL SAFETY NOTES

Be aware of the following general safety notes and the special safety notes in other sections of these directions for use.

Safety alert symbol

This is the safety alert symbol. It is used to alert you to potential personal injury hazards. Obey all safety messages that follow this symbol to avoid possible injury.

MRI Safety Information

See Dentsply Sirona MRI safety information available at <https://www.dentsplysirona.com/ifu>.

2.1 Warnings

- If properly processed and used, adverse effects from this medical device are highly unlikely. However, reactions of the immune system (such as allergies) or localized paresthesia (such as an irritating taste or irritation of the oral mucosa) cannot be completely ruled out. In case of skin sensitization or rash, discontinue use and seek medical attention.
- CEREC Cercon 4D™ Abutment Block restorations are not suitable for patients with clinical symptoms of parafunctional habits or bruxism (see Contraindications).
- Do not inhale dust particles during grinding. Wear suitable protective mask.
- Do not use competitive stains and glazes (except as indicated above) as it may impact the performance of the material.
- Do not use any dental cement, ceramic or other top-half component (and/or hybrid abutment-crown component) materials, scanners, milling units and components, CAD/CAM software, templates, and tools other than those specifically identified in this labeling.
- Choosing the wrong glazing program may lead to damage of the restoration or CEREC SpeedFire®.
- In patients with hypersensitivity to any of the ingredients, this medical device may not be used at all or only under the advice of the dentist or physician in charge.
- Highly angled abutments on small diameter implants are not recommended for use in the posterior region.

2.2 Precautions

- This product is intended to be used only as specifically outlined in these Instructions for Use. Any use of this product inconsistent with the Instructions for Use is at the discretion and is the sole responsibility of the practitioner.
- Wear suitable protective eyewear, clothing and gloves. Protective eyewear is recommended for patients.
- Devices marked "single use" on the labeling are intended for single use only. Discard after use. Do not reuse in other patients in order to prevent cross-contamination.
- CEREC Cercon 4D™ Abutment Block restorations require adequate interocclusal and interproximal thickness. Insufficient wall thickness may lead to premature failure.
- CEREC Cercon 4D™ Abutment Block restorations must be sintered and must be polished or glazed before insertion. Do not allow glaze to enter the screw channel. If you want to do a partial stain and glazing on the vestibular areas then have to polish the unglazed areas after the glazing.
- Insufficient data exist to support use of CEREC Cercon 4D™ Hybrid Abutments in conjunction with multiple-unit superstructure splints, prostheses and restorations.
- Use only in well ventilated areas.
- Do not cement with provisional cements. Use with provisional/temporary cements could cause fracture or separation of restorations.
- Adjustment in the sintered state should be performed with fine diamonds in a highspeed handpiece with water cooling and with low pressure.
- As an alternative it is possible to do adjustments with soft diamond rubber polishers and a handpiece at low speed and low pressure.
- Use of Dentsply Sirona Universal Spray Glaze Fluo with CEREC Cercon 4D™ Abutment Block restorations may result in a brighter shade appearance. Use of Dentsply Sirona Universal Spray Glaze is recommended if this effect is not desired.
- Conventional sintering may lead to higher chroma.

2.3 Interactions

CEREC Cercon 4D™ Abutment Blocks are designed to be fabricated using a Dentsply Sirona CAD/CAM system. Milling Blocks using non-compatible CAD/CAM system may lead to inadequate or unacceptable restorations.

2.4 Adverse reactions

No adverse reactions have been reported for CEREC Cercon 4D™ Abutment Blocks. Should you hear or receive information about any adverse effects, notify Dentsply Sirona.

2.5 Storage conditions

Inadequate storage conditions may shorten the shelf life and may lead to malfunction of the product. Store in a dry place. Protect from moisture. Do not use after the expiration date.

3. STEP-BY-STEP INSTRUCTIONS

3.1 Design

1. Plug a TiBase onto the lab implant of the master model or a CEREC® ScanPost directly onto the implant in the mouth. Plug a scanbody onto this until it comes to rest on the shoulder of the implant without any gaps. The scanbody is scannable without powder or scan spray. See Manufacturer's complete Instructions for Use for TiBase and Scanbody.
2. Make sure to select the CEREC Cercon 4D™ Abutment Block material under Dentsply Sirona in the CEREC® Software prior to scanning as this will ensure the correct firing program will be transferred.
3. Acquire the situation alternatively with e.g., PrimeScan, inEos X5 or CEREC® AC, See Manufacturer's complete Instructions for Use for acquisition equipment and technique.

Design Information

- Maintain a minimum wall thickness of 0.5 mm circularly for abutment (meso) around the screw channel.
- Make sure that no sharp edges or corners are produced.
- Design the outer form of the mesostructure in adherence to the preparation guidelines for the required superstructure.
- If the mesostructure is to be veneered immediately, make sure that this doesn't narrow the screw channel. The connection point to the base and the screw channel should not be coated.

3.2 CAD/CAM processing

CEREC Cercon 4D™ Abutment Blocks for restorations are produced with CEREC® or inLab® CAD/CAM systems by Dentsply Sirona. If you have any questions about these systems, please contact Dentsply Sirona.

3.2.1 Software requirements

CEREC Cercon 4D™ Abutment Blocks are supported by CEREC® version 5.2.7.

3.2.2 Milling Unit and Block sizes

- In the CEREC® software, select the CEREC Cercon 4D™ Abutment Block material. The Dentsply Sirona CAD/CAM milling unit will prompt you to insert a CEREC Cercon 4D™ Abutment Block. For detailed process, please consult the Instructions for Use and technical manuals of the appropriate CAD/CAM systems. Make sure to follow the manufacturers' recommendations.

For CEREC Cercon 4D™ Abutment Blocks dry milling is recommended, but, if desired, wet milling is also possible. CEREC SpeedFire® furnace offers a program with a pre-drying function.

Wet Milling – To avoid contamination and reduced translucency

Using wet milling, you have to be aware that contaminated cooling water (e.g. due to residuals of glass-ceramic particles) can diminish the translucency of the final restoration.

Therefore, renewal of the cooling water and cleaning of the milling chamber and filter is strictly recommended before wet milling.

CAUTION Alternatively, three separate water tanks for glass ceramics, rinsing in-between, and zirconia wet milling could be used.

3.3 Sintering

3.3.1 Preparation for sintering

After the milling process (and prior to sintering) a tungsten carbide tool shall be used to separate the restoration from the blank. Remove the sprue prior sintering.

If further adjustments are needed, it is recommended to do the adjustments before sintering.

Airborne particles – To reduce inhalation risk

- Do not inhale abrasive dusts.
- Use a vacuum system and wear a mask.

3.3.2 Sintering the restoration

A Dentsply Sirona Sintering furnace (CEREC SpeedFire®, inFire HTC speed, or inLab® Profire) is recommended for the sintering process. Before sintering the restoration, it is recommended to free the restorations from dust using compressed air or a ceramic brush.

Do not inhale abrasive dusts. Use a vacuum system and wear a mask.

NOTE: Restorations made from CEREC Cercon 4D™ Abutment Blocks must be sintered in dry conditions. Therefore, wet milled restorations have pre-drying steps included in their sintering programs, prolonging overall sintering time.

Sintering in the CEREC SpeedFire®

When sintering CEREC Cercon 4D™ Abutment Block restorations in the CEREC SpeedFire®, the CEREC® software automatically transfers the job to the CEREC SpeedFire® if the milling machine and the CEREC SpeedFire® are connected to the CEREC® Software.

NOTE: Maximum restoration size

Observe the maximum furnace chamber size when loading the Furnace:

- Diameter: 38 mm
- Height: 20 mm

The restoration (including the Glazing Support Single/Multi Unit) must not exceed the size of the chamber, length (38 mm) and height (20 mm), otherwise the chamber may be damaged.

For further information please refer to the Operating Instructions of the CEREC SpeedFire®.

1. Place the restoration with the occlusal surface facing down directly on the top door insulation.
2. Start the process by touching the start icon. The furnace shuts automatically once the process starts.
3. The furnace opens automatically following successful heat treatment. The process is not yet complete, as a cooling-off phase occurs in an open condition. There will be a signal tone once the cooling process has completed. When the unit's LED status display is green, the furnace can be unloaded after approximately 30 seconds.

Risk of injury

The restoration and parts of the door insulation may still be hot when the unit's LED status display is green.

Always use tweezers to remove the restoration from the furnace. Allow the restoration to cool down for another five minutes before picking it up with hands.

4. Only unload the furnace using metallic or ceramic tweezers. Plastic tweezers are not suitable for removing restorations, as the restoration is still very hot in this phase. Place the restoration on the fireproof tray for further cooling.

NOTE: Sintering times of dry-milled restorations can be reduced by placing the restoration in a 400 °C preheated CEREC SpeedFire®. For this, use the preheating function of CEREC SpeedFire® before starting the sintering process.

inFire HTC speed:

Dry milled restorations

When sintering in the inFire HTC speed use the pre-programmed program "inCoris TZI/ZI speed".

	Heating rate/°C/min	Holding temperature/°C	Holding time/min
S4	99	750	0
S3	99	1510	0
S2	50	1510	30
S1	99	800	5

Wet milled restorations

When sintering in the inFire HTC speed use the pre-programmed program "inCoris TZI/ZI speed wet".

	Heating rate/°C/min	Holding temperature/°C	Holding time/min
S4	99	750	0
S3	50	1510	30
S2	99	800	5
S1	15	80	30

inLab® Profire:

Dry milled restorations

When sintering in the inLab® Profire use the pre-programmed program "CEREC® Zirconia+ speed".

	Heating rate/°C/min	Holding temperature/°C	Holding time/min
S1	99	800	5
S2	50	1510	30
S3	99	750	0

Wet milled restorations

When sintering in the inLab® Profire use the pre-programmed program "CEREC® Zirconia+ speed wet".

	Heating rate/°C/min	Holding temperature/°C	Holding time/min
S1	15	80	30
S2	99	800	5
S3	50	1510	30
S4	99	750	0

Sintering CEREC Cercon 4D™ with other furnaces

For sintering the restorations in furnaces from other manufacturers, follow the firing parameters outlined in Tables below or visit www.dentsplysirona.com.

Dry milled restorations

Start	Heating rate 1	Temp. 1	Holding time 1	Heating rate 2	Temp. 2	Holding time 2	Cooling rate	Temp. 3	Action
°C	°C/min	°C	min	°C/min	°C	min	°C/min	°C	Open door
RT	99	800	5	50	1510	30	99	750	

Wet milled restorations

Start	Heating rate 1	Temp. 1	Holding time 1	Heating rate 2	Temp. 2	Holding time 2	Heating rate 3	Temp. 3	Holding time 3	Cooling rate	Temp. 4	Action
°C	°C/min	°C	min	°C/min	°C	min	°C/min	°C	min	°C/min	°C	Open door
RT	15	80	30	99	800	5	50	1510	30	99	750	

3.3.3 Adjustments of sintered restorations

Surface Condition – To reduce risk of compromised flexural strength

The surface condition of ceramic materials is critical for their flexural strength. Adjustment of the sintered restorations with milling tools must be avoided.

However, if adjustment is necessary, then follow these basic rules:

- Adjustment in the sintered state should be performed with fine diamonds in a highspeed handpiece with water cooling and with low pressure. Adjustments made with diamonds instruments must be followed by polishing.
- As an alternative it is possible to do adjustments with soft diamond rubber polishers and a handpiece at low speed and low pressure. The tool must be applied flat to minimize the chatter.
- Do not adjust or modify the shoulder area to be cemented to the TiBase.

3.4 Polishing and optional Staining and Glazing

CEREC Cercon 4D™ Abutment Block restorations can either be polished or glazed.

3.4.1 Polishing

- CEREC Cercon 4D™ Abutment Block restorations can be polished with standard polishing agents for zirconia ceramics.
- Subsequent heat treatment (depressurization fire) is not necessary/recommended.
- The following Dentsply Sirona Kits for Zirconia are recommended for polishing the restoration.
 - DSZ21 – Dentsply Sirona Zirconia Intraoral Polishing Kit
 - DSU21 – Dentsply Sirona Universal Extraoral Polishing & Shaping Kit

*Note: To be used with little contact pressure!

3.4.2 Optional Staining and Glazing

- Glazing of the abutment crown is optional. If desired, glaze for the meso structure may be applied only on the area of the emergence profile. In both cases prevent infiltration of the glaze into the screw channel and the connection area (bottom of the crown) to the TiBase.

Avoid glaze application to the preparation form in hybrid abutment (meso) constructions.

- Do not apply or allow stain or glaze into the screw channel.
- Before applying the stains and/or glaze, the restoration must be clean and free of oil and other surface contaminants. Clean the restoration with a steam cleaner or in the ultrasonic cleaner with distilled water or a suitable cleaning fluid for 5 minutes. Any contamination after cleaning must be prevented.

Surface Condition – To reduce risk of excessive wear

- Prior to glazing, ensure restorations are properly sintered.
- Restorations must be clean and dry before application.
- Optional staining and/or glazing must be completed before cementation to the TiBase.

Preparation of the restoration for Staining and/or Glazing

Use the Glazing Support Single/Multi Unit holders and CEREC® SpeedPaste. Add a small portion of CEREC® SpeedPaste to the restoration to ensure a safe hold on the pin. Do not fill the entire restoration with the firing paste. Place the Glazing Support Single/ Multi Unit holders in the paste and pick up the tweezers, if required. Observe the operating instructions supplied with the CEREC® SpeedPaste.

Glazing with Dentsply Sirona Universal Spray Glaze or Dentsply Sirona Universal Spray Glaze Fluo

1. Pick up the support with the restoration using tweezers.
2. Shake the spray can vigorously immediately before applying.
3. Maintain a distance of 6 – 10 cm (2.5 – 4.0 in) between the nozzle outlet and the restoration surface.
4. Spray on an even coat of Dentsply Sirona Universal Spray Glaze or Dentsply Sirona Universal Spray Glaze Fluo from all sides directly onto the (stained) restoration. Shake the spray can vigorously between individual sprays.
5. Apply a uniform layer glaze to the surface of the restoration. The applied glaze should be thin and uniform.
6. Ensure the intaglio surface of the restoration is free of glaze. **NOTE: Prevent any glaze from the mating shoulder or from entering the screw channel. For mesostructure parts (custom abutments), do not glaze the outer area which will work as cementing area.** If the glazing material is present on these surfaces of the restoration, remove it with a dry, stiff, short-bristled brush.
7. Keep the can as upright as possible during the spraying process.
8. Wait for a few seconds until the glaze dried and displays an even, whitish coating. Apply more spray as required.
9. If another application of spray is necessary after the firing, the spraying and firing can be performed once more in the same way.

Refer to Instructions for Use (IFU) of Dentsply Sirona Universal Spray Glaze/Dentsply Sirona Universal Spray Glaze Fluo for further details.

NOTE:

- If characterization is desired, the stains must be applied prior to application of Dentsply Sirona Universal Spray Glaze.
- The Use of Dentsply Sirona Universal Spray Glaze Fluo may lead to a brighter shade. Use Dentsply Sirona Universal Spray Glaze if this effect is not desired (see complete Instructions for Use).

Glazing with Dentsply Sirona Universal Stain and Glaze

1. Withdraw desired amount of Dentsply Sirona Universal Glaze from the jar and place it on the mixing palette.

- NOTE:** If the stain or glaze inside the jar has separated, mix thoroughly with a glass or zirconia spatula.
- If a thinner consistency is desired, dilute the material with the Dentsply Sirona Stain and Glaze Liquid (REF 601315/15 ml, 601350/50 ml).
 - Apply a thin layer of the glaze material all over the crown surface in the usual manner using a brush. Make sure not to apply the glaze material too thick or too thin. Too thin of a layer will result in a more matte finish. Too thick layers lead to bubbles on the surface.
 - If a more intensive shade effect is desired, Dentsply Sirona Universal Stains can be applied on the glazed surface of the restoration. Apply a thin layer of the stain material on the crown surface in the usual manner using a brush. Run a stain and glaze cycle in the CEREC SpeedFire® or a ceramic furnace to fix the applied Stains. After that, apply and fire the glaze material. Make sure not to apply the glaze material too thick or too thin.

Refer to Instructions for Use (IFU) of Dentsply Sirona Universal Stain and Glaze for further details.

3.4.3 Glazing the Restoration in the furnace

NOTE: A restoration can be glazed a max of two times in the furnace.

Glazing CEREC Cercon 4D™ Abutment Block restorations in the CEREC SpeedFire®

Glaze restorations individually and do not glaze multiple restorations at the same time. Position the Support with the restoration centrally on the top door insulation and ensure that the Glazing Support Single/Multi Unit or restoration do not protrude out from the door insulation; otherwise, these may collide with the furnace chamber.

Two different glazing programs are available. Select the "GLAZING (Spray)" program for use of Dentsply Sirona Universal Spray Glaze and Dentsply Sirona Universal Spray Glaze Fluo or "GLAZING (Stain & Glaze Paste)" program for use of Dentsply Sirona Universal Stain and Glaze on the CEREC SpeedFire® control panel and allow the unit to cycle.

NOTE: Using Dentsply Sirona Universal Stain and Glaze for crowns and bridges preheating to 400°C is recommended.

Choose the right Glazing program

- Choosing the wrong glazing program may lead to damage of the restoration or CEREC SpeedFire®.

For further information please refer to the Operating Instructions of the CEREC SpeedFire®.

Glazing CEREC Cercon 4D™ Abutment Block restorations with other furnaces

For firing the restoration in furnaces from other manufacturers, follow the firing parameters outlined in Table below.

Drying	Closing	Pre-heating temperature	Pre-heating	Heating rate	Final Temperature	Vacuum	Holding time	Cooling
min	min	°C	min	°C/min	°C	min	min	min
0	2	400	0	55	760	0	2:00	0

Use of glazes and stains by third-party manufacturers

Glazes and stains by third-party manufacturers may require other firing temperature and cycles than those recommended for CEREC Cercon 4D™ Abutment Block. Use only Dentsply Sirona Universal Stain and Glaze System and compatible spray glazes.

3.5 Extraoral Cementation of CEREC Cercon 4D™ structure abutment/crown to TiBase abutment

3.5.1 Preparation of TiBase

- The diameter of the TiBase must not be reduced, modified or adjusted e.g., by grinding.
- Shortening the TiBase is not permitted.
- The underside contact/connection surfaces of the TiBase to the implant should not be sand blasted, polished or otherwise modified. Use of an appropriate analogue to protect connection surfaces is strongly recommended. Wax or impression material may also be placed on exposed emergence profile surfaces for protection.

- Attach the TiBase to the lab implant analogue or polishing tool using screw.
- Cover the hex head of the abutment screw and fill the screw channel with easily removable material, e.g., soft wax, impression material, cotton or plumber's tape.
- Ensure that the CEREC Cercon 4D™ structure can be fully inserted onto the TiBase. The shoulder of the CEREC Cercon 4D™ structure must contact the TiBase flange without perceptible gap. Place small alignment marks on the CEREC Cercon 4D™ structure and the analogue using a felt-tip marker to assist aligning the correct Bu-Li position.
- Sandblast only the surfaces of the TiBase intended for cementation (50 µm aluminum oxide, max. 2.0 bar). Do not sandblast the polished tissue surface of the flange.
- Clean the TiBase with alcohol, ultrasonic or steam. If necessary, reseal the abutment screw and screw channel. Dry with clean, oil-free air.
- Once cleaned, do not touch or contaminate the bonding surface. Calibra Abutment Resin Cement (K240888) should be placed directly on the sandblasted bonding surface according to manufacturer's complete Instructions for use.

3.5.2 Preparation of CEREC Cercon 4D™ structure

- Sandblast the bonding surfaces of the CEREC Cercon 4D™ Abutment Block restoration with 50 µm aluminum oxide at up to 2.0 bar.
- Clean the restoration with an ultrasonic cleaner, steam cleaner or with alcohol.

3.5.3 Cementation of CEREC Cercon 4D™ structure to TiBase

- Prepare and dispense material following the Calibra Abutment Resin Cement (K240888) manufacturer's Instructions for Use.
- Apply the cement to the bonding surfaces of the TiBase and the inside of the restoration. Minimize the cement allowed to enter the restoration screw channel.
- Align the pen marks and push the CEREC Cercon 4D™ structure over the TiBase bonding chimney up to the end stop. Note the locking of the rotation and position lock.
- Immediately seat the restoration with gradual pressure. Verify complete seating. A gentle rocking or vibratory motion may be helpful to insure optimal seating.
- Ensure margin between TiBase and CEREC Cercon 4D™ structure is closed and sealed.

To reduce risk of Incomplete cementation

- Assure restoration is fully seated within cement working time.
- Placement of multiple units may require multiple mixes of cement.

- Protect restoration from contamination and movement until the final set of the cement according to the manufacturer's Instructions for use.

- Debonding – To reduce risk of Inadequate retention due to insufficient stabilization**
- Stabilize restoration during cleanup.
 - Protect restoration from movement during setting.
 - Allow cement to cure without disturbing according to the manufacturer's Instructions for use.

3.5.4 Cleaning Marginal Excess

Self-cure cleanup

Follow the manufacturer's Instructions for Use.

3.5.5 Finishing and Polishing

1. Rinse air barrier material under running water.
 2. Remove wax or other sealing material and any residual cement from the screw channel.
 3. Remove any cement excess at the margin, finish and polish margin using rubber rotary tools.
- Proceed to Hygiene for pre-use appliance preparation.

Residual excess cement – To reduce the risk of soft tissue irritation

- Ensure all excess cement is removed.
- Ensure all surfaces are polished and free of contamination.

3.6 Optional Try-In

- Try-in the restoration for marginal and proximal fit. Make any necessary adjustments as outlined above.

Contamination – To reduce risk of infection

- Restorations should be polished, cleaned, disinfected and sterilized before and after optional try-in. See Hygiene section below.
- Adjustments made during try-in must be re-polished as outlined above.
- Do not glaze restorations cemented to TiBase. Cement bond will be destroyed.

4. HYGIENE, PROCESSING AND DISPOSAL

Cross-contamination – To reduce risk of Infection

- Do not reuse single use products.
- Dispose of in accordance with local regulations.

4.1. Assembled Restoration

Cross-contamination – To reduce risk of Infection

- Restorations must be polished or glazed before cleaning, disinfection and/or sterilization.
- Restorations must be processed according to the validated processing Instructions below prior to delivery.

Ensure all components are able to tolerate processing, e.g., immersion in boiling water for 5 minutes, Automated Washer-Disinfector up to A0≥3.000, and steam autoclaved.

- Dynamic-Air-Removal (Pre-Vacuum) Steam autoclaved at:
 - 132 °C (270 °F) for 4 minutes; Dry Time = 20 minutes
 - 134 °C (273 °F) for 3 minutes; Dry Time = 16 minutes (recommended only for non-US markets)
 - 135 °C (275 °F) for 3 minutes; Dry Time = 16 minutes

without significant effect on physical properties.

- Locally applicable legal regulations and the hygiene standards applicable for a dental practice must be observed.
- Use only the validated processing procedures specified.
- The responsibility for the microbial processing of the Prosthesis Assembly lies with the user.
- The equipment and devices must be properly maintained and serviced at regular intervals.
- The fabricator (laboratory or in-office dental technician) of the assembled Restoration must inform the dentist of the need to process the device before using or inserting it in the patient's mouth. Below are the acceptable methods for cleaning, disinfecting and sterilization. For detailed parameters, see instructions for cleaning, disinfecting and sterilization below.

Process		✓ Acceptable	✗ Unacceptable
Cleaning	Mechanical (Ultrasonic)		✗
	Manual	✓	
	Automated (Washer/Disinfector)	✓	
Disinfection	Manual	✓	
	Automated (Washer/Disinfector)	✓	
	Immersion	✓	
Sterilization	Steam Autoclave	✓	

- The user will select the cleaning, disinfection, and sterilization process according to the equipment available.

Instructions for cleaning, disinfecting and sterilizing TiBase/CEREC Cercon 4D™ Hybrid Abutment and Hybrid Abutment Crown

Warnings	<ul style="list-style-type: none"> High-level disinfection has not been validated as a terminal process for the assembled Hybrid Abutment Restoration prior to delivery. Disinfection is not recommended as a terminal process. Proceed to Sterilization following cleaning and any optional/interim disinfection process. Steam autoclaving sterilization is appropriate and recommended for the assembled Hybrid Abutment Restoration prior to delivery. Cleaning/Disinfection of cemented assembly using an ultrasonic unit is not validated and should not be used. Follow proper infection prevention activities, such as proper hand washing and donning new puncture and chemical resistant gloves at pertinent steps. Do not use hard wire brushes to clean the assembled Hybrid Abutment Restoration as scratching and mechanical degradation can occur. Use only a disinfecting solution which is approved for its efficacy, EPA registered (and/or Health Canada approved) and use in accordance with the IFU of the disinfecting solution manufacturer. Do not allow the assembled Hybrid Abutment Restoration to exceed 137 °C. Always use a pH neutral instrument cleaning solution if authorized. Sterilize metal and plastic components in separate pouches to avoid damage to components. If the integrity of the pouch has been compromised, the device should be processed again before use. In the USA, the assembled Hybrid Abutment Restoration must be cleaned and sterilized before delivery. For the countries that require a three-step process, in which disinfection is required prior to sterilization, the assembled Hybrid Abutment Restoration must be cleaned, disinfected, and sterilized before delivery.
Limitations on Reprocessing	<ul style="list-style-type: none"> The assembled Hybrid Abutment Restoration is intended to be cleaned and sterilized only prior to try-in and/or delivery to the mouth. Repeated processing has not been validated. Cold liquid immersion sterilization, chemical vapor sterilization, and dry heat sterilization methods have not been tested or validated for efficacy and are not recommended for use.
Initial Treatment at the Point of Use	<ul style="list-style-type: none"> If applicable, remove all excess cement, finish and polish assembled restoration. Do not allow residue or any form of contamination to dry on the device. Remove excess soil/gross debris by steam cleaning, wiping with a cloth soaked with 70% Isopropanol, brushing, and/or rinsing under water. Transport to processing area and follow validated processing instructions below. It is recommended that the restoration be processed as soon as is reasonably practical following assembly.
Preparation before cleaning	Process the Hybrid Abutment Restoration fully assembled, adjusted and polished, ready for try-in and/or delivery.
Cleaning: Mechanical	No mechanical cleaning process has been validated.
Cleaning: Manual	<ul style="list-style-type: none"> Cleaning of the product surfaces using a clean, soft nylon brush under fluent tap water for 10 seconds. Pay special attention to crevices and lumens (e.g., of connection points) and to the screw-channel especially for TiBases with captured screws. Brush away all visible contaminations. Rinse the product under tap water. Dry the product completely using a lint free cloth
Disinfection: Manual Optional interim process for countries that require a 3-step process (recommended only for non-US markets)	<ul style="list-style-type: none"> Prepare boiling water (tap water has sufficient quality). Put the product into a heat-resistant container and douse the whole product with boiling water. Remove air bubbles by moving the product e.g. using a tweezer (make sure that you don't have any contact with the boiling water). Pay special attention to crevices and lumens (e.g. of connection points) and to the screw-channel, especially for TiBases with captured screws. Holding time of 5 minutes. Rinse the product thoroughly under tap water. Dry the product completely using a lint free cloth.
Cleaning and disinfection: Automated	<ul style="list-style-type: none"> A thermal disinfectant that complies with ISO 15883 may be used for automatic cleaning and disinfection. Use only properly maintained, calibrated, and approved washer disinfectors according to ISO 15883. Follow the manufacturer's instructions to select the solutions and programs. Remove excess soil with a disposable cloth/paper wipe. In case of gross contamination of the assembled Hybrid Abutment Crown Restoration, remove it by steam cleaning or using a soft brush and running water. Place the assembled Hybrid Abutment Restoration in the washer-disinfectant allowing water and detergent to enter into and drain out of device openings. Run washer-disinfection program with A0 value ≥ 3.000 (e.g., 5 min at ≥ 90 °C) (e.g., Miele Vario TD) using appropriate detergents (e.g., Neodisher MediClean and Neodisher Z, both from Dr. Weigert, Germany or equivalent). After the cycle is complete, visually inspect the assembled Hybrid Abutment Restoration for degradation due to processing. If the restoration or the cement margin is discolored, cracked, distorted, etc., discard and do not use. Proceed to Sterilization following Automated Cleaning and Disinfection.
Drying	The drying instructions are incorporated in the cleaning and disinfection section above.
Maintenance, Inspection and Testing	Visually inspect the assembled Hybrid Abutment Restoration for degradation due to processing. If the restoration or the cement margin is discolored, cracked, distorted, etc., discard and do not use.
Packaging	<ul style="list-style-type: none"> Ensure that the assembled Hybrid Abutment Restoration is completely dry before packaging for sterilization. A lint-free, disposable cloth may be used to dry the assembled Hybrid Abutment Restoration. Place each assembled Hybrid Abutment Restoration in a separate appropriately sized steam-sterilization pouch. Uses of FDA-cleared, ISO 11607 compliant Paper/plastic steam sterilization pouches (e.g., AssurePlus® Sterilization Pouches) are recommended.

Instructions for cleaning, disinfecting and sterilizing TiBase/CEREC Cercon 4D™ Hybrid Abutment and Hybrid Abutment Crown

Sterilization/Steam autoclaving	<ul style="list-style-type: none"> The assembled Hybrid Abutment Restoration must be sterilized utilizing steam autoclave sterilization. The following Dynamic-Air-Removal/Pre-Vacuum steam autoclave cycles may be used: <ul style="list-style-type: none"> • 132 °C (270 °F) for 4 minutes with a minimum drying time of 20 minutes • 134 °C (273 °F) for 3 minutes with a minimum drying time of 16 minutes (recommended only for non-US markets) • 135 °C (275 °F) for 3 minutes with a minimum drying time of 16 minutes <p>An alternate method of sterilization is to place non-bagged instruments into the (SFPP) steam sterilizer and run at one of the above listed cycles.</p> <p>Gravity Displacement steam autoclave cycles have not been validated.</p>
Storage	<ul style="list-style-type: none"> • Do not remove from the pouch until ready for patient delivery to prevent contamination. • Hybrid Abutment Restorations sterilized non-bagged should be used immediately. • Store at room temperature, away from moisture or excessive humidity. • To prevent contamination, store processed Hybrid Abutment Crown Restoration in covered storage such as a drawer or cabinet until delivery. • Prior to delivery, inspect the pouch. If the integrity of the pouch has been compromised, the Hybrid Abutment Restoration should be processed again before delivery.
Additional Information	<ul style="list-style-type: none"> • Prior to use, inspect the Hybrid Abutment Restoration. Discard any Hybrid Abutment Restoration which has become damaged, worn, or distorted. • The instructions provided above have been validated by the manufacturer of the medical device as being capable of preparing a medical device for reuse. It remains the responsibility of the processor to ensure that the processing, as actually performed using equipment, materials and trained personnel in the processing facility, achieves the desired result. This requires verification and/or validation and routine monitoring of the process. Likewise, any deviation by the processor from the instructions provided should be properly evaluated for effectiveness and potential adverse consequences. Use of other reprocessing procedures/methods is at the discretion and sole responsibility of the practitioner.
Manufacturer Contact	Within the United States, call Dentsply Sirona at 1-302-422-4511. For areas outside the United States, contact your local Dentsply Sirona representative.

4.2. Disposal

Dispose of in accordance with local regulations.

5. LOT NUMBER, EXPIRATION DATE AND CORRESPONDENCE

5.1 Do not use after expiration date. ISO standard uses: "YYYY-MM-DD"

5.2 The following numbers should be quoted in all correspondence:

- Reorder number
- Lot number
- Expiration date

5.3 Any serious incident in relation to the product should be reported to the manufacturer and the competent authority according to local regulations.

5.4 EU – A summary of the safety and clinical performance (SSCP) for this product can be found (upon activation) at <https://ec.europa.eu/tools/eudamed> by searching using the Basic UDI-DI number ++EDD1SMCERAMICSZ2 and at <https://www.dentsplysirona.com/ifu> using the reference number (REF).

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