# Patterson® Flowable Composite

## INSTRUCTIONS FOR USE

Patterson Flowable Composite is the ideal versatile radiopaque, flowable, light cured material. The composite is directly injected into the cavity preparation for maximizing adaptation to the preparation.

### COMPOSITION OF PATTERSON FLOWABLE COMPOSITE:

- 44% wt (68% vol) multifunctional methacrylic ester
- 56% wt (32% vol) inorganic filler (0.2 - 1 micron)

### COMPOSITION OF PATTERSON FLOWABLE COMPOSITE

#### LOW VISCOSITY:

- 48% wt (72% vol) multifunctional methacrylic ester
- 52% wt (28% vol) inorganic filler (0.2 - 1 micron)

#### CONTRAINDICATIONS:

- Pulp capping.
- Product may cause skin irritations to some people.
- In such cases, discontinue use and seek medical attention.
- Do not use in conjunction with any eugenol-containing materials.
- Patients with known sensitivity to acrylics or allergic reactions to the materials of this product should not use it.

#### INDICATIONS:

- Conservative Class I (including occlusal surfaces), II, III, IV restorations
- Class V restorations (cervical abrasion, root erosion, abfraction)
- Repair of defects in enamel and porcelain restorations
- Sealants (pit & fissure, implant)
- Radiopaque cavity liner under direct restorations
- Blocking out of undercuts
- Minor core build ups
- Repairing or sealing of temporary crowns
- Cementation of porcelain, ceramic, composite veneers
- Splinting (mobile teeth, fibre bridges)
- Cover stains.

#### STORAGE AND HANDLING:

- Store at temperatures between 2° and 25°C (35°- 77°F).
- Use at room temperature
- Syringe system: when not in use replace the cap tightly.

#### INSTRUCTIONS:

Clean and isolate tooth. A conservative cavity should be prepared, employing standard techniques and instruments, to form a slightly rounded internal form. Margins of the cavity preparation should end in sound and supported enamel with no bevels in stress bearing locations. If bevelling is desired in a low stress location it should extend no further than 1mm at an angle of no greater than 45°. Pre-wedging is also recommended to ensure that the restored tooth will have an adequate contact point. Prophylax all surfaces to be etched including surfaces adjacent to the cavity with an oil free non-fluoride containing paste or a slurry of pumice and water. Rinse thoroughly with water.

#### Shade Selection

1. Do not cover the tooth before shade matching. Select composite shade by matching with a clean and moist tooth using a shade guide.
2. When the adjacent teeth are significantly different in shade, try to match the lighter shade.
3. For visual confirmation of shade selected from the shade guide, place small amounts of composite on a surrounding tooth surface and light cure before bonding.

#### Pulp Protection

4. For deep cavities an appropriate liner or cement should be placed at the deepest point of the cavity.
5. Direct adhesive application: Follow the technique recommended by the adhesive you are using.

#### Self-Etch Technique

Bottle delivery: Dispense one or two drops (depending on restoration size) of self-etch adhesive into a mixing well. Immediately replace the cap.

**NOTE:** Proceed to the next steps quickly as the liquid can evaporate and result in poor bonding.

- a) Apply adhesive into the cavity by scrubbing for 10 seconds with a brush.
- b) Leave on cavity surface for 10 seconds.
- c) Blow with water/oil-free air until you see no movement of the bond (minimum 5 seconds).

**NOTE:** A vacuum aspirator can be used to prevent adhesive liquid from scattering.

- d) Light cure for 10 seconds (460-480nm wavelength, high power LED curing light).

**NOTE:** Enamel subjected to fluoridation should be etched for 90 to 120 seconds.

**IMPORTANT:** If using an etch syringe, ensure that the dispensing tip hub is firmly attached to the syringe by twisting the hub securely onto the syringe. On first usage or after prolonged storage, extrude a small amount onto a mixing pad for familiarity with the etchant's viscosity and rate of extrusion.

#### Selective Etch Technique

- a) Apply phosphoric acid etching gel to tooth enamel only and leave for 10 seconds.

- b) Rinse thoroughly with water and dry with water/oilfree air or with cotton pellets. Do not over dry.
- c) Apply adhesive according to steps in "Self-Etch Technique".

#### Total Etch Technique

- a) Apply phosphoric acid etching gel 37% to tooth structure (enamel & dentine) and leave for 10 seconds.
  - b) Rinse thoroughly with water and dry with water/oilfree air or with cotton pellets. Do not over dry.
  - c) Apply adhesive according to "Self Etch Technique".
6. Directly inject the flowable composite in increments of 2mm or less in:
    - Class V restorations,
    - Pit and fissure sealants,
    - Conservative Class I, II, III and IV restorations
    - Other indications as required.
- CAUTION:** ensure that the dispensing tip is **firmly** attached to the syringe by twisting the tip hub securely onto the syringe.
- WARNING:** Apply restorative at normal room temperature (23°C/74°F).
7. Light cure the composite using a suitable high power LED curing light (460-480nm wavelength) for a minimum of 20 seconds in increments of 2mm.
  8. Polishing of composite
    - (a) Remove excess composite and contour desired shape using a fine diamond or a 12-fluted carbide bur.
 

Tips for polishing using a flexible disc system:

      - (i) The polishing motion **should be constant and unidirectional.**
      - (ii) A back and forth movement over the composite - enamel margin is not recommended.
      - (iii) Keep the surface and polishing disc dry while polishing. A dry surface will produce a smoother, more uniform finish.
      - (iv) Do not use a handpiece with a speed greater than 35,000 rpm.
      - (v) Avoid touching the composite with the mandrel or disc eyelet.
    - (b) For gross reduction, set the speed of the handpiece to approximately 10,000 rpm. Use the coarsest grit disc of the range and polish the restoration starting at the gingival and **moving outward** over the restoration. Replace the disc if cutting efficiency impaired. Wash and dry.
    - (c) For final contouring, use a medium coarse grit at a speed of approximately 10,000 rpm. Wash and dry.
    - (d) For finishing, use a medium fine grit disc at a speed of approximately 30,000 rpm. Wash and dry.
    - (e) Finally, use the finest grit disc at 30,000 rpm. Wash and dry.
 

**NOTE:** After step 8, Polishing Paste may be used to give the composite a lustrous finish, using the following steps:

      - (i) Apply polishing paste to a rubber cup or disc.
      - (ii) Smear a thin layer of paste onto the restoration.
      - (iii) Add a small amount of water to both the tooth and cup. Polish for 30 seconds at low speed and light pressure.
      - (iv) Wash and dry.

#### CEMENTATION OF PORCELAIN, CERAMIC AND COMPOSITE VENEERS

Prior to cementation with Patterson Flowable Composite/Patterson Flowable Composite Low Viscosity, the restoration must be pre-treated according to the manufacturer's instructions for use.

Patterson Flowable Composite/Patterson Flowable Composite Low Viscosity is only suitable for the cementation of restorations that have translucency and are of appropriate thickness to allow polymerization light to pass through, to ensure that the light reaches Patterson Flowable Composite Low Viscosity, so that it can cure completely.

#### PRECAUTIONS:


- Avoid prolonged contact of etchant, adhesive and composite with the skin or oral tissue, as it may cause inflammation of the oral tissues or skin sensitization.
- For professional use only.
- Keep out of reach of children.
- Do not take internally.
- Do not use after expiry date.
- Do not expose material to direct light.
- Caution: Federal Law restricts this device to sale by or on the order of a dentist.

#### FIRST AID

- **Eye (contact):** Wash thoroughly with water and seek medical attention.
- **Skin (contact):** Remove using a cloth or sponge soaked in alcohol. Wash thoroughly with water.
- **Ingestion:** Rinse mouth thoroughly. Drink plenty of water/milk. Seek medical attention if symptoms persist.
- **Inhalation:** No symptoms expected.

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Patterson Dental Supply, Inc.  
1031 Mendota Heights Road  
Saint Paul, MN 55120

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