


# TempoCem ID

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English

### Product description

TempoCem ID is a translucent, dual-cure composite luting cement that is used for the temporary luting of crowns, bridges, veneers, inlays, and onlays, as well as all kinds of temporaries. The material is available in the Smartmix Syringe.

### Indications for use

Temporary and long-term temporary fixing of provisional crowns, bridges, inlays, onlays and veneers and implant-borne temporary prostheses.

### Limitations of use

Do not apply the material directly to the exposed dental pulp.

### Contraindication

Do not use the material where there are known allergies to any of the components or contact allergies.

### Patient target group

Persons receiving treatment as part of a dental procedure, for example:  
 Cementation of indirect restorations and prostheses for the treatment of large dental defects, tooth loss, malocclusions and malformations of the teeth.

### Intended user

Dentist

### Notes for use

- Light-curing units should emit at 450 nm and should be checked regularly. The light intensity should be at least 600 mW/cm<sup>2</sup>. Place the curing light as close as possible to the material. Observe the manufacturer's instructions for the light-curing unit used.
- The material adheres to composites. Apply a thin layer of a separating agent (e.g. vaseline, K-Y Jelly) to existing composite-based fillings, bases or build-up fillings before using the material. In the case of further use of the temporary, apply a separating agent in advance.
- Ultrasound is helpful for the removal of the cement from the restoration.

### Using the Smartmix Syringe (for image, see »Handling«)

1. Before attaching the mixing tip, remove the cap [A] or used mixing tip after turning it 90° counterclockwise, and discard it.
2. Attach a new mixing tip [B].

**Note:** In order to achieve an optimum mixing result, DMG recommends the use of the mixing tips available from DMG. All mixing tips are also available in refill packs.

**Note:** Make sure that the notches on the mixing tip line up with those on the Smartmix Syringe [C].

3. Twist the mixing tip clockwise by 90° to lock it.
4. The material is mixed in the mixing tip when it is expressed and it can be applied directly.

**Note:** When using a Smartmix Syringe for the first time, squeeze out a pea-sized amount of material and discard it.

### Recommended use

1. If required, apply a separating agent to the tooth core, core build-up, or temporary (see also »Notes for use«).
2. Apply the material to the temporary restoration.
3. Insert the restoration within 60 seconds after mixing has begun and attach to the dry, prepared tooth using slight pressure.
4. If necessary, instruct the patient to adopt the habitual occlusion.
5. Remove cement residues within 75 to 120 seconds after inserting the restoration with the aid of a microbrush, brush, foam pellet, or scaler. Excess cement is easier to identify after drying with oil-free and water-free air. Remove interdental excess carefully using dental floss. As an alternative, it is possible to remove excess cement by briefly curing it using light (1 to 2 seconds per side) and then detaching it with a scaler.
6. Allow the material to chemically cure completely for 420 seconds after mixing has begun.
7. Translucent restorations: Additionally cure the material with a suitable light-curing unit for at least 20 seconds per side with an exposure intensity of ≥ 600 mW/cm<sup>2</sup>, as required.

### Timing

Working time (including mixing time)*	≤ 60 s
Removal of excess (without light curing)	75 to 120 seconds after the restoration is inserted
Light-curing time for removal of excess (optional)	1–2 seconds per side
Intraoral curing	≥ 420 s
Light exposure time for curing (optional)	Min. 20 seconds per side with an exposure intensity of ≥ 600 mW/cm <sup>2</sup>

\* The specified times are valid at a room temperature of 23 °C (73 °F) and a normal relative humidity of 50 %. Higher temperatures shorten, lower temperatures lengthen these times.

## Interactions

- Residues of methacrylate-based materials may alter the setting behavior of silicone impression and bite registration materials.
- Materials containing eugenol, moisture and oily air may inhibit polymerization at the contact area.

## Residual risks/side effects

There are no known side effects to date. The residual risk of hypersensitivity to components of the material cannot be ruled out.

## Warnings/precautions

- For dental use only!
- Keep out of the reach of children!
- Avoid contact with the skin! In the event of accidental skin contact, immediately wash the affected area thoroughly with soap and water.
- Avoid contact with eyes! In the event of accidental contact with the eyes, immediately rinse thoroughly with plenty of water and consult a physician if necessary.
- Follow the manufacturer's instructions for other products that are used with the material/materials.
- Serious incidents involving this product must be reported to the manufacturer and to the responsible registration authorities.

## Product performance characteristics

Compressive strength  $\leq 150$  MPa

## Storage/disposal

- Store in a dry place at room temperature 15-25 °C (59-77 °F).
- Protect against exposure to direct sunlight!
- After application, leave the used mixing tip on the Smartmix Syringe as a seal!
- Do not use after the expiration date!
- Disposal must comply with national regulations.

## Composition

Dental glass, EBPADMA, aliphatic urethane methacrylate, unsaturated polyester resin blend, TEDMA, HEMA, SiO<sub>2</sub>, additives.

**Note:** A glossary of all symbols is displayed at: <http://www.dmg-america.com/literature>

**Caution:** US Federal Law restricts this device to sale by, or on the order of a dentist, or other practitioner licensed by law of the state in which he or she practices to use or order the use of this device.

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