High-Speed Handpiece

INSTRUCTIONS FOR USE

043-3904 | 043-3920 | 043-3862 | 043-3938 | 043-3854 | 043-3912

Thank you for selecting the Patterson[®] High-Speed Handpiece. Patterson[®] handpieces are engineered and manufactured to provide exceptional performance, reliability and precision. The advanced turbine design and quality components ensure quiet operation, high torque and years of dependable service. The handpiece is specifically designed to tolerate the rigors of repeated autoclaving. Please take a few minutes to thoroughly read the Instructions for Use to obtain peak performance and longevity from your new handpiece.

GENERAL INFORMATION

Important: *The Food and Drug Administration and The Centers for Disease Control and Prevention strongly recommend that handpieces be sterilized after each patient use.*

- Manually remove all visible soil, such as saliva, from the handpiece prior to the autoclave cycle. See subsequent infection control/reprocessing steps for details.
- Lubricate and maintain your handpiece in accordance with the instructions. Proper care and maintenance will extend the life of your new handpiece.
- The fine wire provided is used to clean the air/water port(s) at the head of the handpiece.

Intended Use: *Used for preparation of dental cavities for restorations such as fillings. High Speed handpieces are reusable, sterilizable medical devices used for the preparation of dental cavities and restorations by trained dental professionals.*

Caution:

- Wear a mask, gloves and protective eyewear when handling contaminated instruments.
- The device provided to you is not sterile. Sterilize the device before your first use.

BUR INSTALLATION/REMOVAL

The handpiece should only be operated with a bur securely in place. **Caution:** Do not use burs that are inappropriate (too long or short) for the handpiece. Burs must be able to be seated and locked in the chuck. **Patterson does not recommend the use of oral surgery burs at any time in our handpieces.**


(Figure 1)

(Figure 2)

(Figure 3)

(Figure 4)

(Figure 5)

For Jacob Chuck model:

1. Place bur changer wrench (084-2401) over handpiece head while pulling back the knob (Figure 1).
2. Lightly push the knob down while turning carefully so the square tip of the tool engages the square tip of the chuck (Figure 2).
3. To loosen the bur, rotate the wrench handle 1/4 to 1/2 turn counter-clockwise (Figure 3). Reverse to tighten.


Important: Do not overtighten. Do not tighten the chuck without a bur or bur blank in place. This will damage the chuck.

For Autochuck model:

1. Press the button at the head of the handpiece firmly while pulling on the bur to remove from the chuck (Figure 4).
2. To insert a bur, push firmly on the button and insert the bur into the chuck until firmly seated (Figure 5).
3. Release pressure on the button to lock the bur in place. Pull on the bur to ensure that it has been fully seated and locked.

Important: Do not press the button while the handpiece is in operation. This will release the chucking mechanism, the bur may be released prematurely and excess heat may be generated. Do not use the handpiece to retract the cheek or other soft tissue.

PORT CONFIGURATIONS

HANDPIECE CONNECTION AND USE

Caution: *The handpiece is supplied non-sterile. Sterilize the handpiece before first use and after each subsequent use.*

1. Apply a disposable sterile barrier sleeve over the handpiece.
2. Following instructions in the “BUR INSTALLATION/REMOVAL” section, install the desired attachment to your handpiece.
3. Prior to connecting the handpiece to the connector line:
 - a. Check that all tubing supplying the handpiece hose is properly connected with the dental unit
 - b. Make sure there are no restrictions in any of the lines
4. Connect the handpiece to the connector line.
 - a. Align the handpiece and port holes on the connector line properly
 - b. Nest the handpiece and connector line and secure
5. Operate the handpiece at 30-35 PSI (measured at the connector). Increasing air pressure will not increase performance and will prove detrimental to the turbine.
6. After use:
 - a. Allow the handpiece to fully stop
 - b. Disconnect the handpiece from the connector line
 - c. Remove and discard the disposable sterile barrier and bur in medical waste
 - d. Reprocess the handpiece as outlined in the “INFECTION CONTROL/STERILIZATION” section

INFECTION CONTROLS/STERILIZATION

Instructions for Steam Sterilizing the Patterson® High-Speed Handpiece.

Warnings:

- These instructions are for use ONLY on the Patterson® High-Speed Handpiece. Any other reusable components or systems should be reprocessed according to individual device's instructions for use.
- Wear a mask, gloves and protective eyewear when handling contaminated instruments.
- The device provided to you is not sterile. Sterilize the device before your first use.
- It is important your handpiece be properly sterilized and lubricated at all times. The life of "O" rings and rotating parts will greatly decrease if your handpiece is not maintained properly.

Limitations on Preprocessing:

- Do not use an automatic washer/disinfector for the Patterson® High-Speed Handpiece. Discoloration of the material may occur.
- Repeated cleaning and steam sterilization cycles have minimum effect on these instruments. End of life is normally determined by wear and damage due to use.
- Do not use chemical disinfectants, or solutions containing phenols, acids or chlorines prior to sterilization or rapid deterioration of the material may occur.
- Never sterilize the Patterson® High-Speed Handpiece in cold liquid disinfection/sterilization, chemical vapor sterilization or dry heat sterilization. These methods have not been tested or validated for efficacy and are not recommended for use.
- Do not immerse the device in an ultrasonic bath.
- Lubricate motor and/or attachments with a dental handpiece lubricant which will withstand autoclave temperatures/pressures.
- Never sterilize with a bur in the handpiece.
- Do not let the handpiece remain in a pouch under wet or damp conditions. Such conditions can result in corrosion of internal parts and shorten the life of the handpiece.

INFECTION CONTROLS/STERILIZATION (CONTINUED)

Point of Use:

- It is recommended that devices are reprocessed as soon as is reasonably practical following use.
- Remove the Patterson® High-Speed Handpiece from the air hose.
- Remove the disposable sterile barrier from the handpiece.
- Manually remove all visible soil, such as saliva, from the handpiece with a disposable cloth or paper wipe.

Containment and Transportation:

- Protect the device from contact with other dental instruments that may cause damage.

Preparation for Decontamination:

- Only process the Patterson® High-Speed Handpiece according to these guidelines. Any other components (e.g. air hose) must be reprocessed according to its product guidelines.
- Spray or drop lubricant/cleaner into the drive air hole at the back of the handpiece until lubricant runs from the handpiece head.
- Connect the handpiece to an air supply. Insert a bur/blank securely into the handpiece.
- Run the handpiece for 30-40 seconds. Dirt and debris will be expelled.
- Wipe exterior surface of handpiece with a dry 2" x 2" gauze or disposable towel.
- Remove rotating instrument from the handpiece.

Cleaning Manual:

- Generously apply running tap water to wet a clean cloth.
- Thoroughly wipe the handpiece with the wetted cloth:
 - Using a soft bristle brush, brush the handpiece until all visible soil is removed, paying attention to areas that may come in more contact with the patient or appear to be stained.
- Rinse the handpiece under portable running tap water for a minimum of 30 seconds to remove solution residue.

INFECTION CONTROLS/STERILIZATION (CONTINUED)

- Dry the handpiece using a clean cloth.
- Visually inspect the handpiece for visible soil. If soil is seen, repeat the cleaning procedure.

Disinfection:

- Do not attempt to disinfect the device. Sterilize only.

Packaging:

- Place each Patterson® High-Speed Handpiece in a separate paper/ plastic steam-sterilizing pouch.

Steam Sterilization:

- Place bagged devices into an ANSI AAMI ST55 or ST8 compliant steam autoclave according to the manufacturer's instructions. The following settings should be used:
Gravity Steam Sterilization
 - Full Cycle: 135 °C (275 °F) for 10 minutes.

Drying:

- To dry, use the drying cycle of the steam autoclave. Set cycle for 60 minutes. Do not exceed 137°C (278 °F).
- After removing the handpiece from the autoclave, allow the handpiece to return to room temperature before use.

Cleaning Manual:

- Visually inspect to ensure that all contamination has been removed.
- Check for distortion, damage or wear. Contact Patterson customer service if devices are damaged or worn.
- ISO 14457 stipulates a durability of at least 250 sterilization cycles. It is recommended that the handpiece be supplied to an authorized service representative for regular service after 600 sterilization cycles or one year of use.

Storage:

- To maintain sterility, instruments should remain bagged until ready for use.

Manufacturer Contact:

- In the United States, contact Patterson customer service at 1-800-873-7683.

HANDPIECE CONNECTION

Prior to connecting the handpiece, check that all tubing supplying the handpiece hose is properly connected with the dental unit. Make sure there are no restrictions in any of the lines. The handpiece and connector lines must align and nest properly. Operate the handpiece at 30-35 psi (measured at the connector). Increasing the air pressure will not increase performance and will prove detrimental to the turbine.

WARRANTY

***Warranty:** Patterson Dental Supply, Inc. warrants this handpiece to the original purchaser against defective workmanship and materials under normal installation, use and service for six months from the date of purchase. For warranty repair, send the product with copy of the original invoice to Patterson. Your product will be repaired promptly and returned to you.*

Manufactured for:
Patterson Dental Supply, Inc.
1031 Mendota Heights Road
Saint Paul, MN 55120

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