LuxaFlow Ultra

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English

Product description

LuxaFlow Ultra is a light-curing, flowable, radiopaque composite specifically for the repair of Luxatemp Ultra and other bis-acrylate materials, as well as for cavity Class III, IV and V restorations.

Indications for use

- Individualization and repair of bis-acrylate temporaries
- Small fillings (cavity Class III, IV and V)
- Minimally-invasive fillings
- Base under other filling materials

Limitations of use

Do not use the material if a dry working area or the recommended application technique is not possible.

Contraindication

- Do not apply the material to the exposed pulp.
- Do not use the material where there are known allergies to any of the components or contact allergies.

Patient target group

Persons receiving treatment as part of a dental procedure, for example:
Loss of tooth substance due to caries, trauma or development disorders

Intended user

Dentist

Notes for use

- Light-curing units should have an output of 450 nm and should be checked regularly.
- The light intensity should be at least 600 mW/cm².
- Position the light-curing unit as close to the material as possible.
- Do not remove the oxygen inhibited layer because this is required for the bond with the next layer.
- When using LuxaFlow Ultra in conjunction with other materials, adhere precisely to the applicable manufacturer's instructions.
- If the material is applied to the patient's mouth with a single use applicator, the single use applicator must be used for this one patient only for hygienic reasons.

Recommended use

Use as customization and repair material

Note: When customizing and repairing a freshly-made provisional, there is no need to roughen and a bonding agent is not required (continue with step 2).

Caution! When processing the temporary, grinding dust can get into the eyes and respiratory tract.

- ▶ Wear a face mask and protective goggles.
 - ▶ Make sure that grinding dust is properly vented – do not inhale.
1. Repairing a provisional that has been worn: Roughen the provisional slightly using a grinder or a sandblaster at the appropriate points. Afterwards apply a bonding agent (e.g. Luxatemp Glaze & Bond) in accordance with the manufacturer's instructions.
 2. Select the desired LuxaFlow Ultra shade. Remove syringe cap and attach a Luer-Lock Tip.
 3. Press the syringe and apply the material with the help of the Luer-Lock Tip, ensuring a maximum layer thickness of 2 mm.
 4. Light-cure each layer for 20 seconds at a light intensity of at least 1000 mW/cm² or 40 seconds at a light intensity of at least 600 mW/cm².
 5. Finish off the repaired region according to the usual methods.

Use as restoration material

1. Clean and dry the preparation.
2. Prepare the cavity in accordance with the general rules for adhesive technique.

Caution! There is a potential risk of damage to the pulp and pain for the patient if the pulp is not sufficiently protected.

- ▶ Protect pulp in deep cavities.
3. Carry out etching according to the preferred technique (e.g. Total Etch). Observe the manufacturer's instructions for the material used!
 4. Apply a bonding agent (e.g. Ecosite Bond). Follow the corresponding instructions for use!
 5. Select the desired LuxaFlow Ultra shade. Remove syringe cap and attach a Luer-Lock Tip.
 6. Press the syringe and apply the material with the help of the Luer-Lock Tip, ensuring a maximum layer thickness of 2 mm.
 7. Light cure each layer for 30 to 40 seconds.
 8. Remove excess, contour the filling surface and finish, e.g. with diamond finishing burs and flexible discs with progressively finer grit sizes.
 9. Use a polishing paste to polish the filling surface and neighboring surfaces. Use a polishing cup or polishing disc for polishing work.

Interactions

- Materials containing eugenol, moisture and oily air may inhibit polymerization at the contact area. Avoid materials containing eugenol, moisture and air containing oil!
- The use of cationic oral rinses, plaque indicators and chlorhexidine may lead to discoloration.
- Leftover methacrylate-based materials may influence the setting behavior of silicone materials.

Residual risks/side effects

There are no known side effects to date. The residual risk of hypersensitivity to components of the material cannot be ruled out.

Warnings/precautions

- For dental use only!
- Keep out of the reach of children!
- Avoid contact with the skin! In the event of accidental skin contact, immediately wash the affected area thoroughly with soap and water.
- Avoid contact with eyes! In the event of accidental contact with the eyes, immediately rinse thoroughly with plenty of water and consult a physician if necessary.
- Follow the manufacturer's instructions for other products that are used with the material/materials.
- Serious incidents involving this product must be reported to the manufacturer and to the responsible registration authorities.
- Wear mouth protection and goggles when working on the temporary, otherwise sanding dust may get into the eyes and respiratory tract!

Product performance characteristics

The product complies with DIN ISO 4049:2019 (Type 2, Class 2, Group 1) and has a polymerization depth of ≥ 2 mm at a curing time of 20 seconds and a light intensity of at least 1000 mW/cm² or with a curing time of 40 seconds and a light intensity of at least 600 mW/cm².

Storage/disposal

- Store in a dry place (15-25 °C/59-77 °F) at room temperature
- Disposal must comply with national regulations.
- Reseal syringes after use!
- Do not use after the expiration date!

Composition

Dental glass, bis-GMA, UDMA, HDDMA, DDDMA, SiO₂, additives. Contains TPO. Inorganic filling materials: approx. 41 vol% (0.02 – 3 µm).

Contains: 2-ethylhexyl 4-(dimethylamino)benzoate

Note: A glossary of all symbols is displayed at: <http://www.dmg-america.com/literature>

Caution: US Federal Law restricts this device to sale by, or on the order of a dentist, or other practitioner licensed by law of the state in which he or she practices to use or order the use of this device.

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