Carbide & Diamond Burs Instructions for Use

SYMBOLS GLOSSARY

	Consult Instructions For Use
	Warning or Precautions.
	Pre-Vacuum or Gravity Displacement Sterilizable in a Steam Sterilizer (Autoclave) at 132°C (270°F)
Rx Only	Caution: U.S. Federal law restricts the sale of this product only to qualified dental professionals
	Date of Manufacture
	Manufacturing Lot Number
	Part Number
	Package contains other recyclable material

GENERAL INFORMATION:

Diamond and carbide burs are rotary instruments used in conjunction with dental handpieces. They are medical devices intended for use on child and adult patients by trained dental professionals.

PRODUCT COMPOSITION:

Diamond dental burs: Bur head contains diamond particles embedded in nickel; shank body made of stainless steel. **Carbide dental burs:** Bur head is made from cemented tungsten carbide; shank is made up of either steel, stainless steel or tungsten carbide. Shank body maybe nickel and/or gold plated.

INDICATIONS FOR USE:

Diamond and carbide burs are used to cut, trim and finish natural and engineered dental materials during various dental procedures.

CONTRAINDICATIONS: None Known.

 **WARNINGS:**

- This product can expose you to nickel and should not be used on individuals with known allergic sensitivity to this metal.
-  **WARNING:** This product can expose you to chemicals including nickel, which is known to the State of California to cause cancer or reproductive harm. (For more information go to www.P65Warnings.ca.gov)
- Never use a bur beyond the manufacturer's recommended RPM.
- Always use water irrigation while using the device. Inadequate use of water and air irrigation may generate undesirable heat and cause patient discomfort, tissue necrosis, or patient burns.

PRECAUTIONS:

- Burs are sold non-sterile and must be sterilized prior to use.
- Use of a dental dam is highly recommended to avoid contamination of the prepared site with bacteria present in the oral cavity and to guard against the potential ingestion of dental materials or other debris.
- Ensure the air pressure to the handpiece does not exceed the manufacturer's recommended setting. Maintain the handpiece in good working order and ensure it is correctly lubricated, as specified by the manufacturer.
- Do not extend the bur from the handpiece chuck. A loose or extended bur could eject from the chuck or break and cause injury.
- Avoid using the bur with a heavy transverse-cutting pressure, a wedging, or levering action, as this may increase the risk of breakage.
- Always wear eye protection and other appropriate personal protective equipment (PPE) when using dental burs.
- Discard burs that show noticeable deterioration in performance.
- Do not use cold disinfectant solutions, as they contain oxidizing agents that will degrade the performance and strength of carbide burs.

ADVERSE EVENTS:

If a serious incident occurs with this medical device, report it to the manufacturer and to the competent authority for the country in which the user/patient is established.

STEP-BY-STEP INSTRUCTIONS:

- Insert the bur fully and securely tighten or latch the device into the handpiece chuck.
- Run the handpiece to ensure the bur is rotating prior to contact with the dental surface and while lifting away from the cutting location.
- Use light pressure and a brushing motion to achieve the best results.
- Use plenty of air and water irrigation while using this device.

CLEANING:

- The use of automated cleaning devices is not recommended for the cleaning of burs.

- Clean burs ultrasonically with the devices inserted in bur blocks or holders to prevent damage from rubbing or vibrating against each other or hard surfaces. An ultrasonic cycle of 5 minutes using a neutral-pH enzymatic ultrasonic detergent is recommended.
- If debris is still present after the ultrasonic cycle, clean the burs with a soft toothbrush and rinse under running tap water until all debris is removed.
- Immediately dry the burs thoroughly with an absorbent towel or paper tissue. Inspect and discard any devices that show signs of damage or corrosion.

STERILIZATION:

- All burs should be sterilized before use.
- After cleaning, dip the devices in an anticorrosion lubricant, such as surgical milk, prior to the sterilization cycle. Place the devices on a sterilization tray. Put the tray in an FDA-cleared or ISO 11607 compliant sterilization pouch and seal.
- Sterilize the devices in the sealed pouch using a pressurized steam sterilizer. For the pre-vacuum steam sterilization method, operate at full cycle with a dwell at 132°C minimum for 4 minutes. For the gravity displacement sterilization method, operate steam sterilizer at full cycle with a dwell at 132°C minimum for 15 minutes. Use distilled water in the steam sterilization process.
- Ensure that the devices go through the full drying cycle before they are removed from the sterilizer. The minimum drying time is 20 minutes for pre-vacuum and 15 minutes for gravity steam sterilization.
- Store sterilized burs, in the sterilization pouch, in a dry environment prior to subsequent use. Before using the devices, check for signs of rust or corrosion. Do not use any burs that are corroded or damaged.
- Use all sterilization devices according to the manufacturer's recommended procedures. It is the user's responsibility to ensure effective sterilization is achieved.

STORAGE AND DISPOSAL:

Store sterilized burs, in the sterilization pouch, in a dry environment prior to subsequent use. Before using the devices, check for signs of rust or corrosion. Do not use any burs that are corroded or damaged. For used burs, dispose of as biohazardous waste following local and national regulations.

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Patterson Dental Supply, Inc.