# Patterson® LC Monoshade Composite / Patterson® LC Opaquer

Instructions for Use

**Patterson® LC Monoshade Composite** is a light cure radiopaque composite for anterior and posterior restorations. **Patterson® LC Monoshade Composite** can be used for all natural tooth shades due to its continuous color matching properties. Besides its permanent gloss and toothlike opalescence, **Patterson® LC Monoshade Composite** shows excellent physical properties and good handling features. High gloss polishability demands only little effort. **Patterson® LC Monoshade Composite** contains dimethacrylate based resins and inorganic fillers < 0.2 µm. The total filler content is 80% by weight and 75% by volume.

**Patterson® LC Opaquer** is a light cure composite that is used in combination with **Patterson® LC Monoshade Composite** to mask discolored teeth. **Patterson® LC Opaquer** contains dimethacrylate based resins and inorganic fillers < 0.2 µm. The total filler content is 80% by weight and 76% by volume. **Patterson® LC Opaquer** is radiopaque.

The material properties of **Patterson® LC Monoshade Composite** and **Patterson® LC Opaquer** meet the requirements of **ISO 4049**, type 1, class 2, group 1.

## Indications Patterson® LC Monoshade Composite

- Restorations of all cavity classes
- Direct veneers
- Diastema closure
- Repair of ceramic restorations (e.g. veneers, crowns)

## Indications Patterson® LC Opaquer

- Masking of the oral cavity in case of extended class III and IV cavities
- Masking of tooth discoloration

## Side Effects

Irritations resulting from direct contact with the pulp cannot be ruled out. Therefore for pulp protection areas close to the pulp should be covered with pulp capping material

## Incompatibility with Other Materials

Do not use in combination with substances containing eugenol because eugenol inhibits the polymerization of the composite. Neither store the composite material in proximity of eugenol containing products, nor let the composite allow coming into contact with materials containing eugenol.

## Preparation of Tooth

Clean the tooth with a fluoride-free polishing paste (e.g. flour of pumice with water) prior to preparation and color matching.

## Cavity Preparation

After isolation (best with a rubber dam) prepare the cavity with minimal tooth reduction.

Pulp protection areas close to the pulp should be covered with a pulp capping material. Refer to the corresponding user instructions.

Apply a common adhesive for enamel and dentine bonding according to the corresponding manufacturer instructions.

The application of the universal adhesive is recommended. This can be used in combination with an etching gel or as a self-etching adhesive. Follow the instructions for product.

## Application and Curing of Patterson® LC Monoshade Composite

In case of deep and extended cavities, best results are obtained with a build-up of a base liner using a light cure resin modified cement or a light cure flowable composite, followed by the application of **Patterson® LC Monoshade Composite**.

In case of small cavities **Patterson® LC Monoshade Composite** can be applied directly.

## Note:

In case of extended class III and class IV cavities use **Patterson® LC Opaquer** palatal to ensure that the restoration may not appear dark.

Use **Patterson® LC Opaquer** also as a first layer for class I, II and V cavities of dark or discolored teeth.

In case of repair of ceramic or composite restorations, use a light cure opaquer as first layer in case of metal or zirconia surfaces.

Apply **Patterson® LC Opaquer** in layers of max. 1 mm. Light cure each increment separately.

By using a polymerization unit (wavelength range 400-500 nm) with a light intensity of at least 1000mW/cm<sup>2</sup>, cure each increment separately for 20 seconds.

Hold the light emission window as close as possible to the filling material.

Apply **Patterson® LC Monoshade Composite** in layers of max. 2 mm. Light cure each increment separately. **Patterson® LC Monoshade Composite** is easy to model. Place it by using a suitable instrument. Matrix strips may be used.

By using a polymerization unit (wavelength range 400-500 nm) with a light intensity of at least 1000mW/cm<sup>2</sup>, cure each increment separately for 20 seconds.

Hold the light emission window as close as possible to the filling material.

## Finishing

Depending on the shape of the area, contour with finishing diamonds, flexible grinding discs or carbide burs.

Polish to high gloss with polishing discs or polishing brushes.

## Storage

Do not store above 25°C (77°F). Protect from direct sunlight. Do not use after expiry date.

Close syringe immediately after use to avoid exposure to light.

## Additional Notes/Warnings

- The ambient light of the dental lamp may start polymerization of the composite.
- Do not use any resin to adjust viscosity of composite restorative material.
- Avoid contact with skin, mucous membrane and eyes.
- Unpolymerized composite may have an irritant effect and can lead to sensitization against methacrylates.
- Commercial medical gloves do not protect against the sensitizing effect of methacrylates.

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**Patterson Dental Supply, Inc.