
V-Print cast

Instructions for use

Product description:

V-Print cast is a light-curing resin for the generative production of burn-out objects for casting and pressing processes.

Indications:

– Production of objects which can be burned out without leaving any residues for casting processes and press ceramics.

Preparation:

V-Print cast is a light-curing resin designed for the generative production of objects using a 385 nm LED DLP printer. Information on approved printers is available from the corresponding device manufacturers and from **VOCO**.

Prepare a print job using slicing software. In addition to the material-dependent construction specifications in these instructions for use, please also observe the dependencies on the positioning, type of support and adaptation from our other documents for your construction.

V-Print cast can be printed in the following layer thickness:

– 50 µm

Processing:

Note: Use separate material containers and cleaning baths for each printing material so as to rule out the possibility of cross-contamination.

The material container should be filled directly before printing is started. When doing so, ensure the material is as free from bubbles as possible and pay attention to the fill level.

Note: Do *not* shake the material before starting printing.

Start the print job observing the previously selected parameters.

We recommend allowing a dripping time of approx. 10 minutes once printing is complete. Then detach the printed objects carefully from the build platform.

Following this, the printed objects must be cleaned, dried and post-exposed in order to guarantee the requisite product characteristics. A detailed description of the two steps out-lined above can be found under Post-processing. Post-process the printed object without delay so as to avoid distortion of its shape and ensure the required precision.

Recommendation: Remove the excess material from the material container following printing. The printer manufacturer’s specifications must be observed and followed, so as to allow you to clean and check the container.

The print material can be returned to its original packaging. If necessary, a separate container can also be used (HDPE, not light-transmissive, airtight).

Please ensure that the excess material intended for further use is free from contaminants and polymerised components. To this end, a stainless steel sieve or similar instrument should always be used when transferring material.

Post-processing:

Cleaning:

We recommend a solvent-resistant cleaning mat for cleaning.

Cleaning should be performed with isopropanol (≥ 98% purity) in an unheated ultrasound bath.

The respective operating manuals and/or instructions for use provided by the manufactureres must be observed.

The printed objects must be cleaned in two or optionally in three stages:

Pre-cleaning:

– Optionally: Pre-clean printed objects carefully for approx. 15 seconds by repeatedly immersing them in a beaker of isopropanol. When doing so, hold the component by a support using tweezers.

– Pre-clean the printed object for 3 minutes in a reusable isopropanol ultrasound bath. Position uncleaned printed objects in the ultrasound bath in such a way that any openings face downwards.

Note: The bath’s cleaning efficacy decreases the more it is used. Replace the bath when the cleaning efficacy decreases.

Main cleaning:

Clean the printed objects for a further 2 minutes in a fresh isopropanol ultrasound bath.

Then dry them carefully with compressed air.

If there are any resin residues on the printed object after the main cleaning or any exits from the undercuts following drying, the printed object can be immersed in the main cleaning bath again briefly, before being dried again.

Preparation for post-exposure:

The support structures do not need to be removed before post-exposure.

Inconvenient support structures can be removed prior to the post-exposure carefully and without exerting pressure using a rotary instrument (e.g., cutting wheel) as close to the printed object as possible. Use a suction device. Remove any remaining acrylic dust carefully using compressed air. Then rinse printed objects for a few seconds with fresh isopropanol and remove any dust using a fine, clean brush if necessary. Dry the printed objects carefully again using compressed air.

We recommend waiting 15 minutes after the last contact with isopropanol before initiating the post-exposure.

Note: Ensure that the printed objects are not exposed to any direct light during this time.

Post-exposure:

No protective gas atmosphere is required. However, ensure that the printed objects do not cast shadows over or touch each other, as this could result in incomplete post-curing.

The post-exposure can be performed with the following devices:

– **Otoflash G171 xenon photoflash unit** **2 × 2,000 flashes**
– **LC-3DPrint Box UV light box** **30 minutes**

In printed objects post-exposed with the Otoflash G171, a cooling-off period of at least 2 minutes with the lid open must be included after the first 2,000 flashes. Once the cooling off period has elapsed, the printed objects should be turned and light-cured with a further 2,000 flashes.

Insufficient cooling can result in distortion of the shape of the printed object.

Finishing:

As a general rule, work with a low contact pressure and speed. This guarantees constant results and also reduces the risk of undesirable signs of processing.

Use a fine-toothed carbide bur or similar to sand the support stubs. The same bur can also be used for subsequent finishing of special structures.

For precise sanding, e.g., between the support strut and the printed object, it is recommended to process the surface in the corresponding areas with sandpaper with different grit sizes. A corresponding result can also be achieved with coarser or finer silicone polishers.

Embedding and casting:

Only use phosphate-bonded bedding compounds for production of the casting mould. A list of recommended bedding compounds can be found at www.voco.dental.

The recommendation is direct mixing of the embedded objects at the final temperature (speed operation).

Follow the bedding compound and device manufacturers’ instructions for use.

Information, precautionary measures:

– Observe the **V-Print cast** safety data sheet!

– Pay attention to the finishing process

V-Print cast contains (meth)acrylates. **V-Print cast** should not be used in case of known hypersensitivities (allergies) to the ingredients.

Contact between uncured **V-Print cast** and the skin/mucous membranes and eyes can cause mild irritation and should be avoided. Wearing protective clothing is recommended. In addition, care must be taken not to inhale any vapours and/or dusts. The wearing of a suitable mask and/or the use of extraction systems is recommended. Further information on the handling of the product can be found on the safety data sheet.

Storage:

Store at between **15°C and 28°C**. Seal bottle again immediately after use. The material will cure if exposed to light. Do not use the product after its expiry date.

Our preparations have been developed for use in dentistry. As far as the application of the products delivered by us is concerned, our verbal and/or written information has been given to the best of our knowledge and without obligation. Our information and/or advice do not relieve you from examining the materials delivered by us as to their suitability for the intended purposes of application. As the application of our preparations is beyond our control, the user is fully responsible for the application. Of course, we guarantee the quality of our preparations in accordance with the existing standards and corresponding to the conditions as stipulated in our general terms of sale and delivery.

This material has been developed solely for use in dentistry. Processing should be done strictly according to the instructions for use. VOCO recognizes its responsibility to replace products if proven to be defective. VOCO does not accept liability for any damage or loss, directly or indirectly, stemming from the use of or inability to use the products described. Before using, it is the responsibility of the user to determine the suitability of the product for its intended use. The user assumes all risk and liability in connection therewith. Descriptions and data constitute no warranty of attributes and are not binding.

CAUTION: U.S. Federal Laws restrict this device to sale by or on the order of a dentist.

No person is authorized to provide any information which deviates from the information provided in the instructions for use.

For questions or comments, please call 1-888-658-2584.

Keep this material out of reach of children.

For dental use only.

An explanation of the symbols used in labeling can be found at www.voco.dental/symbols

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