Dental Zirconia

UTML (Ultra Translucent Multi Layered)

[ENGLISH] INSTRUCTIONS FOR USE

I. Introduction

This IFU is for KATANA Zirconia UTML (Ultra Translucent Multi Layered). KATANA Zirconia UTML is a pre-sintered zirconia disc 98.5 mm in diameter, which contains a plastic ring. This is designed for all milling systems using this generic-type disc. (Please refer to your milling system's technical instructions for correct machine operation.) KATANA Zirconia UTML has 2 available thicknesses: (14mm and 18mm) and has 20 shade variations: (A1, A2, A3, A3.5, A4, B1, B2, B3, B4, C1, C2, C3, C4, D2, D3, D4, EA1, EA2, EA3, ENW). KATANA Zirconia UTML consists of 4 graded shade layers. KATANA Zirconia UTML is recommended for use in fabricating FCZ (Full Contour Zirconia) restorations or the frameworks.

II. Intended Use

KATANA Zirconia is used for the fabrication of the all-ceramic restorations (frameworks, FCZ crowns, FCZ bridges, inlays, onlays and veneers).

III. Sintering Program

Temperature	Programming Rate	Holding Time
Room Temp. — 1550°C (2822°F)	10°C/min. (18°F/min.)	—
1500°C (2822°F)	—	2 hrs
1550°C (2822°F) — Room Temp.	-10°C/min. (-18°F/min.)	—

Temperature	Programming Rate	Holding Time
Room Temp. — 1400°C (2552°F)	50°C/min. (90°F/min.)	—
1400°C (2552°F) — 1500°C (2732°F)	4°C/min. (7°F/min.)	—
1500°C (2732°F) — 1560°C (2840°F)	10°C/min. (18°F/min.)	—
1560°C (2840°F)	—	16 min
1560°C (2840°F) — Room Temp. (*)	-50°C/min. (-90°F/min.)	—

* The restorations may be removed at 800°C (1472°F) or less depending on the circumstances.

Temperature	Programming Rate	Holding Time
Room Temp. — 1450°C (2642°F)	120°C/min. (216°F/min.)	—
1450°C (2642°F) — 1600°C (2912°F)	10°C/min. (18°F/min.)	—
1600°C (2912°F)	—	20 min
1600°C (2912°F) — Room Temp. (*)	-120°C/min. (-216°F/min.)	—

* The restorations may be removed at 800°C (1472°F) or less depending on the circumstances.

IV. Composition

ZrO₂, Y₂O₃ etc.

V. Type and Class (ISO6872:2015)

Type:II/ Classe:3

VI. Physical Properties

Coefficient of Thermal Expansion (25-500°C (77-932°F)): 9,7×10⁻⁶/K

VII. Directions for Use

- Take the disc from the packaging and confirm that the disc does not have a crack or other damage.
- Place the disc into the milling machine; then begin the milling process following the milling systems technical instructions.
- After milling, remove the restorations from the disc with a diamond bur, etc.
- The cutting waste or dust, which is attached to the restorations, can be removed with a gentle air stream.

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- The margins should be prepared with a deep chamfer and rounded shoulders, with cutting edges and corners rounded to eliminate sharp preparation corners. The angle of the axial surface should be within the range of 5 to 15 degrees.

- When preparing teeth, avoid the following: deep shoulders, J-margins, knife edges, serrated margins, non-tapered abutments, undercuts, guide grooves, the formation of retentive holes, and sharp corners.

- Keep the following thickness of this product for fabricating prosthetics:

Location & indication	Wall thickness
Anterior crown or bridge	0.8 mm or more
Veneer	0.4 mm or more*
Posterior crown or bridge	1.0 mm or more
Inlay or onlay	1.0 mm or more

* 0.4 mm or more of this product is for full zirconia veneers. Keep thickness 0.8 mm or more, if it is used for combination with the porcelain.

- Use the following cross-sectional areas for connectors when fabricating bridges.

Location & indication	Connector cross section
Anterior 2- or 3-unit bridges	12 mm² or more
Premolar 2- or 3-unit bridges	16 mm² or more

- Choose a shade color that is brighter than the intended color for a thick restoration, as it may look duller depending on the thickness of the restorations.
- When using a sintering furnace for the first time and changing a sintering condition, colors after sintering may vary. Sinter a small piece of Zirconia beforehand and confirm the color.

Caution in conjunction with usage:

- Do not use the disc if there is a crack noticed after removing it from the package.
- Dry milling is recommended. If wet milling is used, the translucency of the disc may be reduced.
- If you find a crack in the restorations, do not use.
- The above sintering recommendation is only a guideline; some adjustments may be required depending on each individual furnace.
- When milling the disc, use caution when approaching the milling bur with the plastic ring to prevent detachment of the disc. Mill the disc as if leaving a zirconia part contacting the internal side of the plastic ring at 2 mm.
- Cut and remove the plastic ring prior to sintering the whole disc.
- Dispose of this product as a medical waste to prevent infection.

Storage:

- Store in a cool and dry place. Keep away from direct sunlight.
- The product should be stored at 10-30°C (50-86°F)
- Do not remove the disc from its packaging during storage.
- The disc is fragile, and requires care when handling.
- The product must be stored in an appropriate place where only dental personnel have access.
- The product must be used by the expiration date indicated on the package.

[WARRANTY]

Kuraray Noritake Dental Inc. will replace any product that is proven to be defective. Kuraray Noritake Dental Inc. does not accept liability for any loss or damage, direct, consequential or special, arising out of the application or use of or the inability to use these products. Before using, the user shall determine the suitability of the products for the intended use and the user assumes all risk and liability whatsoever in connection therewith.

[NOTE]

If a serious accident attributable to this product occurs, report it to the manufacturer's authorized representative shown below and the regulatory authorities of the country in which the user/patient resides.

*KATANA™ and ™CERABIEN™ are trademarks of NORITAKE CO., LIMITED.

	Medical Device
---	----------------

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