

# Cetylcide-G<sup>®</sup>

High-Level Disinfectant/Sterilant Concentrate and Diluent



Cetylcide-G is a pleasantly-scented glutaraldehyde-based sterilant/high-level disinfectant. The tip-and-pour bottle measures either 1 quart or 1/2 gallon of use-solution for a total yield of 15 quarts (3.75 gallons). This uniquely designed system allows you to mix only what you need; therefore eliminating waste as with products requiring you to activate an entire quart or gallon at a time. The 3.2% strength use-solution may be reused for up to 28 days at 20°C (68°F). An 8 oz bottle of Diluent is supplied with each quart of Cetylcide-G.



Manufacturers Suggested Retail Price per Quart	\$60
Yields of Use Dilution per Quart	15 quarts (3.75 gallons)
Retail Cost per Unit of Use-Dilution	\$4.00 per quart or \$16.00 per Gallon
Retail Cost per Unit for Leading Brand	\$19.00 per quart or \$30.99 per Gallon
Retail Cost per Unit for Private Label	\$11.49 per quart or \$26.99 per Gallon

Description	Cetylcide-G High Level Disinfectant /Sterilant & Diluent
Size	1Qt/8oz
Item No.	0122
MSRP	\$60.00

## Key Features

- Activate small quantities as needed
- Effective at room temperature (20°C/68°F)
- Approx 8 oz of Cetylcide-G concentrate makes 1 Gallon of ready to use solution
- 1 quart makes 3.75 Gallons of use solution
- Tip-and-pour bottle
- Up to 28 day use-reuse
- Can be formulated with ordinary tap water
- Anti-corrosion formula
- Dilutes to 3.2% glutaraldehyde when mixed as directed

## Benefits

- Economical - only \$4.00 for a quart of use solution
- No waste - only use what you need
- Save inventory space
- No need to elevate solution temperature
- No mess dispensing
- Simplifies disinfection
- Allows for manual processing
- Compatibility
- Compliant
- Very little goes a long way



[www.cetylite.com](http://www.cetylite.com)

### Contact:

Cetylite Industries, Inc.  
9051 River Road  
Pennsauken, NJ 08110  
800.257.7740  
marketingdept@cetylite.com

# How to Prepare Cetylcide-G Use Solution

(See Insert for Detailed Instructions)



Fig. 1  
Fill C-Tub or other appropriate container with 1 quart, 1/2 gallon, or 1 gallon of water, depending on your needs.



Fig. 2  
With both caps closed, tilt the 8 oz Cetylcide-G Diluent bottle forward and fill the measuring chamber.



Fig. 3  
Using the angular pour-back scale, adjust to the amount needed by tilting the bottle back.



Fig. 4  
Hold bottle upright and check the fluid level for accuracy. See Fig. 5 for amounts needed.

Fig 5

Cetylcide-G Diluent	Amount of Water
15 mL	1 quart
30 mL	1/2 gallon
60 mL	1 gallon



Fig. 6  
Remove the Pour Cap and pour Cetylcide-G Diluent into the C-Tub or other appropriate container.



Fig. 7  
Using the same technique as demonstrated above, fill the Cetylcide-G Concentrate measuring chamber. See Fig. 8 for amounts needed.

Fig 8

Cetylcide-G	Water/Diluent Mixture
63 mL	1 quart
126 mL	1/2 gallon
252 mL	1 gallon



Fig. 9  
Remove the Pour Cap. Pour Cetylcide-G Concentrate into the Water/Cetylcide-G Diluent mixture. Fill-in and apply the supplied green Use-Solution label to the container.

## How to Use Cetylcide-G Use Solution

### For use as a high-level disinfectant

Immerse instruments in Cetylcide-G Use Solution for 40 minutes at 20°C.

### For use as a Sterilant

Immerse instruments in Cetylcide-G Use solution for 10 hours at 20°C. Remove instruments from the solution and rinse thoroughly with sterile water following the rinsing instructions.

### Rinsing

Immerse instruments completely in three separate copious volumes of water (sterile water if taken from sterilant). Each rinse should have a minimum of 1 minute in duration unless otherwise noted by the device manufacturer. Use fresh portions of water for each rinse. Discard the water following each rinse. Do not reuse the water for rinsing or any other purpose as it will be contaminated with glutaraldehyde. If your state requires neutralizing of glutaraldehyde before disposing, please ask about our Hyde-Out neutralizer product.

### When using Cetylcide-G Use Solution



DO NOT TOUCH



DO NOT SPRAY



DO NOT WIPE



Member