Patterson® DC Nano-reinforced Self-etch Bond

Carefully read instructions prior to use

Patterson® DC Nano-reinforced Self-etch Bond is a dual-curing, self-etching bond reinforced with nanoparticles for creating a durable bond between tooth substance and light-, self- and dual-curing restoratives without marginal leakage. This all-in-one adhesive consists of 2 components, which eliminates the need for a separate etching step, yet still achieves the same adhesive properties as total-etch bonds. **Patterson® DC Nano-reinforced Self-etch Bond** tolerates residual moisture and can be used universally with all light-, self- and dual-curing composites compomers and Omrocers®.

Indications:

- direct self- or dual-curing composite restorations and core build-ups
- direct light-curing composite/compomer-/Omocer®-based restorations
- fixing root pins/posts with dual- or self-curing composite cements
- indirect restorations; with the application of dual-curing/self-curing composite cements for bonding inlays, onlays, crowns and bridges

Contraindications:

Patterson® DC Nano-reinforced Self-etch Bond contains organic acids, BIS-GMA, HEMA, TMPTMA, campherchinon, amines (DABE), BHT, catalysts, fluorides and ethanol. Do not use in patients with known allergies to these ingredients.

Warnings and precautionary measures:

- Phenolic substances (e.g. temporary materials containing eugenol) disturb polymerization and therefore should not be used with **Patterson® DC Nano-reinforced Self-etch Bond**.
- Avoid contact with the oral mucosa. Contact of the oral mucosa with **Patterson® DC Nano-reinforced Self-etch Bond** may result in a whitish discoloration caused by the coagulation of proteins. This is a reversible occurrence that subsides within a few days. In case of contact with the eyes, rinse with plenty of clean water and consult a physician. Application brushes provide precise and thin application of **Patterson® DC Nano-reinforced Self-etch Bond**.

Application:

1. Preventing contamination

A dry working field is essential due to the subsequent application of the restorative. Avoid contamination with blood or saliva. Using a rubber dam is recommended.

2. Preparation

Prepare the cavity according to the principles of the adhesive restoration technique. Rinse cavity thoroughly with water. Remove excess moisture with a faint air jet. Do not overdry dentin. The dentinal surface of the cavity should remain slightly moist. If the tooth has not been prepared (e.g. cervical defect) the surface must be thoroughly cleaned and enamel margin roughened with a finishing diamond. Areas in close proximity to the pulp should be protected with a suitable calcium hydroxide-based liner.

3. Mixing

Patterson® DC Nano-reinforced Self-etch Bond is dual-curing and therefore must be applied immediately after mixing. Mix 1 drop of Liquid 1 and 1 drop of Liquid 2 on a mixing palette with an application brush for approximately 2 sec (produces a self-etching adhesive).

4. **Patterson® DC Nano-reinforced Self-etch Bond** application

Apply the adhesive in a layer of medium thickness to the enamel/dentin and rub into the tooth surface for 20 sec.

Caution: The mixed bonding is light-curing and thus an intensive exposure to ambient light should be avoided. The OP-Light during the application should be dimmed.

5. Drying

Dry the adhesive layer for at least 5 sec with an air syringe.

6. Polymerization

6.1. Direct restoration

Light-, self- and dual-curing materials

The bonding layer must be polymerized with blue light (halogen or LED lights) for 10 sec when using light-, self- and dual-curing restoratives.

6.2. Indirect restoration

Light-curing luting composites

The bonding layer must be polymerized with blue light (halogen or LED lights) for 10 sec when using light-curing luting materials.

Self- or dual-curing luting materials

When using dual- or self-curing luting materials, **Patterson® DC Nano-reinforced Self-etch Bond** also provides excellent adhesion without light-curing. **Patterson® DC Nano-reinforced Self-etch Bond** can also be light-cured for 10 sec with bluelight (halogen or LED lights) as an optional step.

Luting root posts

- Massage the mixed **Patterson® DC Nano-reinforced Self-etch Bond** (liquid 1 and 2) into the occlusal surface of the core and root canal for 20 sec with an applicator. Disperse solvent with an oil-free air-jet for 5 sec. Remove excess material by dabbing it off with an absorbent strip of paper. Do not light-cure the bonding layer.
- Having prepared the root post in accordance with the manufacturer's instructions, insert it into the root canal using a suitable luting material. (see respective instructions for use). The aim is to obtain a slight excess of cement.
- Light-cure for at least 40 sec to fix the root post and cure the **Patterson® DC Nano-reinforced Self-etch Bond** applied to the occlusal surface.

7. Advice on luting Maryland Bridges

For luting Maryland bridges, the enamel must either be roughened or conditioned with a phosphoric acid etching gel before applying **Patterson® DC Nano-reinforced Self-etch Bond**. Self-etching bonding materials are less effective on unprepared enamel. Overapplication of bonding material to enamel that has not been ground or bevelled can lead to margin discoloration. For optimum adhesion results with **Patterson® DC Nano-reinforced Self-etch Bond**, the enamel should be prepared in advance in class III, IV, V cavities and with diastema closures.

Advice for storage and usage:

Tightly close **Patterson® DC Nano-reinforced Self-etch Bond** bottles after use. Store in an upright position to guarantee flow-back of the liquid. Store at temperatures between 39°F - 73°F (4°C - 23°C). Allow **Patterson® DC Nano-reinforced Self-etch Bond** to reach room temperature before use. Do not expose **Patterson® DC Nano-reinforced Self-etch Bond** to direct sunlight or the chair light. Do not use **Patterson® DC Nano-reinforced Self-etch Bond** after the expiration.

This material has been developed solely for use in dentistry. Processing should be done strictly according to the instructions for use. Patterson® recognizes its responsibility to replace products if proven to be defective. Patterson® does not accept liability for any damage or loss, directly or indirectly, stemming from the use of or inability to use the products described. Before using, it is the responsibility of the user to determine the suitability of the product for its intended use. The user assumes all risk and liability in connection therewith. Descriptions and data constitute no warranty of attributes and are not binding.

CAUTION:

U.S. Federal Laws restrict this device to sale by or on the order of a dentist.

No person is authorized to provide any information which deviates from the information provided in the instructions for use.

**Keep this material out of reach of children.
For dental use only.**

Manufactured for:
Patterson Dental Supply, Inc.
1031 Mendota Heights Road
Saint Paul, MN 55120

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