
LED curing light plus
INSTRUCTIONS FOR USE

**E
N
G
L
I
S
H**

Thank you for purchasing the Patterson Dental Inc., LED curing light. This curing light is a powerful, cordless, lightweight curing light using an LED for the polymerization of composite restorative materials that contain the photoinitiator Camphorquinone. It produces a strong blue light in the wavelength range of 440 – 480 nm, the relevant range for Camphorquinone containing products. It features a 'NO RAMP' mode and 'RAMP' mode for convenient polymerization shrinkage minimization of composite materials. In addition to this, the LED attachments can be changed depending on the clinical application.

PRODUCT CONTENTS LIST

- Handpiece (battery pack, control section and standard LED attachment)
- Charger
- Plug pack
- 100 barrier sleeves
- 3 lens caps
- Radii Plus light shields (5 Pack)

UNPACKING AND INITIAL CHARGING

The curing light has been carefully packed to withstand any damage during transit. Ensure all parts contained within the box match the product contents list above.

- 1) Remove all parts from the box.
- 2) Place the plug pack cable into the charger.
- 3) Plug the plug pack into an available power outlet and turn power outlet on. The charger will conduct a brief LED test: the charge light LED should momentarily show red/green and the 5 blue radiometer LEDs should turn on/off.
- 4) Place the curing light handpiece into the charger. The handpiece can be rotated 360° in the charger. Depending on the length of time from manufacture to package opening, the charge remaining in the battery can vary. The charge light will show the below three conditions.

Charge Light Color	Battery condition
Red	discharged
Flashing Green	charging (OK to use)
Green	charged

To ensure the highest level of performance, allow the curing light to fully charge prior to first use. On average, it will take approximately 1-3 hours to initially charge.

OPERATION

• Turning unit on and off

Simply pick up the curing light handpiece, press the on/off button once to activate and de-activate the light. Audible beeps facilitate accurate time measurement of curing. Replace the handpiece back in the charger when not in use.

• NO RAMP mode and RAMP mode.

The curing light has two curing modes, NO RAMP and RAMP mode. NO RAMP mode will output 100% intensity immediately after pressing the on/off button whereas RAMP mode will output a reduced intensity for the first 5 seconds before irradiating at 100% intensity.

• Switching between the NO RAMP mode and RAMP mode.

The curing light has been set at the factory to the NO RAMP mode. To switch between the NO RAMP mode and RAMP mode, when the unit is off, hold down the on/off button for 3 seconds. Two audible beeps will be heard. The unit will now be programmed for RAMP mode. To switch back to the NO RAMP mode, hold down the on/off button for 3 seconds, a single audible beep will be heard indicating the unit has switched back to NO RAMP mode. NOTE: In RAMP mode the visible light starts slightly dimmer and progressively gets stronger, as opposed to in NO RAMP mode, where the visible light is constant from the start.

RAMP mode indicator:

This unit has been fitted with a blue indicator light above the on/off button. When in RAMP mode, the indicator light will show blue. The indicator light is off in the NO RAMP mode.

• Audible beeps for NO RAMP and RAMP modes:
Depending on the mode set (NO RAMP or RAMP), the beeping cycle will differ. Below is a table explaining the audible beeps for each mode.

Audible beeps	Time between beeps (seconds)		Total elapsed time (seconds)	
	NO RAMP	RAMP	NO RAMP	RAMP
1 short beep	10 sec	15 sec	10 sec	15 sec
2 short beeps	10 sec	10 sec	20 sec	25 sec
3 short beeps	10 sec	10 sec	30 sec	35 sec
4 short beeps	10 sec	10 sec	40 sec	45 sec
5 short beeps	10 sec	10 sec	50 sec	55 sec
1 long beep	10 sec	10 sec	60 sec	65 sec

At the end of 60 seconds (or 65 seconds in case of "RAMP" mode) one long beep will be heard before the blue LED is automatically switched off. NOTE: The above cycle is based on the Standard LED attachment supplied in this kit. Other access types of LED attachments will have different cycles. Refer to instruction sheets associated with those products for beeping cycle.

• Rotating and changing the LED attachment

The LED attachment, can be rotated 360° to ensure accurate and detailed polymerization in any area of the mouth. WARNING: Do not rotate the LED attachment at the LED end, instead turn at the end closer to the control section as shown below:

• Lens Cap

The curing light has been supplied with a lens cap fitted. The lens cap should be inspected prior to each use for scratches, cracks or foreign matter that may reduce the effectiveness of the LED light source. Generally lens caps should be replaced every 2 weeks. To remove the lens cap carefully unscrew the old lens cap counter-clockwise and fit the new lens cap by tightening clockwise.

• Built in Radiometer

The curing light's built-in radiometer provides a simple test to check the functionality of the unit. With the curing light on, carefully position the lens cap on the built-in radiometer. 4 or 5 lights should illuminate indicating a fully functional unit. If only 1, 2 or 3 lights illuminate, it could be due to the following reasons:

- 1) There could be some material on the lens cap such as composite: - the lens cap should be replaced.
- 2) The lens cap could be damaged or cracked: - it should be replaced.
- 3) The battery could be extremely low - the unit should be recharged.
- 4) In a severe case the light could be malfunctioning: - the unit should be serviced.

Note: Portable radiometers can give inaccurate and inconsistent results. These results can vary by as much as 300 mW/cm² or more.

• Disposable barrier sleeves

Transparent, disposable barrier sleeves are supplied to protect the handpiece and prevent cross contamination. Testing showed no significant difference in spectral and energy output when operating with or without a barrier sleeve in place. Slide the barrier sleeve over the end of the handpiece before placing the orange light shield on the end of the handpiece.

It is essential to use disposable barrier sleeves to prevent any liquids getting on the handpiece. Evidence of liquids entering handpiece causing damage due to not using barrier sleeves will void warranty.

• Automatic shut off function using the standard LED attachment

Should the standard LED attachment become too hot due to excessive use (after about 4 minutes of continuous irradiation), the unit will automatically shut off and be unusable for about 1 minute until it cools down. If using a different LED attachment, please refer to the instruction sheet for that attachment's automatic shut-off sequence as it may vary.

• Changing the LED Attachment

The curing light allows you to change the LED attachment depending on the curing situation. To change the LED attachment, pull off the LED attachment carefully by placing fingers as shown below. Please refer to LED attachment instructions for use for operation guidelines.

• Changing the battery pack

Carefully unscrew the battery pack counter-clockwise to remove and replace. A single beep will be heard when the battery pack is successfully connected.

• Cleaning the curing light:

The ideal way to protect patients from cross-infection is to utilize curing light barrier sleeves. Barrier sleeves are also essential to keep the curing light clean. When cleaning the curing light handpiece, ensure lens cap, battery pack, control section and LED attachment are connected as one unit.

- **LIQUID CONTACT WITH TERMINALS OF HANDPIECE AND CHARGER WILL CAUSE DAMAGE VOIDING WARRANTY.**
- **LIQUID CONTACT WITH LED WILL CAUSE DAMAGE VOIDING WARRANTY.**
- **LIQUIDS MUST NOT ENTER ANY PART OF THE HANDPIECE OR CHARGER.**
- **DO NOT SPRAY LIQUIDS DIRECTLY ONTO, OR APPLY COPIOUS AMOUNTS OF LIQUIDS TO ANY PART OF THE CURING LIGHT.**
- **DO NOT AUTOCLAVE.**

General Surface Cleaning

General surface cleaning can be conducted with anti-microbial surface disinfectant such as Glutaraldehyde, Chlorhexidine gluconate, and 70% Isopropyl alcohol. When applying the disinfectant to any component of the device, spray the disinfectant agent onto a piece of cloth and wipe over the external areas to be cleaned.

• Curing Times

Curing times will differ for different formulations of composite restorative materials. Follow the manufacturers instructions for recommended curing times. Always bench test new materials before use in-vivo.

• Disposal of battery pack

As a means of protecting the environment, the curing light battery pack is fitted with a Lithium Ion battery. Dispose of battery pack in accordance with local legal regulations.

TROUBLE SHOOTING

Problem	Causes	Solution
• Charge light flashes red	Possible battery fault	Replace battery
• Handpiece blue LED flashes then goes off	Battery charged too low	Re-charge battery
• audible beep heard but no light irradiates	LED attachment is not correctly inserted	Correctly insert LED attachment (see below)

- No lights appear on the charger when the hand piece is placed on the charger.
 1. Check terminals on the battery are clean and free of debris.
 2. Check power is getting to charger, turn power outlet on. The charge light LED should momentarily show red/green and the 5 blue radiometer LEDs should turn on/off.
 3. Check the power outlet is turned on, if so use a different power outlet.
 4. Check the plug pack is correctly connected to the charger.
 5. Check that the handpiece is placed into the charger correctly.
 - 6 If the problem still exists after checking 1, 2, 3 and 4 above the charger may be faulty. Return the unit to the supplier for service.

When returning faulty units please make sure that all parts are returned together in the original packaging.

SAFETY AND PRECAUTIONS

The safe operation of this unit depends on strictly following the operating instructions in this manual. The Manufacturer accepts no liability for any damage resulting from the use of this unit for any other purpose than the polymerization of dental composite materials.

CAUTION! Irradiation of the eyes bears an inherent health risk. Hence, the light must not be directed towards the eyes. Exposure must be restricted to the area of the oral cavity in which clinical treatment is intended. Suitable blue-light filtering safety goggles should be used. The unit emits a high intensity light and intensive light exposure of soft tissues (gingiva, oral mucosa and skin) should be avoided as this exposure may cause damage or irritation. The emitted light should be placed directly above the material to be cured. If applicable cover soft tissue areas.

Do not use the curing light in patients or by users with:

- heart pacemaker implants who have been advised to be cautious in regards to their exposure to small electrical devices.
- a history of photo biological reactions (including individuals with urticaria solaris or erythropoietic protoporphyria) or who are currently on photo-sensitizing medication (including 8-methoxypsoralen or dimethylchlorotetracycline).
- a history of cataract surgery. These people may be particularly sensitive to the exposure to light and should be discouraged from curing light treatment unless adequate safety measures, such as the use of protective goggles to remove blue light, are undertaken.
- a history of retinal disease. These people should seek advice from their ophthalmologist prior to operating the unit. In operating the curing light unit, this group of individuals must take extreme care and comply with all safety precautions (including the use of suitable light-filtering safety goggles).

WARRANTY PERIOD

Patterson Dental Supply, Inc. extends a five (5) year warranty against defects in materials or workmanship to the original purchaser of this product, except the battery which has a two (2) year warranty. The Manufacturer agrees to correct any defects which develop within the warranty period, either by repair or replacement, at its option. This warranty is valid providing factory inspection indicates that any such defect developed during normal and proper use subject to the conditions below.

WARRANTY CONDITIONS - PLEASE READ CAREFULLY

- Please fill out the warranty card included in kit and send promptly back to Patterson Dental Supply, Inc. with proof of purchase. Failure to do so may void your warranty.
- a) Claims for damage in shipment should be filed promptly with the transportation company.
- b) All shipments claimed defective can only be returned to the Manufacturer with the written consent of the Manufacturer. All returned products must be accompanied by a full description of discrepancy or malfunction.
- c) Only authorized service personnel are allowed to carry out repairs to the curing light. Manufacturer shall be released from all obligations under this warranty if repairs or modifications are made by persons other than Manufacturer's own or authorized service personnel.
- d) Use only authorized parts exclusively to replace defective components. The product warranty does not cover any damage resulting from the use of third-party replacement parts.
- e) **Manufacturer shall be released from all obligations under this warranty in the event of improper installation; damage due to casualty; or obvious misuse including but not limited to damage from liquids and autoclaving.**
- f) Manufacturer shall in no event be liable for any consequential damages arising from the misuse of the product.
- g) Only the above stated warranty shall be applicable, and all other warranties, expressed or implied, are hereby excluded.

Any enquiries concerning details of operation should be addressed to Patterson Dental Supply, Inc. Corporate Office, 1031 Mendota Heights Road, Saint Paul, MN 55120. Please specify the model and serial number.

TECHNICAL DATA

Handpiece (battery, control section and standard LED attachment)	
- wavelength range:	440 nm – 480 nm
- peak wavelength range:	460 nm
- RAMP mode elapsed time to full intensity:	5 seconds
- light intensity:	1500 mW/cm ² (peak)
- duration of continuous use NO RAMP mode:	60 seconds
- RAMP mode:	65 seconds
- total continuous run time with fully charged battery:	3 hours, 20 minutes
- dimensions:	22mm diameter x 243mm long
- weight:	176g / 6.2 ounces
- battery:	2 x 4.2V Lithium Ion – 1400 mAh

Plug Pack

- input voltage:	90 – 264V – 50/60 Hz
- output voltage:	12V DC
- output current:	750 mA
- dimensions:	70mm high x 45mm wide x 50mm long
- weight:	80g / 2.82 ounces

Charger

- time to charge battery:	60 – 180 minutes
- operating temperature:	10°C - 40°C
- total height with handpiece inserted:	255mm
- dimensions:	70mm wide x 40mm high x 130mm long
- weight:	84g / 3 ounces

Distributed by:
Distribué par:

**PATTERSON
DENTAL**
Corporate Office
1031 Mendota Heights Road
Saint Paul, MN 55120

**PATTERSON
DENTAL**
Montréal Québec
H3M 3E6

Transport and storage environment:
Ambient temperature: 0°C - 40°C (32°F - 104°F)
Relative humidity: 10% - 85%
Atmospheric pressure: 0.5 atm - 1.0 atm (500 hPa to 1060 hPa)

CAUTION
(refer to accompanying instructions)

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