High Performance Alginate Substitute

Patterson® High Performance Alginate Substitute is an addition cured elastomeric vinyl impression material specifically formulated to be used as an alternative to traditional alginate materials. Intended to be used primarily for taking impressions for case study models, temporary crown and bridge fabrication, partial denture impressions, Orthodontic models and counter impressions. This mint scented product offers excellent detail reproduction, quick setting time, dimensional accuracy and stability in both an easy to use 50 ml 1:1 ratio HP cartridge and 380ml 5:1 ratio cartridge. Combined, it's non-runny consistency, ideal tear strength and hardness makes **Patterson® High Performance Alginate Substitute** the logical choice when it comes to looking for a replacement to messy alginate impression materials.

Preparation for Use 50 ml Cartridges:

Important: Before attaching a mixing tip to the cartridge follow these steps:

- 1) Load the cartridge into the dispensing gun.
- 2) After the cartridge is properly positioned in the dispensing gun, remove and discard the shipping cap. Extrude a small amount of material to make certain that the product is exiting the barrels evenly and that there is no plugging of material.
- 3) After aligning the point on the mixer hub with the notch in the cartridges, push on the mixing tip an then rotate the hub 90° clockwise locking the mixer to the cartridge.
- 4) Squeeze the dispensing gun trigger to extrude material. As the material flows through the mixing tip it will be mixed. The mixing tip can be used directly in the mouth.
- 5) After using, do not remove the used mixing tip. Leave it in place on the cartridge where it acts as a seal for the remaining material in the cartridge. The next time the cartridge is used, remove and discard the old mixing tip, extrude a small amount of material and begin again with step #3.

Preparation for Use 380 ml Cartridges:

- 1) Position the cartridge nozzle side up with the nose cap tab facing you. While firmly holding the cartridge with both hands, using your thumb gently push on the tab until the nose cap unlocks. Using your thumb and forefinger slowly lift and remove the nose cap. Once free discard the nose cap, **DO NOT REUSE**.
- 2) Load the 380ml cartridge into the dispenser following the manufacturer's direction.
- 3) Following the operating instructions supplied by the machine manufacture, bleed the cartridge by extruding a small amount of material until the base and catalyst are evenly dispensing.

Wipe the extruded materials from the end of the cartridge using an upward movement to avoid cross contamination. This is recommended for first time use and when the cartridge has been stored upright for a prolong period of time.

- 4) Attach the dynamic mixing tip to the extrusion nozzle by aligning the holes of the cartridge with the corresponding mixer ports. Secure the tip to the cartridge using a bayonet locking ring (to lock rotate ring clockwise approximately a ¼ turn). Note: only use Sulzer/MixPac dynamic mixers as other manufacturer's tips are not compatible with this cartridge system.
- 5) Following the operating instructions supplied by the machine manufacture, start dispensing material discarding the first half inch of product. Immediately after dispensing the required amount of material, if the dispensing machine does not automatically release pressure from the pistons, manually move the plunger back relieving pressure on the pistons thus avoiding excess discharge of material. Do not remove the mixing tip when finished.
- 6) During the workday, the cartridge may be stored in the machine with the used mixer attached. To start a new impression, remove the locking ring (to unlock rotate ring counter clockwise approximately ¼ turn), then remove and discard the used mixer. Attach the new dynamic mixer starting at step #4.
- 7) At the end of the workday or when the cartridges will not be used for an extended period of time, it is recommend that the cartridge be removed from the machine and stored in an upright position.

Instruction for Use 50 ml Cartridges/380 ml Cartridges:

- 1) Prepare the teeth in the usual manner for impression taking.
- 2) In any technique where an impression tray is required, the first

step is to coat the tissue side of the pre-selected tray with vinyl adhesive. The optimum bond occurs when the adhesive film is applied as thin and even as possible. You can use either a perforated or non-perforated metal or plastic impression tray.

- 3) To load the impression tray follow the operating steps in the "preparation for use" section above. Fill the tray with the desired amount of material holding the tray in a tilted position keeping the tip immersed in the material to avoid air bubbles. Seat the filled impression tray into the mouth passively maintaining the position until set.
- 4) When fully set remove the impression in an axial direction without any rocking motion. The impression can now be gently rinsed and disinfected.

NOTES:

Patterson® High Performance Alginate Substitute is a non-toxic material that when set, is highly chemically resistant. Before the material is cured, spots on clothing can be removed using a commercial dry cleaning fluid.

Patterson® High Performance Alginate Substitute will not work properly and may never set when sulfur, eugenol or amines are present. Do not use with polyether, polysulfide or condensation curing silicone impression materials. Certain gingival retraction cords which contain sulfides can result in the sulcus areas remaining uncured.

Patterson® High Performance Alginate Substitute:

- Compatible with dental gypsum and epoxy materials. Impressions can be poured after 2 hours and up to 30 days later.
- Impressions can be disinfected with any standard disinfecting solution.

- Impressions can be silver or copper plated.
- Multiple pours are acceptable if care is taken.
- Shelf life is two years from date of manufacture when stored at room temperature (74°F or 23°C).

Alginate Substitute - 50ml - Fast Set
Work 1:00 Min | Work/Set 2:00 min

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Patterson Dental Supply, Inc.