

RIVA STAR™



RIVA STAR TREATMENT KIT
 Reorder code 8800503
 10 x Riva Step 1 capsules
 10 x Riva Step 2 capsules
 10 x Microbrush Silver
 10 x Microbrush Green



GINGIVAL BARRIER 2 SYRINGE REFILL
 Reorder code 7750006
 2 x 1g Gingival Barrier Syringes

riva star
 TOOTH DESENSITISING AND ANTI-CARIES AGENT
 AGENT DÉSENSIBILISANT ET ANTI-CARIE

SDI

KEEP REFRIGERATED
 GARDER REFRIGÉRÉ

Riva Star (Step 1 and Step 2)
 10 x Silver Capsules Step 1
 (with silver applicator brush)
 10 x Green Capsules Step 2
 (with green applicator brush)
 Riva Star (Étape 1 et Étape 2)
 10 x Capsules argentées - Étape 1
 (avec brosse d'application argentée)
 10 x Capsules vertes - Étape 2
 (avec brosse d'application verte)
 NPN 80085846

REFERENCES

1. Heo J, Tompkins GR, Love RM; Evaluation of the Antimicrobial Activity of Combined Silver Diamine Fluoride (SDF) and Potassium Iodide (KI) as an Endodontic Medicament; University of Otago, New Zealand, 2010.
2. Knight GM, McIntyre JM; The effect of silver fluoride and potassium iodide on the bond strength of auto cure glass ionomer cement to dentine. Australian Dental Journal 2006;51:42-45.
3. Knight GM; Silver Bullet. Australian Dental Journal 2010, Aesthetic update.

SDI YOUR SMILE.
 OUR VISION.



NEXT GENERATION SDF
SDF + KI
 SMART DESENSITIZING AND
 CARIES ARRESTING AGENT

CREATED BY SDI
RIVA STAR™



- Clinically tested & proven for decades
- Esthetically pleasing for the patient - reducing staining
- Approved for Caries Arrest
- Non-invasive - suitable for MID
- Patented Technology
- Immediate dentinal hypersensitivity relief in adults
- Increased bond strengths when used with auto cure GIC 1 such as Riva SC HV



COMPOSITES | GLASS IONOMERS
 TOOTH WHITENING | EQUIPMENT
 CEMENTS | ADHESIVES | RIVA STAR
 ETCHANTS | SEALANTS
 AMALGAMS | ACCESSORIES



YOUR SMILE.
 OUR VISION.

SDI (North America) Inc.
 1279 Hamilton Parkway
 Itasca, IL 60143
 Phone 1 800 228 5166
 Fax: 630-361-9222
 usa.canada@sdi.com.au
 www.sdi.com.au



RIVA STAR IS A NEW GENERATION CARRIES ARRESTING AGENT AND DESENSITIZER WHICH IMMEDIATELY REMOVES TOOTH SENSITIVITY

Riva Star - Actions

- Silver ●●● antimicrobial
- Iodide ●●● antimicrobial + reduces staining
- Fluoride ●●● bacteriostatic + antimicrobial effect
- High pH ●●● bacteriostatic + antimicrobial + encourages precipitation of minerals within the caries



A **STAR** product containing three powerful anticaries components:



Together, these have a powerful effect to provide instantaneous control of progressing caries and relieving hypersensitivity. Riva Star is the latest generation of a silver diamine fluoride (SDF) system that takes dentistry to another level.

Triple-innovators


Riva Star was designed by three renowned Australian clinicians, bringing together their clinical knowledge and experience of using these three exceptional components, and combining them into one STAR product - Riva Star!

THE ULTIMATE ANSWER TO PHARMACOLOGICAL CARRIES MANAGEMENT

Riva Star - Effective Caries Management

Riva Star - Antimicrobial activity¹

Riva Star is an effective biofilm inhibitor. This is supported by study results from Otago University (NZ), that Riva Star has a higher inhibition zone against four bacterial species compared to sodium hypochlorite.



Riva Star - Indications


- **Step 1 is for children 3 years to adults**
- **Step 2 is for adult use only**
- both step 1 and 2 are indicated for anti-caries / anti-cavity
- prevents, fights and /or protects against cavities or caries
- helps reverse the tooth decay process at the earliest stage before it can become a cavity
- helps reduce the sensitivity in cervical hypersensitive teeth in adults

Riva Star - Ideal for use with GICs

Conventional glass-ionomers bond very well to surfaces treated with Riva Star.²

Riva Star - A desensitizing star offering immediate relief

The silver fluoride and potassium iodide action in Riva Star blocks the microscopic tubules that make up dentine. A low-solubility precipitate is formed that gives instant relief. It acts by preventing the flow of fluid that produces sensitivity.²



Exposed dentinal tubules are well known for causing dentine hypersensitivity. Insoluble silver iodide precipitates blocking dentinal tubules.

Riva Star - Reduces the risks of endodontic therapy³

Pretreating non-symptomatic deep carious lesions with Riva Star can be the difference between retaining pulpal vitality or having the tooth proceed to endodontic intervention. Riva Star is a perfect solution to enhancing patient care and services of modern day practices.

Riva Star - Patented procedure

Unlike other silver fluoride systems, Riva Star's two step patented procedure minimises the risk of staining. By applying the potassium iodide solution over the silver fluoride, a silver iodide creamy white precipitate is formed and washed off.

Riva Star - Two years and beyond desensitizing effect

Clinical studies support the immediate effect of Riva Star and its long term effect of over two years.