Patterson® Flowable Composite

INSTRUCTION FOR USE

Patterson® Flowable Composite is the ideal versatile fluoride releasing, flowable, radiopaque, light cured material. The composite is directly injected into the cavity preparation for maximizing adaptation to the preparation. The shade range consists of thirteen popular shades - enamel A1, A2, A3, A3.5, B2, B3, C2, C3, D3, incisal, gingival, bleach; and OA2 dentin shade.

COMPOSITION:

35% wt (50% vol) multifunctional methacrylic ester
65% wt (50% vol) inorganic filler (0.4-10 micron)

INDICATIONS:

Class V restorations; tunnel / air abrasion preparations; pit and fissure sealants; minimal class I,II,III and IV restorations; radiopaque cavity liners; cervical abrasion/erosion lesions; abfraction lesions; implant seals; incisal abrasions; repair enamel defects; minor core build-ups; temporary crowns; cement porcelain/ceramic/ composite veneers; splint mobile teeth; attach fibre bridges; restorations in deciduous teeth; repair porcelain restorations; block out undercuts in inlay/onlay/ crown preparations; cover stains.

INSTRUCTIONS:

Clean and isolate tooth. A conservative cavity should be prepared, employing standard techniques and instruments, to form a slightly rounded internal form. Margins of the cavity preparation should end in sound and supported enamel with no bevels in stress bearing locations. If bevelling is desired in a low stress location it should extend no further than 1mm at an angle of no greater than 45°. Pre-wedging is also recommended to ensure that the restored tooth will have an adequate contact point. Prophyl all surfaces to be etched including surfaces adjacent to the cavity with an oil free non-fluoride containing paste or a slurry of pumice and water. Rinse thoroughly with water.

Shade Selection

Shade selection should be made with a just cleaned and moist tooth.

NOTE: Do not use uncured paste for color matching as there is a slight color change upon curing.

Isolation

Isolation techniques must be used to prevent contamination. Rubber dam is the preferred mode of isolation.

Pulp Protection

For deep cavities an appropriate liner or cement should be placed at the deepest point of the cavity.

1. Acid Etching

Thoroughly dry the surface to be etched with dry, oil-free air. Etch tooth surface with 37% phosphoric acid.

(a) Enamel only

Etch surface for at least 20 seconds.

(b) Dentin and enamel

Using the "total etch" technique etch the surface including any glass ionomer for at least 20 seconds.

Enamel subjected to fluoridation should be etched for 90 to 120 seconds.

Etching Precautions: Ensure that the dispensing tip hub is **firmly** attached to the syringe by twisting the hub securely onto the syringe. Avoid acid contact with oral tissues, eyes and skin. If accidental contact occurs wash thoroughly with water. In the case of eye contact, wash eye for 15 minutes and seek medical attention. Use matrix strips to protect adjacent tooth surfaces during etching. On first usage or after prolonged storage, extrude a small amount onto a mixing pad for familiarity with the etchant's viscosity and rate of extrusion.

2. Wash thoroughly with water.

3. Remove excess water. Keep moist. Avoid contamination e.g. saliva.

4. Apply bonding agent to saturate all internal surfaces according to manufacturer's instructions.

5. Directly inject the flowable composite in increments of 2mm or less in:

5.1 Class V restorations,

5.2 Tunnel preparations,

5.3 Pit and fissures sealants,

5.4 Minimal Class I, Class II, Class III and Class IV restorations.

Or other indication as required.

CAUTION: ensure that the dispensing tip is **firmly** attached to the syringe by twisting the tip hub securely onto the syringe.

WARNING: Apply restorative at normal room temperature (23°C/74°F).

8. Light cure the composite using a suitable visible light source for a minimum of 20 seconds in increments of 2mm.

9. Polishing of composite

(a) Remove excess composite and contour desired shape using a fine diamond or a 12-fluted carbide bur.

Tips for polishing using a flexible disc system:

(i) The polishing motion should be constant and unidirectional.

(ii) A back and forth movement over the composite-enamel margin is not recommended.

(iii) Keep the surface and polishing disc dry while polishing. A dry surface will produce a smoother, more uniform finish.

(iv) Do not use a handpiece with a speed greater than 35,000 rpm.

(v) Avoid touching the composite with the mandrel or disc eyelet.

(b) For gross reduction, set the speed of the handpiece to approximately 10,000 rpm. Use the coarsest grit disc of the range and polish the restoration starting at the gingival and moving outward over the restoration. Replace the disc if cutting efficiency impaired. Wash and dry.

(c) For final contouring, use a medium coarse grit at a speed of approximately 10,000 rpm. Wash and dry.

(d) For finishing, use a medium fine grit disc at a speed of approximately 30,000 rpm.

Wash and dry.

(e) Finally, use the finest grit disc at 30,000 rpm. Wash and dry.

Note: After step 9, polishing paste may be used to give the composite a lustrous finish, using the following steps.

(i) Apply polishing paste to a rubber cup or disc.

(ii) Smear a thin layer of paste onto the restoration.

(iii) Add a small amount of water to both the tooth and cup. Polish for 30 seconds at low speed and light pressure.

(iv) Wash and dry.

NOTES ON DIGITAL X-RAY UNITS: The glass filler in composites is what makes composites radiopaque. Flowable composites have less glass filler than standard composites, making them less radiopaque. Flowable composites are more radiopaque than dentin but less radiopaque than enamel. Typically, dentists that are new to digital x-ray are the dentists that have difficulty seeing the flowable composite. There are several factors to consider if the flowable composite cannot be seen:

- The position of the x-ray head may need to be adjusted.
- The radiation may need to be adjusted.
- The type of x-ray head may affect the result (Newer is better)
- The grey scale on the screen may need to be adjusted.

Doctors are also encouraged to contact their digital x-ray manufacture representative to help them understand all of the options and features to improve their radiographs. Doctors should be able to learn additional tips from their x-ray representative. This flowable composite meets the ISO standard (ISO 4049:2000(E).) for radiopacity.

PRECAUTIONS:

Avoid prolonged contact of bonding agents and composite with the skin or oral tissue, as it may cause inflammation of the oral tissues or skin sensitization.

Any persons having known resin allergies should immediately discontinue the use of bonding agents and composite.

Keep out of reach of children.

Do not take internally.

Do not refrigerate.

Store at temperatures between 10° and 25°C (50°-77°F).

When not in use replace the cap tightly.

Use at room temperature.

Do not use after expiry date.

Do not expose material to direct light.

Caution: Federal Law restricts this device to sale by on the order of a dentist.

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Patterson Dental Supply, Inc.