Application :

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PATTERSON®
DENTAL

USA

Patterson® Glass Ionomer Luting Cement

corresponds to DIN EN ISO 9917
Carefully read instructions prior to use

Product description:

Patterson® Glass Ionomer Luting Cement is a radiopaque glass ionomer cement for luting in application capsules for direct application. **Patterson® Glass Ionomer Luting Cement** is characterized by good rheological properties, a high level of compressive strength, fluoride release and a low film-thickness.

Each capsule contains a minimum of 0.19 ml of usable cement. Shades: universal white.

Indications:

- Luting of metal-based crowns, bridges, inlays and onlays
- Luting of high-strength all-ceramic crowns and bridges (such as aluminum oxide, lithium disilicate and zirconium oxide ceramics)
- Luting of metal-based and high-strength all-ceramic crowns and bridges (such as aluminum oxide, lithium disilicate and zirconium oxide ceramics) on implant abutments
- Luting of metal root posts or indirect metal post core build-ups
- Luting of prefabricated steel crowns
- Luting of orthodontic bands

Contraindications:

Patterson® Glass Ionomer Luting Cement contains polyacrylic acid. **Patterson® Glass Ionomer Luting Cement** should not be used in case of known hypersensitivities (allergies) to this ingredient.

Patterson® Glass Ionomer Luting Cement is not intended for non-reinforced all-ceramic restorations such as silicate ceramics, nor for indirect composite restorations, orthodontic brackets, ceramic inlays, ceramic onlays or veneers.

Method of Use:

Preparation:

Pre-treatment of the restoration according to the manufacturer's instructions (e.g. sandblasting of the inside of metal restorations). If necessary, degrease using alcohol or acetone. Mask the dentin near to the pulp with calcium hydroxide.

Clean the prepared cores (e.g. with a cleaning paste), rinse with water and dry. Do not over dry the dentin. The surface should appear slightly moist.

Activating the application capsule and mixing the material:

Only extract the application capsule from the tubular bag immediately before use. Shake the application capsule briefly by hand before use, in order to loosen the powder. During activation, the liquid from the plunger must be pressed into the inside of the capsule as fully as possible. To achieve this, the end of the application capsule where the plunger of liquid protrudes is pressed down onto a hard surface (e.g. a table) as far as the plunger will move (see fig. 1). Only activate the application capsule immediately before use. Load the activated application capsule immediately into a high-frequency mixer with approx. 4,000 oscillations per minute and allow it to mix for 10 seconds. Remove the application capsule from the mixer immediately and apply the material. Working time and consistency may be influenced by the respective type of capsule mixer.


Fig. 1: Activating the application capsule

Application:

Carefully open the outlet by lifting the application tip until you feel a stop (see fig. 2). Place the opened application capsule into the holding bracket of an application instrument. The cannula can be aligned as required by turning the whole application capsule in the holder. By operating the lever, dispense the material directly onto the surfaces which have been prepared for luting and distribute it using, for example, the application tip of the application capsule. Always apply a slight excess of material. Insert the restoration immediately applying light pressure and hold it in position continuing to apply light pressure while curing takes place. Avoid unnecessary waiting times between activation, mixing and application, since the incipient setting may make dispensing and the positioning of the restoration more difficult.

Fig. 2: Opening the application capsule

At room temperature, **Patterson® Glass Ionomer Luting Cement** allows a window of approx. 3 minutes working time, before setting begins after 5 to 7 minutes (3 - 5 mins. at 99 °F (37 °C)).

Higher temperatures shorten and lower temperatures extend these working times.

Prevent the ingress of water and saliva during application and setting phase.

Excess material can be removed particularly easily during the rubber-elastic phase (before complete curing) using a scaler or a curette, for example.

Notes and precautions:

Avoid direct application to dentin near the pulp or to exposed pulp. In case of contact with the eyes, rinse immediately with plenty of water and seek medical advice.

Instruments used need to be rinsed immediately with cold water before residues of material cure. Following the removal of orthodontic bands, residues of material can be removed quickly and gently using an ultrasonic scaler.

Storage:

Store **Patterson® Glass Ionomer Luting Cement** in a dry place at 59 °F to 73 °F (15 °C to 23 °C).

Only remove it from the foil immediately before use.

Do not use the product after its expiry date.

Product ordering information:

111-0626

This material has been developed solely for use in dentistry. Processing should be done strictly according to the instructions for use. Patterson® recognizes its responsibility to replace products if proven to be defective. Patterson® does not accept liability for any damage or loss, directly or indirectly, stemming from the use of or inability to use the products described. Before using, it is the responsibility of the user to determine the suitability of the product for its intended use. The user assumes all risk and liability in connection therewith. Descriptions and data constitute no warranty of attributes and are not binding.

CAUTION: U.S. Federal Laws restrict this device to sale by or on the order of a dentist.

No person is authorized to provide any information which deviates from the information provided in the instructions for use.

**Keep this material out of reach of children.
For dental use only.**

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