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PATTERSON®
DENTAL

Patterson® SE Adhesive Luting

Carefully read instructions prior to use

Patterson® SE Adhesive Luting is a universal, dual-curing, self-adhesive, resin cement. It is used for cementing indirect restorations made from ceramic, composite and metal as well as for cementing posts. Employing an additional adhesive or conditioner is not required when using **Patterson® SE Adhesive Luting**. **Patterson® SE Adhesive Luting** contains bi-functional methacrylate, acid methacrylate and inorganic fillers. The filler percentage is 70% w/w.

Patterson® SE Adhesive Luting can be used with any commercial coupling silane. Please observe the manufacturer's instructions.

The double-barrel syringe automatically mixes the material and prevents errors and bubbles. **Patterson® SE Adhesive Luting** comes with tapered cannulae for easily accessible cavities and restorations as well as blunt end cannulae. Attaching the intraoral tip to the blunt cannulae facilitates accessing difficult to reach areas and post drillings.

Indications:

- Definitive cementing of inlays, onlays, crowns and bridges (not Maryland bridges) made from ceramic, composite and/or metal
- Definitive cementing of metal, ceramic and fibre-reinforced posts

Contraindications: Patterson® SE Adhesive Luting contains Bis-GMA, aliphatic, aromatic and acid methacrylate as well as benzoyl peroxide, amines and BHT. Allergies to components cannot be excluded.

Do not use in patients with known sensitivities.

Pulp protection: Areas in close proximity to the pulp should be covered with a hard-setting calcium hydroxide preparation to avoid pulp irritations.

Method of use:

1. Removing the temporary and cavity preparation

Remove the temporary restoration and excess temporary cement by using appropriate methods. Clean the cavity or stump with a eugenol-free polishing paste or pumice paste. Any use of substances such as sealants, desensitizers, EDTA solutions e.g. after cleaning the cavity or stump is not recommended because their residues may have a detrimental effect on the polymerization reaction of **Patterson® SE Adhesive Luting**.

2. Preparation

Check the fit and contact points of the definitive restoration. If silicone has been used for verification, thoroughly remove silicone residue afterwards. Checking the occlusion on glass-ceramic restorations should only be carried out after definitive cementing to reduce the risk of fracture. Remove any residue from the surface of the tooth and restoration after completing the try-in. Prepare the restoration according to the manufacturer's instructions and carefully dry. Dry the surface of the tooth preferably with either a foam pellet or gentle puff of air, leaving the dentin slightly moist.

3. Working and setting times

Working time: 90 s (73°F (23°C))

Intra-oral working time after insertion: 30 - 60 s

Setting time: 4 min (99°F (37°C))

These times may vary based on temperature, humidity, storage conditions and age of product.

4. Applying and seating of the restoration

Due to technical reasons slight differences in the filling level between the syringe chambers can occur prior to the first use. Therefore remove syringe cap and squeeze the material out until the material is delivered evenly from both chambers. Attach a mixing tip. The notches on the syringe and the mixing tip must line up exactly. Turn the mixing tip 90° clockwise to lock the tip into place. To guarantee homogeneous mixing of both pastes, please discard the first 2 mm of the extracted material. After that, apply the material immediately and directly onto the prepared contact surfaces. An intraoral tip can also be attached to the tip of the mixing tip with blunt end, if necessary. Seat restoration gently onto tooth prep, allowing the cement to flow from all margins in order to leave slight excess amount. Either replace transport cap or use the attached mixing tip as a cap.

5. Excess removal

Remove excess material from difficult to reach areas first. Clean contact points and interproximal areas with dental floss. Excess cement can be removed with a suitable instrument (scaler) either in the gel phase (ca. 60 - 90 s after insertion) or after brief polymerization (1 - 2 s). Large-volume excess is easier to remove. The restoration must be held in position during the polymerization phase.

As an alternative, the coarse excess of the resin cement can be removed with a foam pellet or a similar object immediately after insertion. Leave a slight excess and cover it with a glycerine gel (air-block) or light-cure it.

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6. Curing

Polymerization time can be controlled with additional light-curing. Cure translucent restorations with light for 20 s on each side. Light-curing the cementing gap of opaque restorations for 10 s on each side is recommended. The auto-polymerization of **Patterson® SE Adhesive Luting** is completed 4 min after insertion.

7. Finishing

Excess can be processed immediately after polymerization. Remove cured excess material with a finishing diamond or flexible polishing disk. Finish approximal areas with finishing/polishing strips.

Cementing endodontic posts:

1. Prepare root canal for treatment. Prepare the endodontic post according to the manufacturer's instructions.
2. Disinfect the root canal, rinse with plenty of water and avoid hydrogen peroxide.
3. Remove the transport cap from the double-barrel syringe. Subsequently attach the blunt mixing tip and turn to adjust its position. Attach the intraoral tip to the mixing tip.
4. Dry the root canal and dentin with paper points, leaving the dentin slightly moist.
5. Applying **Patterson® SE Adhesive Luting** directly into the root canal is recommended (begin applying from apical, allowing the intraoral tip to be immersed in the resin cement to avoid air bubbles during the application). Do not use a lentulo! The endodontic post can also be wetted with **Patterson® SE Adhesive Luting**. Insert the endodontic post with a turning motion with the aim of leaving a slight excess of cement. Remove the rough excess and light-cure **Patterson® SE Adhesive Luting** in the direction of the endodontic post to anchor it.

Warnings and precautionary measures:

In case of contact with the eyes, immediately rinse with water and consult a doctor.

Do not use phenolic substances (e.g. substances containing eugenol and temporary cements), as they can disturb polymerization of Patterson® SE Adhesive Luting. Only use eugenol-free temporary cements for temporary luting.

For Maryland bridges please use conventional luting systems in combination with total etch or self-etch bondings.

Uncured methacrylate resin may cause contact dermatitis and damage the pulp. Avoid contact with skin, eyes, and soft tissue. Wash thoroughly with water after contact. In case of contact with skin or eyes, flush immediately with water. Consult a doctor for eyes. Do not take internally.

The warranted product characteristics can only be achieved when using the mixing tips as approved by Patterson®. All the mixing tips for use with our respective products are available from Patterson®.

Storage:

Store at temperatures between 36°F - 46°F (2°C - 8°C). Do not use after the expiration date.

Allow **Patterson® SE Adhesive Luting** to reach room temperature before use.

Product ordering information:

111-0634

This material has been developed solely for use in dentistry. Processing should be done strictly according to the instructions for use. Patterson® recognizes its responsibility to replace products if proven to be defective. Patterson® does not accept liability for any damage or loss, directly or indirectly, stemming from the use of or inability to use the products described. Before using, it is the responsibility of the user to determine the suitability of the product for its intended use. The user assumes all risk and liability in connection therewith. Descriptions and data constitute no warranty of attributes and are not binding.

CAUTION: U.S. Federal Laws restrict this device to sale by or on the order of a dentist.

No person is authorized to provide any information which deviates from the information provided in the instructions for use.

Keep this material out of reach of children.

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