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Patterson® Temporary Luting Cement

corresponds to DIN EN ISO 3107 type 1

Carefully read instructions prior to use

Patterson® Temporary Luting Cement is a eugenol-free cement with calcium hydroxide, is easy to process and does not influence the setting of resin-based filling and luting materials.

Patterson® Temporary Luting Cement is recommended for patients with known allergies to eugenol.

Patterson® Temporary Luting Cement contains calcium hydroxide to support tooth vitality.

Indications:

Patterson® Temporary Luting Cement is used for cementing temporary crowns and bridges, inlays and onlays.

Contraindications:

Patterson® Temporary Luting Cement contains zinc oxide, magnesium oxide, calcium hydroxide, fatty acids, modified rosin and vegetable oils. Do not use in cases of known sensitivities or allergies.

Method of use:

Patterson® Temporary Luting Cement

Due to technical reasons slight differences in the filling level between the syringe chambers can occur prior to the first use. Therefore remove syringe cap and squeeze the material out until the material is delivered evenly from both chambers. Place and lock the mixing tip with a 90° turn clockwise. The material is automatically mixed in the mixing tip.

All double-barrel syringes come with a detachable piston. The piston can simply be pulled from the double-barrel syringe and replaced without causing any material leakage. The plugs in the syringe chambers remain in place, preventing cross-contamination of the materials.

Application:

Apply a thin layer of **Patterson® Temporary Luting Cement** inside the temporary crown or bridge and affix it to the prepared tooth.

Patterson® Temporary Luting Cement – in tubes remains workable for about 1.5 - 2 min, **Patterson® Temporary Luting Cement** – in cartridges/ in double-barrel syringes about 1 min. All materials set in the mouth within 4 - 6 min.

Place the temporary into the mouth. Remove any material excess after 3 - 5 min.

Precautionary measures:

Pulp protection (e.g. calcium hydroxide cement) is recommended in cases of exposed pulp or inadequate dentin thickness.

The sensitivity of cores can be reduced by using desensitizing preparations prior to temporary luting. Remove these materials completely before permanent luting (Do not use preparations containing eugenol!).

Extensive superficial and multi-surfaced cavities as well as temporary long-term sealing work should be treated with a harder cement.

The warranted product characteristics can only be achieved when using the original mixing tips as approved by Patterson®. All the mixing tips for use with our respective products are available from Patterson®.

Supplementary instructions/information:

Moisture (humid cores) further facilitates removal of the temporary.

Before removing **Patterson® Temporary Luting Cement** residues, spray with a cooling spray and wait for a short time.

Store **Patterson® Temporary Luting Cement** at temperatures between 39 - 73°F (4 - 23°C). Close tubes tightly after use.

Patterson® Temporary Luting Cement corresponds to EN 23107 and ISO 3107/1988.

Product ordering information

111-0618

This material has been developed solely for use in dentistry. Processing should be done strictly according to the instructions for use. Patterson® recognizes its responsibility to replace products if proven to be defective. Patterson® does not accept liability for any damage or loss, directly or indirectly, stemming from the use of or inability to use the products described. Before using, it is the responsibility of the user to determine the suitability of the product for its intended use. The user assumes all risk and liability in connection therewith. Descriptions and data constitute no warranty of attributes and are not binding.

CAUTION: U.S. Federal Laws restrict this device to sale by or on the order of a dentist.

No person is authorized to provide any information which deviates from the information provided in the instructions for use.

Keep this material out of reach of children.**For dental use only.****Ciment pour scellements provisoires Patterson®**

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