Patterson® Bite Registration
Vinyl Polysiloxane (VPS) Bite Registration Material DIN 13903
Vanilla Fast, Vanilla Super Fast, Unflavored Fast,
Unflavored Super Fast, Rigid Unflavored Fast

Manufactured for:
 Patterson Dental Supply, Inc.
 1031 Mendota Heights Road
 Saint Paul, MN 55120

37362/0819

Product Descriptions / Indication

Patterson® bite registration materials (Vanilla Fast, Vanilla Super Fast, Unflavored Fast, Unflavored Super Fast, and Rigid Unflavored Fast), are syringeable addition-curing elastomeric bite registration materials intended for making accurate occlusal records.

- Patterson® bite registration materials (Vanilla Fast, Vanilla Super Fast, Unflavored Fast, Unflavored Super Fast, and Rigid Unflavored Fast), can be used for making accurate occlusal registrations.
- Patterson® bite registration materials (Vanilla Fast, Vanilla Super Fast, Unflavored Fast, Unflavored Super Fast, and Rigid Unflavored Fast), can be used for applications, which require a hard bite registration silicone.
- Patterson® bite registration materials (Vanilla Fast, Vanilla Super Fast, Unflavored Fast, Unflavored Super Fast, and Rigid Unflavored Fast), Fast can be used for "small applications" because of their shorter working time.

Warnings

Do not use Patterson® bite registration materials for conventional impressions as their high end hardness values severely inhibit withdrawal from the mouth and may result in inadvertent extraction of teeth.

Cautions

Do not use after expiration date.
 Do not leave any residual bite registration material in the sulcus or oral cavity.
 Do not swallow bite registration material! If swallowed: consult a medical doctor if problems arise or persist.
 Avoid contact with the eyes. If there is accidental contact with the eyes: rinse immediately and thoroughly with an eyewash or water. Consult a medical doctor if problems arise or persist.
 The product does not normally cause allergic reactions; however, for sensitive persons, an allergic reaction to the product cannot be ruled out.

Instructions for use

Use of products containing active sulfur, aluminum chloride or nitrogen compounds (retraction cords containing ferric sulfate, polysulfide impression materials etc.) in conjunction with this product will interfere with the setting reaction of Patterson® bite registration material. Use of these materials requires the area to be rinsed thoroughly to remove all residues.
 If gloves are worn when handling the material, they should be tested with a sample of the bite registration material beforehand to ensure compatibility. Do not use latex gloves.
 Do not allow the material to enter the sewer or water system to avoid environmental hazard.
 Avoid contact with clothing, as the material cannot be removed by dry cleaning.

Note

The total working time and intraoral setting time should be strictly adhered to in order to ensure optimal bite registrations. Should the temperature of the material deviate noticeably from 23 °C (74 °F), the total working time and setting times will be affected.

Instructions for use

Insert the cartridge into a dispensing gun DS-50 1:1/2:1.
 Ensure that the notches on the cartridge base are pointing downwards. The clasp will not close if the cartridge is not properly inserted into the dispensing gun.
 Remove the cartridge cap. The cap can be replaced after initial use.
 Advance plunger into cartridge and dispense a small amount of Patterson® bite registration material until equal amounts are extruded at the same rate.
 Install a mixing tip on the cartridge and turn cap 90 degrees in a clockwise direction to lock in place.
 Inject material directly onto the patient's occlusal surfaces or occlusal rims.
 Replace the plug after use. Before using the cartridge again, remove the cartridge cap or the filled mixing tip and check the cartridge outlet openings to ensure that no polymerized materials have sealed the openings.
 Clear any clogged openings and dispense a small amount of impression material until equal amounts are extruded at the same rate.
 Install a new mixing tip and continue as described above.

Cleaning and disinfection of the Dispensing Gun:

- Remove the cartridge.
- Remove any uncured paste using an alcohol wipe.
- The dispensing gun can be used several times if properly disinfected. Patterson® recommends exchanging the gun if there are visible traces of use. Disinfection can be carried out using a commercially available dipping disinfectant. Patterson® recommends the use of glutaraldehyde-based solutions.

Storage Conditions

Store in a dry place at room temperature; do not expose to sunlight and avoid exposure to extreme temperature fluctuations.

Disinfection

Bite registrations can be disinfected using for example, a 2% acidic glutaraldehyde solution. Use disinfection solutions. See manufacturer instructions.
 For professional use only. Caution: Federal (U.S.A.) Law restricts this device to sale by or on the order of a dentist.

Limitation of liability

Except where prohibited by law, Patterson® will not be liable for any loss or damage arising from this product, whether direct, indirect, special, incidental or consequential, regardless of the theory asserted, including warranty, contract, negligence, or strict liability.

Product Overview and Technical Data *Total setting time (removal time from mouth) from start of mix

Product name	Working time at 23° C / 74° F	Intraoral setting time at 35° C/95°F ≥	Total setting time* ≥	Hardness (approx. Shore)	Linear dimensional change (approx.) %	Mixing technique	Mixing element
Vanilla Fast	30 seconds	1 minute 30 seconds	2 minutes	A 90/HD 58	-0.2	Dispensing gun DS-50 1:1/2:1	Mixing tip, green, Ø 6.5 mm
Vanilla Super Fast	15 seconds	45 seconds	1 minute	A 90/HD 58	-0.2	Dispensing gun DS-50 1:1/2:1	Mixing tip, green, Ø 6.5 mm
Fast Unflavored	30 seconds	1 minute 30 seconds	2 minutes	A 90/HD 58	-0.2	Dispensing gun DS-50 1:1/2:1	Mixing tip, green, Ø 6.5 mm
Super Fast Unflavored	15 seconds	45 seconds	1 minute	A 90/HD 58	-0.2	Dispensing gun DS-50 1:1/2:1	Mixing tip, green, Ø 6.5 mm
Fast Rigid	30 seconds	1 minute 30 seconds	2 minutes	D 43/HD 76	-0.1	Dispensing gun DS-50 1:1/2:1	Mixing tip, green, Ø 6.5 mm

Legend	
Temperature limit	Use-by date
Keep away from sunlight	Keep dry
Batch code	Consult instructions for use
Dual system	Diameter
Cartridge	Static mixing tip MB (do not re-use)
dispensing gun DS 50 1:1/2:1	Mixing tip Type B MB
Millimeter	Milliliter
Less than or equal to	Greater than or equal to
Rx only For professional use only. Caution: Federal (USA) Law restricts this device to sale by or on the order of a dentist.	

The information provided for Patterson® products is based on comprehensive research and experience in application technology. Results are furnished to the best of our knowledge, subject to technical changes within the framework of product development. However, users must comply with and consider all recommendations and information in connection with any use. Updated: Feb 11th, 2019

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