GOCCLES®

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PRELIMINARY NOTES

- Read all instructions before using the product
- Understand how a healthy oral cavity looks like through Goccles in order to properly appreciate what may be abnormal
- Collect the patient's relevant medical and dental history before proceeding with the screening.


EXAMINATION

1. Conduct an extra-oral and intra-oral examination to the naked eye, also palpating all the structures of the head and neck
2. Repeat the intra-oral examination using Goccles eyewear in combination with a dental curing light
 - Reduce the room light for better contrast
 - The device has to be worn in adherence to the face, avoiding external light to reach the operator's eyes. Use the provided neck cord to optimize the fitting
 - Characteristics of the curing light: emission spectrum between 380 and 500 nm, emission peak between 440 and 480 nm. **IMPORTANT** Provide the patient a proper device for eye protection:

foto 1 - picture 1 - fotografia 1

foto 2 - picture 2 - fotografia 2

orange filters with UV525 protection are the best solution.

- While performing the examination, maintain the curing light at a distance of approximately 20-50 cm from the oral cavity to optimize the visualization of the natural tissue fluorescence
3. Abnormal tissue typically should appear as an irregular, dark area that stands out against the green fluorescence pattern of surrounding healthy tissues (picture 1)
 4. If a suspicious area is discovered, re-evaluate it trying to identify any causes for the region to appear abnormal. Take into consideration its consistency, its appearance both with naked eye and through Goccles and any relevant patient history information
 5. Record all significant findings and inform the patient about the appropriate course of action

FOLLOW UP

- If lesions or alterations of the oral mucosa are found, it's crucial that the clinician examines the patient after a couple of weeks
- Proceed with a new investigation of the oral cavity with Goccles to assess the presence of the abnormal area, after that, evaluate any change in the area, also checking if the presumed causative agent has been removed
- If the abnormal area has not disappeared, proceed the examination recommending the patient to perform a specialistic visit for further investigations

Characteristics of suspicious lesions and oral cancer

- Leukoplakia: a white patch that cannot be wiped off with gauze and with no obvious explanation
- Erythroplakia: a fiery red patch that cannot be characterized either clinically or pathologically
- Erythroleukoplakia: an abnormal patch of red and white tissue on mucous membranes in the mouth
- Evident loss of fluorescence
- The lateral tongue and floor of mouth are the regions at greatest risk for cancer development
- Asymmetry and/or irregular shape
- The lesion doesn't spontaneously resolve and generally increases in size over time

Other phenomenon evident through fluorescence visualization

- The presence of dark pigments (such as liquorice or coffee) usually appears as dark spots both with naked eye inspection and autofluorescence examination with Goccles (picture 2)
- Inflammation typically appears as darker areas characterized by poorly-defined borders due to the excess of blood content. If the normal fluorescence returns after a light pressure on the region, then the lesion may have an inflammatory origin
- Hyperkeratosis may often appear bright through Goccles because of strong keratin fluorescence. Goccles is not suitable for the study of hyperkeratotic lesions, however it can help to identify their margins which appear as a dark halo surrounding the main lesion

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

Mucosa of a healthy individual appears as the examination of auto-fluorescence; a: b: vestibular gingival; c: lingual belly; d: dorsum of the tongue; e: buccal mucosa; f: hard palate; g: floor of the mouth.

Loss of fluorescence in patients with mild dysplasia; a: patient with lichen oral planus; b: the same patient with loss of fluorescence of buccal mucosa (arrows); patient who has already undergone previous biopsies for leukoplakia, left buccal mucosa; d: halo intralesional loss of fluorescence within the apparently normal mucosa (arrows)

Hyper-fluorescence in a patient with hyperkeratotic lesion (base of tongue and floor of the mouth)

PIERREL Pharma srl
Strada Statale Appia 46/48 81043
Capua (CE) Italy

goccles@pierrelgroup.com

www.goccles.com