

VITA Easyshade[®] Compact

VITA Easyshade[®] Compact



Exclusive
Distributor


PATTERSON
DENTAL

VITA[®]

Total Shade Control in the Palm of Your Hand.


The all-new VITA Easyshade Compact delivers the fastest, easiest way to measure tooth shades under any lighting condition. Easyshade Compact gives you the ability to conquer the trickiest shade challenges all day, every day, instantly and accurately. Now you can measure natural teeth, check restoration shades and track bleaching progress with the push of a button.



VITA Easyshade® Compact

Cordless, Sleek and Incredibly Fast.


The pressure of getting a perfect shade match has never been higher due to patients' increased attention to esthetic appearance. The VITA Easyshade Compact quickly and precisely measures 16 VITA Classical Shades, 29 VITA System 3D-Master® shades and even those hard-to-see intermediate shades. High-esthetics with VITA materials has never been easier to achieve.



VITA Easyshade® Compact

Powerful Plug-and-Play Features.

The amazing design of Easyshade Compact hints at the powerful technology working underneath its sleek shell. Spectrophotometer technology, the latest OLED display and white-LED color-corrected lighting just for starters.



OLED display with easy-to-read icons

Versatility and Flexibility.

Measure a single central, check crown accuracy before seating, get a more detailed shade analysis with a detailed map – it's all here with Easyshade Compact. And with built-in memory that stores the last 25 measurements, you never have to worry about forgetting the result.


Select the appropriate measurement setting easily using the navigation buttons and easy-to-follow icons:

- Tooth single mode
- Tooth areas mode: maps the cervical, middle and incisal shades
- Verify restoration mode: check a restoration for shade accuracy
- Shade tab mode: practice proper placement of probe on VITA shade tab




Quality Control at the Push of a Button.

For the dental laboratory, Easyshade Compact verifies shade accuracy during the buildup process to determine where adjustment may be required.



By displaying which color parameters need attention (value, chroma or hue), Easyshade Compact makes it a simple matter to know precisely which VITA 3D-Master or VITA Classical shades should be used to make shade corrections.



VITA Classical Shade Guide



VITA 3D-Master® Toothguide



VITA Bleachedguide®

VITA Easyshade® Compact

Welcome to the Evolution.

Today, millions of satisfied patients around the world enjoy beautiful smiles from their dentists, aided by the superior esthetics of VITA shades and materials. Explore the entire system of VITA materials, including CAD/CAM blocks, all-ceramic and metal porcelains, denture teeth and composites.



VITABLOCS® Guide 3D-Master



VITA Linearguide™

VITA Easyshade® Compact



VITA Easyshade® Compact

After Shade Measurement,
Prescribe VITA Materials for
Superior Esthetic Results!

**VITA SHADE,
VITA MADE.**

Ordering Information

Vident Part Number	Product Description
E75-201	VITA Easyshade® Compact
Patterson Part Number	
706-3605	

To order call Patterson Dental at
800-873-7683

Exclusive
Distributor


PATTERSON
DENTAL

www.vident.com • 800-828-3839 USA • 800-263-4778 CAN

VITA®

© Vident 2011. All rights reserved. L-9404V Rev.0111. Easyshade, 3D-Master, VITABLOCKS, Bleachedguide, and Linearguide are trademarks or registered trademarks of VITA Zahnfabrik.