



A PERFECT PAIR

Quality glove options
for the task at hand



FOR THOSE WHO ENSURE HEALTH AND CREATE BEAUTIFUL SMILES

Our proven, high quality, single-use infection prevention and control solutions for PPE, sterilization, office essentials, and preventive care are thoughtfully designed for simplicity and comfort without compromise to safety and clinical efficacy.

Making each day in your practice a little easier.

Medicom was founded in 1988 in response to the urgent need for medical gloves for healthcare professionals during the global HIV crisis.

Understanding the critical role that quality gloves play in safety and protection against cross-contamination, we have continued to develop a diverse portfolio of premium medical examination and surgical glove options designed to meet the individual clinical needs and preferences.

Backed by over 35 years of expertise, we invite you to begin your journey here to finding your perfect pair.



CHOOSING THE RIGHT GLOVES

Protection and comfort

When it comes to hand protection, it is important to choose the right type of glove for each specific task or application. Our broad selection of powder-free medical examination and surgical gloves is innovatively designed to ensure superior protection and the comfort you deserve.

Our examination gloves feature cutting-edge technologies and special formulations to meet the evolving needs of dental professionals and their patients.

Our surgical offering includes synthetic alternatives that offer the benefits of latex gloves without the risks to patients and staff who suffer from Type I (latex) allergies.



Which material best fits your needs?

Each glove material has unique strengths and limitations. The following descriptions are representative of basic material characteristics.

Glove Attribute	Nitrile (Exam)	Polychloroprene (Exam/Surgical)	Polyisoprene (Surgical)	Latex (Exam/Surgical)	Vinyl (Exam)
Level of Barrier Protection	Excellent	Excellent	Excellent	Excellent	Poor
Strength and Durability	Excellent	Excellent	Very Good	Excellent	Poor
Puncture Resistance	Excellent	Good	Good	Very Good	Poor
Elasticity	Very Good	Excellent	Excellent	Excellent	Fair to Poor
Comfort and Fit	Very Good	Very Good	Excellent	Excellent	Fair
Allergy Potential	Low	Low	Low	High	Very Low

Individual glove performance may vary.

FIT AND FEEL

Grip

Gloves with textured fingertips ensure a secure grip, even when wet. Fully textured gloves provide ultimate grip.

Weight/Thickness


Glove weight and thickness should match the level of protection and tactile sensitivity required for specific tasks.

Size

The right size of glove will ensure a comfortable fit that is not too snug or too loose. Examination gloves are ambidextrous for added convenience. Surgical gloves are hand-specific (right and left) and made in more precise sizing.

Cuffs

Beaded cuffs on gloves prevent roll-down and protect against fluid drips. Extended cuffs offer additional coverage for wrists and forearms.



CHEMO SAFE

Tested and approved for use with chemotherapy drugs.



FENTANYL SAFE

Tested and approved for use with fentanyl.



CHEMICAL RESISTANT

Tested and approved for use with certain chemicals.



HAND HEALTH

Infused with moisturizing agents to soothe and rejuvenate dry hands.



NO ACCELERATORS ADDED

Helps reduce the risk of Type IV allergies caused by chemical accelerators commonly used in the manufacturing process.



ANTIMICROBIAL TECHNOLOGY

Helps prevent cross-contamination by eliminating microbes within five minutes of contact with the exterior of the glove.

Consider the intended use of your gloves

Medical examination gloves are ideal for patient examinations, blood tests, and patient care. Surgical gloves must be used for sterile procedures.

For handling drugs, chemical components or food, you must use gloves that comply with industry standards for these tasks.

The importance of fit

Wearing properly fitting gloves is crucial to working safely and efficiently. Gloves that are too tight can tear easily, contribute to hand fatigue and cause injuries like carpal tunnel syndrome. Gloves that are too loose decrease tactile sensitivity and manual dexterity and increase the risk of cross-contamination.

NON-LATEX STERILE SURGICAL GLOVES

Offering enhanced comfort, safety and performance for surgical procedures

Our synthetic polychloroprene and polyisoprene gloves deliver the performance of natural rubber latex without the risk of Type I allergies or latex sensitization.

- Produced on ergonomically-designed formers for enhanced comfort and fit
- Polymer coating for easy damp and dry-hand donning
- Available in several thickness options to suit your needs

Our non-latex sterile surgical line up can be used in a variety of procedures.*

	General surgery	Cardiovascular	Dental/Maxillofacial	Ear, Nose and Throat	Endovascular	Labour and Delivery	Laparoscopic/Robotics	Neuro	Obstetrics	Ophthalmology	Orthopedics	Pediatrics	Plastics	Thoracic	Urology	Vascular
1144 AssureTouch™ Guard Micro Polychloroprene		•	•	•	•		•	•		•		•	•			•
1145 AssureTouch™ Sense Polychloroprene	•	•		•	•	•	•	•	•	•	•	•	•	•	•	•
1146 AssureTouch™ Defend Polychloroprene	•	•		•	•	•	•	•	•	•	•	•	•	•	•	•
1147 AssureTouch™ React Polyisoprene	•	•		•	•	•	•	•	•	•	•	•	•	•	•	•

*Reference based on publicly available data and information of common practice. Not intended as recommendation of use.

Need a higher level of protection?

Wearing two pairs of surgical gloves, or "double gloving," can help reduce the risk of cross-contamination in the operating room. The underglove provides an extra layer of protection, especially if the outer glove is punctured. Choosing undergloves in contrasting colors also makes it easier to detect perforations.

Do you have a latex allergy/sensitivity, chemical allergy or contact dermatitis?

Choosing synthetic gloves can help avoid potentially severe allergic reactions to latex. Those experiencing adverse reactions to synthetic gloves may benefit from switching to gloves not made with chemical accelerators. Try Medicom nitrile medical examination gloves with no accelerators added.

GET TO KNOW YOUR GLOVES

Non-Sterile Medical Examination Gloves



POLYCHLOROPRENE Ocean Pacific® NeoNatural

Get the best of both worlds

Combines the softness of latex with the strength of nitrile for reliable protection without sacrificing comfort or dexterity.



Item	1228	
Colour	Emerald ●	
Weight	6 g	
Thickness: mm (mil)	Finger 0.11 (4.33)	Palm 0.08 (3.15)
Size	XS, S, M, L, XL	
Packaging	100 gloves/box 10 boxes/case	
Grip	Textured fingertips	



NITRILE Ocean Pacific® Eclipse

Precise fit for precise tasks

Thin and stretchy glove that conforms to hands to help minimize fatigue and maximize tactility.


Item	OP-EC	
Colour	Blue ●	
Weight	3 g	
Thickness: mm (mil)	Finger 0.05 (1.97)	Palm 0.05 (1.97)
Size	XS, S, M, L, XL	
Packaging	200 gloves/box 10 boxes/case	
Grip	Textured fingertips	



NITRILE Ocean Pacific® Lab Master

We get your sensitive side

Manufactured without chemical accelerators known to cause skin reactions, which makes them safe for professionals with Type IV chemical allergies or sensitivities.



Item	OP-LM2	
Colour	Indigo ●	
Weight	3.2 g	
Thickness: mm (mil)	Finger 0.08 (3.15)	Palm 0.05 (1.97)
Size	XS, S, M, L, XL	
Packaging	200 gloves/box 10 boxes/case	
Grip	Textured fingertips	



NITRILE Ocean Pacific® TruAloe



Tough protection that puts a smile on your hands

Aloe vera-infused to soothe and moisturize hands during use. Manufactured without chemical accelerators known to cause skin reactions.




Item	1226	
Colour	Aloe Green ●	
Weight	3.5 g	
Thickness: mm (mil)	Finger 0.1 (3.94)	Palm 0.06 (2.36)
Size	XS, S, M, L, XL	
Packaging	200 gloves/box 10 boxes/case	
Grip	Textured fingertips	



NITRILE Ocean Pacific® Elements

81% biodegradable in less than 1.5 years

Specially formulated for accelerated biodegradation in landfills, without compromising product quality. Tested as per ASTM D5511.



Item	OP-BIO	
Colour	Blue ●	
Weight	3.5 g	
Thickness: mm (mil)	Finger 0.07 (2.76)	Palm 0.05 (1.97)
Size	XS, S, M, L, XL	
Packaging	200 gloves/box 10 boxes/case	
Grip	Textured fingertips	



NITRILE Ocean Pacific® Blue Lite

Ideal balance of comfort and strength

A robust glove providing enhanced dexterity and tactile sensitivity, coupled with an excellent fit at palms and thumbs.



Item	OP-BL2	
Colour	Blue ●	
Weight	3.5 g	
Thickness: mm (mil)	Finger 0.08 (3.15)	Palm 0.06 (2.36)
Size	XS, S, M, L, XL	
Packaging	200 gloves/box 10 boxes/case	
Grip	Textured fingertips	



NITRILE Ocean Pacific® Intuition

Dexterity, softness and style

A strong glove offering excellent dexterity and tactility in a beautiful shade of lilac. Soft as silk.



Item	OP-IN	
Colour	Lilac ●	
Weight	3.5 g	
Thickness: mm (mil)	Finger	Palm
	0.08 (3.15)	0.06 (2.36)
Size	XS, S, M, L, XL	
Packaging	200 gloves/box 10 boxes/case	
Grip	Textured fingertips	



NITRILE Ocean Pacific® Ninja

Superior puncture resistance and grip

Optimal protection and ultimate grip. Approved by Health Canada for food handling.



Item	OP-NI	
Colour	Black ●	
Weight	5.3 g	
Thickness: mm (mil)	Finger	Palm
	0.12 (4.72)	0.1 (3.94)
Size	S, M, L, XL	
Packaging	200 gloves/box 10 boxes/case	
Grip	Fully textured	



NITRILE AssureTouch™ Basics

Value without compromise

Affordable protection coupled with excellent manual dexterity and tactile sensitivity.



Item	9985	
Colour	Blue ●	
Weight	3.5 g	
Thickness: mm (mil)	Finger	Palm
	0.07 (2.76)	0.05 (1.97)
Size	XS, S, M, L, XL	
Packaging	100 gloves/box 10 boxes/case	
Grip	Textured fingertips	



NITRILE AssureTouch™ Noir NEW

Safety meets sophistication

The strength and durability of nitrile in a stylish, light-weight glove. Provides exceptional protection and tactility while minimizing the visibility of blood.



NITRILE SafeBasics® True Fit Thin

Versatile all-purpose glove

The unrivaled strength and durability of nitrile for reliable protection and exceptional fit.



Item	9955	
Colour	Black ●	
Weight	3.5 g	
Thickness: mm (mil)	Finger	Palm
	0.11 (4.3)	0.06 (2.36)
Size	XS, S, M, L, XL	
Packaging	200 gloves/box 10 boxes/case	
Grip	Textured fingertips	

Item	1185		1186	
Colour	Blue 		Pink 	
Weight	3.5 g		3.2 g	
Thickness: mm (mil)	Finger	Palm	Finger	Palm
	0.08 (3.15)	0.05 (1.97)	0.07 (2.76)	0.05 (1.97)
Size	XS, S, M, L, XL		XS, S, M, L	
Packaging	300 gloves/box, 8 boxes/case			
Grip	Textured fingertips			



NITRILE SafeTouch® Advanced Comfort

Ultimate comfort and tactility

Offers excellent tactile sensitivity and dexterity while remaining highly resistant to punctures and tears.

Item	1196	
Colour	White ○	
Weight	2.5 g	
Thickness: mm (mil)	Finger	Palm
	0.05 (1.97)	0.05 (1.97)
Size	XS, S, M, L, XL	
Packaging	300 gloves/box 10 boxes/case	
Grip	Textured fingertips	



NITRILE SafeTouch® MicroDefense

An added level of protection

Kill up to 99.999% of tested microbes within just five minutes of contact. Safe for skin and oral use.



Item	1230	
Colour	Violet Blue ●	
Weight	3.2 g	
Thickness: mm (mil)	Finger 0.07 (2.76)	Palm 0.05 (1.97)
Size	XS, S, M, L, XL	
Packaging	SX to L: 200 gloves/box XL: 180 gloves/box 10 boxes/case	
Grip	Textured fingertips	



NITRILE SafeTouch® Advanced Form

Optimal protection and tactile sensitivity

Provides an optimal balance between protection and tactile sensitivity.

Item	1177	
Colour	Blue ●	
Weight	3.5 g	
Thickness: mm (mil)	Finger 0.11 (4.33)	Palm 0.07 (2.76)
Size	XS, S, M, L, XL	
Packaging	200 gloves/box 10 boxes/case	
Grip	Textured fingertips	



NITRILE SafeTouch® Advanced Rejuvenate

Suffer from dry hands?

Infused with lanolin and vitamin E to soothe and moisturize dry hands. Versatile and highly resistant with a comfortable latex-like feel.



Item	1176	
Colour	White ○	
Weight	4 g	
Thickness: mm (mil)	Finger 0.1 (3.94)	Palm 0.08 (3.15)
Size	XS, S, M, L, XL	
Packaging	100 gloves/box 10 boxes/case	
Grip	Fully textured	



LATEX SafeTouch® Connect

Optimal comfort and grip

Versatile, low protein content latex glove.

Item	1124	
Colour	Natural ○	
Weight	6 g	
Thickness: mm (mil)	Finger 0.08 (3.15)	Palm 0.08 (3.15)
Size	XS, S, M, L, XL	
Packaging	100 gloves/box 10 boxes/case	
Grip	Fully textured	



LATEX SafeTouch® Rejuvenate

All you love in latex and more

Infused with lanolin and vitamin E to soothe and moisturize dry hands. Soft latex with an ergonomic fit.



Item	1163	
Colour	Natural ○	
Weight	6 g	
Thickness: mm (mil)	Finger 0.12 (4.72)	Palm 0.11 (4.33)
Size	XS, S, M, L, XL	
Packaging	100 gloves/box 10 boxes/case	
Grip	Textured fingertips	



VINYL SafeTouch® EverStrong

The perfect fit for low-risk activities

Versatile PVC (vinyl) with a smooth texture for optimal tactile sensitivity.

Item	1129	
Colour	Clear ○	
Weight	5.3 g	
Thickness: mm (mil)	Finger 0.06 (2.36)	Palm 0.06 (2.36)
Size	S, M, L, XL	
Packaging	100 gloves/box 10 boxes/case	
Grip	Smooth texture	

Sterile Surgical Gloves



POLYCHLOROPRENE AssureTouch™ Guard Micro

Polychloroprene combines nitrile-like strength with latex-like comfort, without the risks of Type I allergies or latex sensitization. Ideal for procedures requiring tactile sensitivity. Micro thickness allows for double gloving.

Item	1144	
Colour	Natural White ○	
Weight	9.8 g	
Thickness: mm (mil)	Finger 0.13 (5.12)	Palm 0.11 (4.33)
Size	5.5, 6, 6.5, 7, 7.5, 8, 8.5, 9	
Packaging	1 pair/pouch 50 pairs/box 4 boxes/case	
Grip	Textured palm	



POLYCHLOROPRENE AssureTouch™ Sense

Polychloroprene combines nitrile-like strength with latex-like comfort, without the risks of Type I allergies or latex sensitization. Standard thickness makes them suitable for a wide range of surgical procedures.

Item	1145	
Colour	Natural White ○	
Weight	13 g	
Thickness: mm (mil)	Finger 0.17 (6.69)	Palm 0.15 (5.91)
Size	5.5, 6, 6.5, 7, 7.5, 8, 8.5, 9	
Packaging	1 pair/pouch 50 pairs/box 4 boxes/case	
Grip	Textured palm	



POLYCHLOROPRENE AssureTouch™ Defend

Polychloroprene combines nitrile-like strength with latex-like comfort, without the risks of Type I allergies or latex sensitization. Mint green colour aids in instant identification as non-latex and quicker detection of tears when double gloving. Standard thickness makes them suitable for a wide range of surgical procedures.

Item	1146	
Colour	Mint Green ○	
Weight	13.5 g	
Thickness: mm (mil)	Finger 0.17 (6.69)	Palm 0.15 (5.91)
Size	5.5, 6, 6.5, 7, 7.5, 8, 8.5, 9	
Packaging	1 pair/pouch 50 pairs/box 4 boxes/case	
Grip	Textured palm	



POLYISOPRENE AssureTouch™ React

Polyisoprene offers latex-like elasticity and comfort without the risks of Type I allergies or latex sensitization. Thiuram and DPG-free for a reduced risk of Type IV chemical allergies. Standard thickness makes them suitable for a wide range of surgical procedures.

Item	1147	
Colour	Natural White ○	
Weight	11.9 g	
Thickness: mm (mil)	Finger 0.18 (7.09)	Palm 0.17 (6.69)
Size	5.5, 6, 6.5, 7, 7.5, 8, 8.5, 9	
Packaging	1 pair/pouch 50 pairs/box 4 boxes/case	
Grip	Textured palm	

FIND YOUR MATCH

Which criteria are most important to you? Use the chart below to find the glove that suits you best.



Non-Sterile Medical Examination Gloves

	Item	Product Name	Weight	g	Drug/ Chemical Resistance	Speciality	Gloves/ Box	Colour
POLYCHLOROPRENE	1228	Ocean Pacific® NeoNatural	Heavy	6	Chemo tested / Chemical resistant		100	Emerald
	OP-EC	Ocean Pacific® Eclipse	Light	3			200	Blue
NITRILE	OP-LM2	Ocean Pacific® Lab Master	Light	3.2	Chemo Tested	No accelerators added	200	Indigo
	1226	Ocean Pacific® TruAloe	Light	3.5	Chemo Tested	Infused with Aloe Vera/No accelerators added	200	Aloe Green
	OP-BIO	Ocean Pacific® Elements	Light	3.5	Chemo Tested / Fentanyl Safe	Anaerobic biodegradation	200	Blue
	OP-BL2	Ocean Pacific® Blue Lite	Light	3.5	Chemo Tested		200	Blue
	OP-IN	Ocean Pacific® Intuition	Light	3.5	Chemo Tested		200	Lilac
	OP-NI	Ocean Pacific® Ninja	Heavy	5.3	Chemo Tested / Fentanyl Safe		200	Black
	9985	AssureTouch™ Basics	Light	3.5	Chemo Tested		100	Blue
	9955	AssureTouch™ Noir NEW	Light	3.5			200	Black
	1186	SafeBasics® True Fit Thin	Light	3.2	Chemo Tested		300	Pink
	1185	SafeBasics® True Fit Thin	Light	3.5	Chemo Tested / Fentanyl Safe		300	Blue
	1196	SafeTouch® Advanced Comfort	Light	2.5			300	White
	1230	SafeTouch® MicroDefense	Light	3.2	Chemo Tested	Antimicrobial	XS-L: 200 XL: 180	Violet Blue
	1177	SafeTouch® Advanced Form	Light	3.5			200	Blue
	1176	SafeTouch® Advanced Rejuvenate	Medium	4		Infused with lanolin and vitamin E	100	White
LATEX	1124	SafeTouch® Connect	Heavy	6			100	Natural
	1163	SafeTouch® Rejuvenate	Heavy	6		Infused with lanolin and vitamin E	100	Natural
VINYL	1129	SafeTouch® EverStrong	Heavy	5.3			100	Clear

Light = under 3.6 g | Medium = 3.6 g - 4.5 g | Heavy = over 4.5 g

Sterile Surgical Gloves

	Item	Product Name	Weight	g	Drug/Chemical Resistance	Speciality	Gloves/Box	Colour
POLYCHLOROPRENE	1144	AssureTouch™ Guard Micro	Light	9.8			50 pairs	○ Natural White
	1145	AssureTouch™ Sense	Medium	13			50 pairs	○ Natural White
	1146	AssureTouch™ Defend	Medium	13.5			50 pairs	● Mint Green
POLYISOPRENE	1147	AssureTouch™ React	Medium	11.9			50 pairs	○ Natural White

Light = under 11 g | Medium = 11 g - 14 g | Heavy = over 14 g



Order your
FREE samples

NOTES

This image shows a single sheet of white paper with horizontal ruling lines. The lines are evenly spaced and run across the width of the page. There are no margins, text, or other markings on the paper.



Proudly headquartered in Canada.
Globally connected.

Dental CAN

medicom.com | 1-800-361-2862

Rev.05 March 2025

©2025 AMD Medicom Inc. All rights reserved.