




























<p><u>PROG. NO. 10</u> Contour 400mg FAST</p>  <p>4000 CPM 11.5 SEC.</p>	<p><u>PROG. NO. 11</u> Contour 600mg FAST</p>  <p>4000 CPM 10.5 SEC.</p>	<p><u>PROG. NO. 12</u> Contour 800mg FAST</p>  <p>4000 CPM 9.5 SEC.</p>	<p><u>PROG. NO. 13</u> Contour 400mg REG.</p>  <p>4000 CPM 10.5 SEC.</p>
<p><u>PROG. NO. 14</u> Contour 600mg REG.</p>  <p>4000 CPM 10.0 SEC.</p>	<p><u>PROG. NO. 15</u> Contour 800mg REG.</p>  <p>4000 CPM 9.5 SEC.</p>	<p><u>PROG. NO. 16</u> Sybraloy 400mg REG.</p>  <p>4000 CPM 8.5 SEC.</p>	<p><u>PROG. NO. 17</u> Sybraloy 600mg REG.</p>  <p>4000 CPM 7.5 SEC.</p>
<p><u>PROG. NO. 18</u> Sybraloy 800mg REG.</p>  <p>4000 CPM 7.0 SEC.</p>	<p><u>PROG. NO. 19</u> Sybraloy 400mg FAST</p>  <p>4000 CPM 8.5 SEC.</p>	<p><u>PROG. NO. 20</u> Sybraloy 600mg FAST</p>  <p>4000 CPM 7.5 SEC.</p>	<p><u>PROG. NO. 21</u> Sybraloy 800mg FAST</p>  <p>4000 CPM 7.0 SEC.</p>
<p><u>PROG. NO. 22</u> Tytin FC 400mg FAST</p>  <p>4000 CPM 7.0 SEC.</p>	<p><u>PROG. NO. 23</u> Tytin FC 600mg FAST</p>  <p>4000 CPM 7.0 SEC.</p>	<p><u>PROG. NO. 24</u> Tytin FC 800mg FAST</p>  <p>4000 CPM 7.0 SEC.</p>	<p><u>PROG. NO. 25</u> Tytin FC 400mg REG.</p>  <p>4000 CPM 7.0 SEC.</p>

<p><u>PROG. NO. 26</u> Tytin FC 600mg REG.</p>  <p>4000 CPM 7.0 SEC.</p>	<p><u>PROG. NO. 27</u> Tytin FC 800mg REG.</p>  <p>4000 CPM 7.0 SEC.</p>	<p><u>PROG. NO. 28</u> Tytin 400mg SLOW</p>  <p>4000 CPM 8.5 SEC.</p>	<p><u>PROG. NO. 29</u> Tytin 600mg SLOW</p>  <p>4000 CPM 7.5 SEC.</p>
<p><u>PROG. NO. 30</u> Tytin 800mg SLOW</p>  <p>4000 CPM 7.0 SEC.</p>	<p><u>PROG. NO. 31</u> Tytin 400mg REG.</p>  <p>4000 CPM 8.5 SEC.</p>	<p><u>PROG. NO. 32</u> Tytin 600mg REG.</p>  <p>4000 CPM 7.5 SEC.</p>	<p><u>PROG. NO. 33</u> Tytin 800mg REG.</p>  <p>4000 CPM 7.0 SEC.</p>
<p><u>PROG. NO. 34</u> Contour Tablets</p>  <p>3700 CPM 10.0 SEC.</p>	<p><u>PROG. NO. 35</u> Tytin Tablets</p>  <p>3400 CPM 9.0 SEC.</p>	<p><u>PROG. NO. 36</u> Sybraloy Tablets</p>  <p>3650 CPM 10.0 SEC.</p>	<p><u>PROG. NO. 37</u> Glass Ionomer Capsules</p> <p>3900 CPM 13.0 SEC.</p>
<p><u>PROG. NO. 38</u> Glass Ionomer Capsules</p> <p>4000 CPM 10.0 SEC.</p>	<p><u>PROG. NO. 39</u> Amalgam Capsules Spherical-Particle</p> <p>3700 CPM 8.0 SEC.</p>	<p><u>PROG. NO. 40</u> Amalgam Capsules Dispersed-Particle</p> <p>3800 CPM 10.0 SEC.</p>	