


# VITABLOCS® RealLife—ingeniously three-dimensional.

Truly natural anterior esthetics with a simple mouse click



Introducing the only machinable block that mimics the tooth's natural enamel-layered-over-dentin design for anterior esthetics that look just like natural teeth. Custom positioning in the CEREC® software allows a nearly infinite choice of 3-dimensional placement for amazing individual

**VITA SHADE,  
VITA MADE. VITA®**

effects. And while the extreme esthetic possibilities may be new, VITABLOCS RealLife benefit from being composed of the same fine-structure feldspar ceramic found in highly successful VITABLOCS Mark II (proven over 25 years in millions of cases) making RealLife the perfect combination of proven longevity and new technology.



**Contact Patterson Dental at 800-873-7683**

## VITABLOCS® RealLife for CEREC®/inLab MC XL

VITABLOCS RealLife are the first CAD/CAM machinable ceramic blocks to offer 3-dimensional esthetics, utilizing the latest Sirona software solution.



### Benefits

- Fast, simple and natural chairside anterior restorations
- More accurate and detailed simulation of natural tooth structure and esthetics
- Design flexibility with movement in all three spatial directions
- Position of the restoration in the block can be changed with the software at will to reproduce the shade effect in regard to translucency, chroma, and lightness
- Eliminates the need for characterizing with stains or individualizing with layered veneering materials
- Maximizes the potential of CAD/CAM technology
- Made of VITABLOCS fine-structure feldspar ceramic, clinically proven millions of times over 25 years

### Indications

Highly esthetic full and partial anterior crowns and veneers; up to and including premolars

### Properties

Flexural strength values of approx. 150 MPa and modulus of elasticity of approx. 45 GPa

### Product type/shape

VITABLOCS feldspar ceramic in size RL-14/14 (14 x 14 x 18mm) and 6 shades: 0M1C\*, 1M2C, 2M1C,\* 2M2C, 2M3C\* and 3M2C.

\*Available January, 2011

VITA SHADE,  
VITA MADE. **VITA®**



### Special Introductory Offer

Buy 2 Packages of VITABLOCS®  
RealLife 14/14 and receive  
**FREE VITA Patient Mirror\***

To redeem offer: Purchased product must be reflected on one Patterson Dental invoice. Fax your invoice to VITA Marketing at 714-961-6209 by 12/31/10. This offer may not be combined with any other offer. Vident, a VITA company, reserves the right to change or cancel this offer at any time.



| Vident Item Number | Description                 | Shade | Quantity |
|--------------------|-----------------------------|-------|----------|
| VX55-4480          | RealLife 14/14 (14x14x18mm) | 1M2C  | 5-pack   |
| VX55-4484          | RealLife 14/14 (14x14x18mm) | 2M2C  | 5-pack   |
| VX55-4491          | RealLife 14/14 (14x14x18mm) | 3M2C  | 5-pack   |



**Contact Patterson Dental at 800-873-7683**