


## VITABLOCS® TriLuxe *forte* All-New Finer Layering



Is Esthetic CAD/CAM Dentistry Your *forte*?



**VITA®**

The all-new VITABLOCS TriLuxe *forte* provides excellent CAD/CAM esthetics. Based on the highly successful TriLuxe bloc, *forte* features even finer nuances of shade transition from the enamel to the neck. With more intense color saturation and increased fluorescence in

the cervical area, *forte* blocs provide unsurpassed natural shade effects and a chameleon blend with natural dentition. Contact your Patterson representative today and discover the high esthetics of VITABLOCS TriLuxe *forte* for yourself. [www.vident.com](http://www.vident.com)



**Contact Patterson Dental at 800-873-7683**

# VITABLOCS® TriLuxe *forte*

## Features/Benefits

- Shades in the VITA 3D-Master® system provide consistent, accurate shade matching.
- Based on VITABLOCS Mark II fine-structure feldspar ceramic, with proven clinical reliability in more than 14 million restorations.
- Abrasion behavior like natural tooth enamel, reducing wear on opposing dentition.
- Fast and simple characterization with Vita Akzent stains or shading pastes for lifelike restorations.
- Enhanced esthetics with cutback and layering using VITAVM®9 veneering porcelain.
- Superior milling properties and protection of milling instruments in the CEREC and inLab system.
- Increased color saturation and fluorescence.
- Available in large block size to accommodate large crowns.



Patterson No.	Description	Shade	Quantity
707-5781	TF-12 (10x12x15mm)	1M2C	10-pack
707-5799	TF-12 (10x12x15mm)	2M2C	10-pack
707-5807	TF-12 (10x12x15mm)	3M2C	10-pack
707-5716	TF-14 (12x14x18mm)	1M2C	5-pack
707-5724	TF-14 (12x14x18mm)	2M2C	5-pack
707-5732	TF-14 (12x14x18mm)	3M2C	5-pack
707-5633	TF-14/14 (14x14x18mm)	1M2C	5-pack
707-5641	TF-14/14 (14x14x18mm)	2M2C	5-pack
707-5658	TF-14/14 (14x14x18mm)	3M2C	5-pack



**Contact Patterson Dental at 800-873-7683**

©Vident 2007. VITABLOCS, VM and 3D-Master are registered trademarks of VITA Zahnfabrik. \*Offer valid through April 30, 2008. Fax your Patterson Dental invoice to Vident at 714-961-6209 and we will ship your Vitablocs guide directly. Allow 2-4 weeks for delivery. Not valid with any other offer.