



Richmond

CATALOG



THROUGHOUT THE HEALTH CARE WORLD, WHEN THE NAME RICHMOND IS MENTIONED, QUALITY COTTON PRODUCTS IMMEDIATELY COME TO MIND.

It is no wonder – we have been manufacturing them since 1895.

What began as cotton pellets and cotton balls expanded into braided cotton rolls, sponges, pound rolls and header rolls; today we offer a full line of premium quality, single use, disposable products and dispensers designed to make life easier for the health care professional. Our production is vertically integrated which allows us to control the manufacturing quality process from purifying the cotton through to the finished product. Innovation has long been the hallmark for the Richmond brand for over 100 years.

It has been well over a century since Richmond started manufacturing and we still hold true to our core values: quality in product, integrity in service, and loyalty to customers.

Richmond is a division of Barnhardt Manufacturing Company. Barnhardt is a diversified company consisting of Barnhardt Natural Fibers Group, Carolina Absorbent Cotton, Intrinsic, and North Carolina Foam Industries.



RichmondDentalAndMedical.net

855-841-0743



Braided Cotton Rolls



DIAMETERS:

small | medium | large

LENGTHS:

3/4" | 1-1/2" | 4" | 6"

Over the years, Richmond Dental & Medical has continually provided the highest quality Braided Cotton Roll; this superb product has become the preferred choice of dental professionals. Our attention to creating the best dental roll has resulted in a product with superior absorbency, flexibility and comfort.

- **Does not stick to mucosa.** The 100% cotton braided roll is held together with silky yarn instead of chemicals (which cause the roll to stick to mucosa), resulting in ultimate comfort for the patient.
- **Unique wicking action** provides a drier field and outstanding durability as our Braided Rolls absorb up to 30% faster.
- **Made with purified medical-grade cotton**
- **Innovative flexible design** for a secure fit into the buccal curvature and the sublingual space.
- **Infection control roll dispensers are available;** they are fully autoclavable and will keep rolls sterile until use; available in blue and white.



CODE	PRODUCT	DESCRIPTION	PACKAGING	CASE COUNT	ORIGIN
□ 200215	1-1/2" ECONOMY PAK	MED. DIA. (3/8") - NON-STERILE	2500/BAG	4 BAGS (10,000 ROLLS)	USA
□ 200204	1 1/2" JUNIOR PAK	MED. DIA. (3/8") - NON-STERILE	2000/BOX	12 BOXES	USA
□ 201201	1 1/2"	LG. DIA. (1/2") - NON-STERILE	1000/BOX	12 BOXES	USA
□ 200206	1 1/2"	MED. DIA. (3/8") - STERILE	2000/BOX	12 BOXES	USA
□ 201210	1 1/2"	SM. DIA. (5/16") - NON-STERILE	2000/BOX	12 BOXES	USA
□ 201205	6"	LG. DIA. (1/2") - NON-STERILE	100/BOX	4 BOXES	USA
□ 201208	6"	MED. DIA. (3/8") - NON-STERILE	200/BOX	4 BOXES	USA
□ 201213	6"	SM. DIA. (5/16") - NON-STERILE	200/BOX	4 BOXES	USA
□ 201226	4"	MED. DIA. (3/8") - NON-STERILE	250/BOX	4 BOXES	USA
□ 201218	3/4" PEDO ROLLS	MED. DIA. (3/8") - NON-STERILE	1400/BOX	4 BOXES	USA
□ 201219	3/4" PEDO ROLLS	SM. DIA. (5/16") - NON-STERILE	1400/BOX	4 BOXES	USA
□ 200245	VARIETY PAK	ASSORTED SIZES - NON-STERILE	1325/BOX	2 BOXES	USA
□ 200224	INFECTION CONTROL ROLL DISPENSER	BLUE - AUTOCLAVABLE	1 DISP + 200 ROLLS	6 DISPENSERS	USA
□ 200424	INFECTION CONTROL ROLL DISPENSER	WHITE - AUTOCLAVABLE	1 DISP + 200 ROLLS	6 DISPENSERS	USA

All products must be ordered in full case quantities. Please contact Richmond customer service for information about quantity discounts and product category discounts.

 **Richmond**

RichmondDentalAndMedical.net
855-841-0743

 **PURIFIED COTTON**

 **cotton**

 **Made in USA**

Manufactured where grown, in the USA for Richmond, Charlotte, NC

6/15/17

Wrapped Rolls



AVAILABLE IN:
Sterile
Non-sterile

DIAMETER:
Medium

LENGTH:
1-1/2"

Our Wrapped Rolls are a great, reliable alternative among dental rolls. Their soft, nonwoven absorbent casing contains no starch which means they will not stick to mucosa.

- 100% purified medical-grade cotton interior wrapped in nonwoven fabric
- Sealed with adhesive approved by the FDA
- Starch-free; will not stick to mucosa
- Infection control roll dispensers are available; they are fully autoclavable and will keep rolls sterile until use; available in blue and white.



CODE	PRODUCT	DESCRIPTION	PACKAGING	CASE COUNT	ORIGIN
□ 200404	1-1/2"	MED. DIA. (3/8") - NON-STERILE	2000/BOX	12 BOXES	USA
□ 200406	1-1/2"	MED. DIA. (3/8") - STERILE	2000/BOX	12 BOXES	USA
□ 200224	INFECTION CONTROL ROLL DISPENSER	BLUE - AUTOCLAVABLE	1 DISP + 200 ROLLS*	6 DISPENSERS	USA
□ 200424	INFECTION CONTROL ROLL DISPENSER	WHITE - AUTOCLAVABLE	1 DISP + 200 ROLLS*	6 DISPENSERS	USA

*Dispensers are shipped with Braided Rolls only.

All products must be ordered in full case quantities. Please contact Richmond Dental customer service for information about quantity discounts and product category discounts.

Econo Rolls



DIAMETER:
Medium

LENGTH:
1-1/2"

Satisfaction in savings. Richmond Dental & Medical's Econo Rolls are for practices whose first consideration is price. Made with 100% cotton with a methylcellulose binder, the absorbent roll is smooth and rigid.

- Economical
- 100% cotton
- Smooth and rigid
- Non-Sterile
- Durable infection control roll dispensers are available; they are fully autoclavable and will keep rolls sterile until use; available in blue and white.



CODE	PRODUCT	DESCRIPTION	PACKAGING	CASE COUNT	ORIGIN
□ 216206	1-1/2"	MED. DIA. (3/8") - NON-STERILE	2000/BOX	12 BOXES	CHINA
□ 200224	INFECTION CONTROL ROLL DISPENSER	BLUE - AUTOCLAVABLE	1 DISP + 200 ROLLS*	6 DISPENSERS	USA
□ 200424	INFECTION CONTROL ROLL DISPENSER	WHITE - AUTOCLAVABLE	1 DISP + 200 ROLLS*	6 DISPENSERS	USA

*Dispensers are shipped with Braided Rolls only.

All products must be ordered in full case quantities. Please contact Richmond Dental customer service for information about quantity discounts and product category discounts.

MultiPly®

AVAILABLE IN:

2x2 (8-Ply)

3x3 (8-Ply)

4x4 (8-Ply)

Richmond Dental & Medical 8 layer MultiPly® Rayon/Polyester Nonwoven Sponges and are perfect for oral surgery and periodontic applications. These additional layers provide twice the absorption of traditional sponges.

- Most absorbent nonwoven sponge available
- 8-plys provide thickness perfect for bite pressure
- Maintains dry surface
- Lint-free
- Autoclavable, leak proof dispensers available for 2x2 and 4x4 sizes, and come in two colors: blue and white
- MultiPly® sponges are also available in pre-sterilized SteriPocket®. These sponges are pre-packaged sterile sponges and are ideal to distribute to patients after extractions.



CODE	PRODUCT	DESCRIPTION	PACKAGING	CASE COUNT	ORIGIN
□ 300637	2"x2"	RAYON/POLY 8-PLY - NON-STERILE	100 SPONGES/SLEEVE	30 SLVS	ISRAEL
□ 300638	3"x3"	RAYON/POLY 8-PLY - NON-STERILE	100 SPONGES/SLEEVE	25 SLVS	ISRAEL
□ 300639	4"x4"	RAYON/POLY 8-PLY - NON-STERILE	100 SPONGES/SLEEVE	20 SLVS	ISRAEL
□ 300101	2"x2" SPONGE DISPENSER	BLUE - LEAKPROOF	1 DISPENSER	N/A	USA
□ 300103	2"x2" SPONGE DISPENSER	WHITE - LEAKPROOF	1 DISPENSER	N/A	USA
□ 300111	4"x4" SPONGE DISPENSER	BLUE - LEAKPROOF	1 DISPENSER	N/A	USA
□ 300113	4"x4" SPONGE DISPENSER	WHITE - LEAKPROOF	1 DISPENSER	N/A	USA
□ 303120	2"x2", STERIPOCKET	2 RAYON/POLY SPONGES/PK - STERILE	200 PKS/BX	4 BOXES	ISRAEL
□ 303122	4"x4", STERIPOCKET	2 RAYON/POLY SPONGES/PK - STERILE	100 PKS/BX	4 BOXES	ISRAEL

All products must be ordered in full case quantities. Please contact Richmond Dental customer service for information about quantity discounts and product category discounts.

100% Cotton Nonwoven Sponges

AVAILABLE IN:

2x2 (4-Ply)

3x3 (4-Ply)

4x4 (4-Ply)

The 100% Cotton Nonwoven Sponge is the perfect general use sponge combining the most desirable characteristics of sponges into one product. It features a smooth, textured and lint-free surface that provides exceptional patient comfort and maximum absorbency.

- Lint-free so no fibers are left in your patient's mouth
- All-natural and biodegradable
- Larger sizes are ideal for surface and instrument cleaning
- Autoclavable, leak proof dispensers available for 2"x2" and 4"x4" sizes, and come in two colors: blue and white
- High absorbency

CODE	PRODUCT	DESCRIPTION	PACKAGING	CASE COUNT	ORIGIN
□ 300631	2"x2"	100% COTTON - NON-STERILE	200 SPONGES/SLEEVE	25 SLVS	USA
□ 300632	3"x3"	100% COTTON - NON-STERILE	200 SPONGES/SLEEVE	20 SLVS	USA
□ 300633	4"x4"	100% COTTON - NON-STERILE	200 SPONGES/SLEEVE	10 SLVS	USA
□ 300101	2"x2" SPONGE DISPENSER	BLUE - LEAKPROOF	1 DISPENSER	N/A	USA
□ 300103	2"x2" SPONGE DISPENSER	WHITE - LEAKPROOF	1 DISPENSER	N/A	USA
□ 300111	4"x4" SPONGE DISPENSER	BLUE - LEAKPROOF	1 DISPENSER	N/A	USA
□ 300113	4"x4" SPONGE DISPENSER	WHITE - LEAKPROOF	1 DISPENSER	N/A	USA

All products must be ordered in full case quantities. Please contact Richmond Dental customer service for information about quantity discounts and product category discounts.

Rayon/Poly Nonwoven Sponges

AVAILABLE IN:

2x2 (4-Ply)

3x3 (4-Ply)

4x4 (4-Ply)

Richmond Dental & Medical's Nonwoven Rayon/Polyester Sponges are an economical and lint-free alternative. They will not leave fibers in dental work, nor will they leave lint residue on instruments when used for cleaning.

- Absorbent and strong
- Perfect for wipes
- Great for disinfecting or cleaning instruments
- Handkerchief folded
- Autoclavable, leak proof dispensers available for 2"x2" and 4"x4" sizes, and come in two colors: blue and white



CODE	PRODUCT	DESCRIPTION	PACKAGING	CASE COUNT	ORIGIN
□ 300634	2"x2"	RAYON/POLY - NON-STERILE	200 SPONGES/SLEEVE	25 SLVS	ISRAEL
□ 300635	3"x3"	RAYON/POLY - NON-STERILE	200 SPONGES/SLEEVE	20 SLVS	ISRAEL
□ 300636	4"x4"	RAYON/POLY - NON-STERILE	200 SPONGES/SLEEVE	10 SLVS	ISRAEL
□ 300101	2"x2" SPONGE DISPENSER	BLUE - LEAKPROOF	1 DISPENSER	N/A	USA
□ 300103	2"x2" SPONGE DISPENSER	WHITE - LEAKPROOF	1 DISPENSER	N/A	USA
□ 300111	4"x4" SPONGE DISPENSER	BLUE - LEAKPROOF	1 DISPENSER	N/A	USA
□ 300113	4"x4" SPONGE DISPENSER	WHITE - LEAKPROOF	1 DISPENSER	N/A	USA

All products must be ordered in full case quantities. Please contact Richmond Dental customer service for information about quantity discounts and product category discounts.

All Gauze Sponges



AVAILABLE IN:

2x2 (8-ply or 12-ply)

3x3 (12-ply)

4x4 (8-ply or 12-ply)

Proof that a product that works well can also be simple and inexpensive. All Gauze Sponges have a loose weave, giving them a textured, cheesecloth-like appearance and a mildly abrasive surface, which means they provide a good grip.

- Simple yet effective
- For use in a wide variety of applications
- Mildly abrasive texture ideal for cleaning instruments and surfaces
- Autoclavable, leak proof dispensers available for 2"x2" and 4"x4" sizes, and come in two colors: blue and white



CODE	PRODUCT	DESCRIPTION	PACKAGING	CASE COUNT	ORIGIN
□ 300420	2"x2"	ALL-GAUZE 12-PLY - NON-STERILE	200 SPONGES/SLEEVE	25 SLVS/CASE	CHINA
□ 300430	2"x2"	ALL-GAUZE 8-PLY - NON-STERILE	200 SPONGES/SLEEVE	25 SLVS/CASE	CHINA
□ 300421	3"x3"	ALL-GAUZE 12-PLY - NON-STERILE	200 SPONGES/SLEEVE	20 SLVS/CASE	CHINA
□ 300422	4"x4"	ALL-GAUZE 8-PLY - NON-STERILE	200 SPONGES/SLEEVE	20 SLVS/CASE	CHINA
□ 300423	4"x4"	ALL-GAUZE 12-PLY - NON-STERILE	200 SPONGES/SLEEVE	10 SLVS/CASE	CHINA
□ 300101	2"x2" SPONGE DISPENSER	BLUE - LEAKPROOF	1 DISPENSER	N/A	USA
□ 300103	2"x2" SPONGE DISPENSER	WHITE - LEAKPROOF	1 DISPENSER	N/A	USA
□ 300111	4"x4" SPONGE DISPENSER	BLUE - LEAKPROOF	1 DISPENSER	N/A	USA
□ 300113	4"x4" SPONGE DISPENSER	WHITE - LEAKPROOF	1 DISPENSER	N/A	USA

All products must be ordered in full case quantities. Please contact Richmond Dental customer service for information about quantity discounts and product category discounts.

Cotton-Filled Sponges



AVAILABLE IN:

2x2 (8-Ply)
(Sterile, Non-sterile)

3x3 (4-Ply)
(Non-sterile)

4x4 (4-Ply)
(Non-sterile)

This variety of sponge looks similar to the all-gauze variety with its cheesecloth-like outer layer, but contains an inner layer of 100% cotton. This provides extra absorbency, as well as additional patient comfort. The cotton-filled sponge is dentists' traditional, reliable choice.

- All gauze with an inner layer of 100% cotton
- Cotton center provides greater absorbency
- Autoclavable, leak proof dispensers available for 2x2 and 4x4 sizes, and come in two colors: blue and white
- 2"x2" size available in pre-sterilized SteriPocket®. These sponges are pre-packaged sterile sponges and are ideal to distribute to patients after extractions.



CODE	PRODUCT	DESCRIPTION	PACKAGING	CASE COUNT	ORIGIN
□ 300405	2"x2"	COTTON FILLED 8-PLY - STERILE	200 SPONGES/SLEEVE	25 SLVS	CHINA
□ 300410	2"x2"	COTTON FILLED 8-PLY - NON-STERILE	200 SPONGES/SLEEVE	25 SLVS	CHINA
□ 300311	3"x3"	COTTON FILLED 4-PLY - NON-STERILE	200 SPONGES/SLEEVE	20 SLVS	CHINA
□ 300312	4"x4"	COTTON FILLED 4-PLY - NON-STERILE	200 SPONGES/SLEEVE	10 SLVS	CHINA
□ 300101	2"x2" SPONGE DISPENSER	BLUE - LEAKPROOF	1 DISPENSER	N/A	USA
□ 300103	2"x2" SPONGE DISPENSER	WHITE - LEAKPROOF	1 DISPENSER	N/A	USA
□ 300111	4"x4" SPONGE DISPENSER	BLUE - LEAKPROOF	1 DISPENSER	N/A	USA
□ 300113	4"x4" SPONGE DISPENSER	WHITE - LEAKPROOF	1 DISPENSER	N/A	USA
□ 300408	2"x2", COTTON FILLED SPONGES	2 SPONGES/PK - STERILE	300 PK/BX	8 BOXES	CHINA

All products must be ordered in full case quantities. Please contact Richmond Dental customer service for information about quantity discounts and product category discounts.



THE DENTAL
ADVISOR
++++

SteriPocket®



AVAILABLE IN:
2x2 & 4x4 Rayon/poly
nonwoven (8-Ply)
2x2 Cotton-filled (4-Ply)



As hard as dental professionals work, any effort to make life easier and more efficient is always welcome. That's why Richmond Dental & Medical's SteriPocket® Sponges continue to be as popular as they are effective. SteriPocket® Sponges are individual, pre-sterilized packets of 2 cotton-filled or 2 rayon/poly nonwoven sponges. Easy, economical and professional, they are always ready when you need them.



- Perfect for dispensing to patients after an extraction
- Patients leave office confident that the sponges remain sterile until opened
- Saves time Assistants do not have to sterilize and package sponges for patients.
- Saves money Assistants can easily control the number of sponges distributed.

CODE	PRODUCT	DESCRIPTION	PACKAGING	CASE COUNT	ORIGIN
☐ 303120	2"x2", RAYON/POLY SPONGES	2 SPONGES/PK - STERILE	200 PKS/BX	4 BOXES	ISRAEL
☐ 303122	4"x4", RAYON/POLY SPONGES	2 SPONGES/PK - STERILE	100 PKS/BX	4 BOXES	ISRAEL
☐ 300408	2"x2", COTTON FILLED SPONGES	2 SPONGES/PK - STERILE	300 PKS/BX	8 BOXES	CHINA

All products must be ordered in full case quantities. Please contact Richmond Dental customer service for information about quantity discounts and product category discounts.

Earloop Face Mask

LEVEL I

AVAILABLE IN:

Anti-fluid
Anti-fluid, Anti-fog
Anti-fluid, Anti-fog
with shield

COLORS:

Pink
Blue
Lavender

They are as integral to a dental procedure as the most expensive equipment in the room. After all, you can't work without protection from airborne germs and fluids. Fortunately, our earloop face masks give you just such protection without sacrificing breathability. For added protection, we offer a face mask with an eye protection shield.

- 100% latex-free, including the earloops
- Face mask pleats cover nose, mouth and chin for full protection
- 3-ply spunbound fabric provides comfort, filtration and excellent breathability
- Particle filtration efficiency of $\geq 95\%$ at 1.0 micron
- Soft aluminum nosepiece maintains shape when formed to face and is fully enclosed for comfort

CODE	PRODUCT	DESCRIPTION	PACKAGING	CASE COUNT	ORIGIN
□ 400501	ANTI-FLUID	PINK	50 MASKS/BOX	6 BOXES	USA
□ 400502	ANTI-FLUID	BLUE	50 MASKS/BOX	6 BOXES	USA
□ 400509	ANTI-FLUID	LAVENDER	50 MASKS/BOX	6 BOXES	USA
□ 400503	ANTI-FLUID, ANTI-FOG	PINK	50 MASKS/BOX	6 BOXES	USA
□ 400504	ANTI-FLUID, ANTI-FOG	BLUE	50 MASKS/BOX	6 BOXES	USA
□ 400510	ANTI-FLUID, ANTI-FOG	LAVENDER	50 MASKS/BOX	6 BOXES	USA
□ 400505	ANTI-FLUID, ANTI-FOG W/SHIELD	PINK	25 MASKS/BOX	4 BOXES	USA
□ 400506	ANTI-FLUID, ANTI-FOG W/SHIELD	BLUE	25 MASKS/BOX	4 BOXES	USA
□ 400511	ANTI-FLUID, ANTI-FOG W/SHIELD	LAVENDER	25 MASKS/BOX	4 BOXES	USA

All products must be ordered in full case quantities. Please contact Richmond Dental customer service for information about quantity discounts and product category discounts.

Reflective Shields Plus®



AVAILABLE IN:

- Large
- Small
- Dual Pack

Absorbs 2x's faster than the leading competitors & absorbs up to 30% more!



The Reflective Shields Plus® are exceptionally popular with both dentists and patients. And it's no wonder: they are uniquely crafted to provide absorption and superior illumination during common procedures. The cheek retractor specially designed heart shape and distinctive reflective film are both innovative and practical.



- Does not swell in mouth and compromise the work area
- Greater absorbency providing a dry field within seconds
- Added comfort; Special heart shape easily accommodates buccal curvatures and keeps shield in check with the parotid gland
- Unique reflective film allows illumination of the oral cavity
- Durable nonwoven absorbent material for moisture retention
- Ideal for various procedures such as sealants, composite fillings, restorations and crown preps

CODE	PRODUCT	DESCRIPTION	PACKAGING	CASE COUNT	ORIGIN
600700	REFLECTIVE SHIELDS PLUS®	LARGE REFLECTIVE MOISTURE ABATEMENT PADS	50/BOX	20 BOXES	USA
600710	REFLECTIVE SHIELDS PLUS®	SMALL REFLECTIVE MOISTURE ABATEMENT PADS	50/BOX	20 BOXES	USA
600720	REFLECTIVE SHIELDS PLUS®	DUAL PACK	25-LG, 25-SM/BOX	20 BOXES	USA

All products must be ordered in full case quantities. Please contact Richmond Dental customer service for information about quantity discounts and product category discounts.



Cotton Pellets

SIZES:

#0 - #5

REFILLS:

#1 - #4

They are small, but that doesn't keep them from doing their job well. Richmond Dental & Medical's 100% Cotton Pellets are designed to provide outstanding quality and absorbency in a variety of uses.

- 100% cotton - ideal for applying medicaments
- Highly absorbent
- Chrome tilt-top dispensers available for sizes #1 through #4 they hold one size each
- Plastic temple dispensers for use with sizes #1 through #4 they hold four sizes at a time
- EZ UP® Pellet Dispenser for sizes #1 through #5



CODE	PRODUCT	DESCRIPTION	PACKAGING	CASE COUNT	ORIGIN
□ 100105	SIZE 0	11/32" DIAMETER	400 PELLETS/BOX	6 BOXES	USA
□ 100106	SIZE 1	9/32" DIAMETER	2000 PELLETS/BOX	6 BOXES	USA
□ 100107	SIZE 2	7/32" DIAMETER	2550 PELLETS/BOX	6 BOXES	USA
□ 100108	SIZE 3	5/32" DIAMETER	2500 PELLETS/BOX	12 BOXES	USA
□ 100109	SIZE 4	1/8" DIAMETER	3000 PELLETS/BOX	12 BOXES	USA
□ 100110	SIZE 5	3/32" DIAMETER	1300 PELLETS/BOX	12 BOXES	USA
□ 100900	ASSORTED REFILLS	2 - #1'S, 4 - #2'S, 4 - #3'S, 2 - #4'S		1 BOX	USA
□ 100901	#1 REFILLS	9/32" DIAMETER	12 REFILLS/BOX	1 BOX	USA
□ 100902	#2 REFILLS	7/32" DIAMETER	12 REFILLS/BOX	1 BOX	USA
□ 100903	#3 REFILLS	5/32" DIAMETER	12 REFILLS/BOX	1 BOX	USA
□ 100904	#4 REFILLS	1/8" DIAMETER	12 REFILLS/BOX	1 BOX	USA
□ 600091	TILTTOP #1 PELLETS	PELLET DISPENSER - CHROME	1 DISP. W/PELLETS	6 DISPENSERS	SWITZ.
□ 600092	TILTTOP #2 PELLETS	PELLET DISPENSER - CHROME	1 DISP. W/PELLETS	6 DISPENSERS	SWITZ.
□ 600093	TILTTOP #3 PELLETS	PELLET DISPENSER - CHROME	1 DISP. W/PELLETS	6 DISPENSERS	SWITZ.
□ 600094	TILTTOP #4 PELLETS	PELLET DISPENSER - CHROME	1 DISP. W/PELLETS	6 DISPENSERS	SWITZ.
□ 600601	TEMPLE DISP.	PELLET DISPENSER - PLASTIC	1 DISP. + 10 REFILLS	1 DISPENSER	USA
□ 600095	PELLET DISPENSER	1 DISPENSER, 1 LARGE GRID, 1 SMALL GRID	1	1	SWITZ

All products must be ordered in full case quantities. Please contact Richmond Dental customer service for information about quantity discounts and product category discounts.

Richmond

RichmondDentalAndMedical.net
855-841-0743



Manufactured where grown, in the USA for Richmond, Charlotte, NC

6/15/17

Cotton Tip Applicators



LENGTHS:
3"
6"

Cotton tipped applicators are useful when applying medicaments or in general wound care. Three and six inch length sizes are available.



- 100% purified medical grade cotton tip
- Secure wooden stem
- 100 per package; 10 packages per case

CODE	PRODUCT	DESCRIPTION	PACKAGING	CASE COUNT	ORIGIN
400401	COTTON-TIPPED APPLICATORS	6" 100 PER PAK NON-STERILE	10 PACKS/BOX	10 BOXES	USA
400402	COTTON-TIPPED APPLICATORS	3" 100 PER PAK NON-STERILE	10 PACKS/BOX	10 BOXES	USA

All products must be ordered in full case quantities. Please contact Richmond Dental customer service for information about quantity discounts and product category discounts.

Medical Balls



AVAILABLE IN:

Large (Non-sterile)

Medium (Non-sterile)

Even products that seem as simple as this can make a difference to your patients. Their versatility and ease of application will make a difference to your practice.

- 100% cotton, rayon and polyester
- Ideal for general use
- High absorbency
- Soft and compact

CODE	PRODUCT	DESCRIPTION	CASE COUNT	ORIGIN
□ 187124P	COTTON BALLS	MEDIUM, NON-STERILE	100 BALLS/BAG, 36 BAGS/CASE	USA
□ 187123P	COTTON BALLS	MEDIUM, NON-STERILE	300 BALLS/BAG, 36 BAGS/CASE	USA
□ 187106P	COTTON BALLS	MEDIUM, NON-STERILE	2,000 BALLS/BAG, 2 BAGS/CASE	USA
□ 187100P	COTTON BALLS	LARGE, NON-STERILE	1,000 BALLS/BAG, 2 BAGS/CASE	USA
□ 187404P	RAYON BALLS	MEDIUM, NON-STERILE	2,000 BALLS/BAG, 2 BAGS/CASE	USA
□ 187407P	RAYON BALLS	LARGE, NON-STERILE	1,000 BALLS/BAG, 2 BAGS/CASE	USA
□ 187301P	POLY BALLS	LARGE, NON-STERILE	1,000 BALLS/BAG, 2 BAGS/CASE	USA
□ 187300P	POLY BALLS	MEDIUM, NON-STERILE	2,000 BALLS/BAG, 4 BAGS/CASE	USA

All products must be ordered in full case quantities. Please contact Richmond Dental customer service for information about quantity discounts and product category discounts.

Cotton Pound Roll



AVAILABLE IN:
Sterile
Non-sterile

WIDTH:
12 inches

WEIGHT:
1 pound

The cotton pound roll by definition is a 12" wide cotton roll weighing 16 ounces. The support and comfort will be felt by your patient.

- 100% purified medical-grade cotton
- Multi-purpose use
- High absorbency
- Perfect for cushioning or dressing

CODE	PRODUCT	DESCRIPTION	PACKAGING	CASE COUNT	CASE DIM (INCH)	CASE WEIGHT	ORIGIN
1098401	COTTON POUND ROLL	12" WIDTH, NON-STERILE	1 ROLL	25/CASE	22L X 22W X 13L	28 LBS	USA
1098403	COTTON POUND ROLL	12" WIDTH, STERILE	1 ROLL	12/CASE	14L X 18W X 14H	16 LBS	USA

All products must be ordered in full case quantities. Please contact Richmond Dental customer service for information about quantity discounts and product category discounts.

Infection Control Roll Dispenser



COLORS:
blue
white

Richmond Dental & Medical's durable IC roll dispensers are designed for infection control, dispensing one sterile cotton roll at a time. Our standard size works with any 1.5" medium cotton roll.

- Sterilize in autoclave or chemclave
- Keeps rolls sterile until use
- Dispenser comes packaged with 200 braided cotton rolls
- Available in blue and white



CODE	PRODUCT	DESCRIPTION	PACKAGING	CASE COUNT	ORIGIN
□ 200224	INFECTION CONTROL ROLL DISPENSER	BLUE - AUTOCLAVABLE	1 DISP + 200 ROLLS	6 DISPENSERS	USA
□ 200424	INFECTION CONTROL ROLL DISPENSER	WHITE - AUTOCLAVABLE	1 DISP + 200 ROLLS	6 DISPENSERS	USA

All products must be ordered in full case quantities. Please contact Richmond customer service for information about quantity discounts and product category discounts.

Solution Trays

PERFECT FOR DISPENSING THE FOLLOWING MEDICAMENTS:

Chlorhexidine
Alcohol
Formocresol*
Liquid Etch
Tubulicid*
Hemodent
Caries Indicator
Primer and Bond*
Irrigating Solution

The Richmond Dental & Medical's Solution Tray is an elegant chair-side tray, complete with dropper bottles and dishes, designed to dispense a variety of solutions accurately, conveniently and economically, while reducing the risk of cross contamination. The Richmond Dental & Medical Solution Tray is a "must have" for every dental practice.

- The anti-reflux pump offers drop-by-drop control of the flow, allowing for economical dispensing of solutions
- Color-coded dropper bottles and dishes allow immediate, accurate identification
- Variety of solutions on one tray saves time and is convenient
- Each tray arrives with two Dappen dishes of each color
- Available with a choice of #1, #2, #3 or #4 Tiltop dispenser, specify on order
- *Opaque glass bottles for light-sensitive solutions also are available upon request
- Autoclave sterilization, except for the pumping device, which should be sterilized, with Ethylene Oxide or with a 3% Sodium Hypochlorite solution.

CODE	PRODUCT	DESCRIPTION	PACKAGING	CASE COUNT	ORIGIN
□ 440720	SOLUTION TRAY W/TILTOP	3 BOTTLE TRAY W/CHOICE OF TILTOP	EACH	1 BOX	SWITZERLAND
□ 440730	SOLUTION TRAY W/4 BOTTLES	4 BOTTLE TRAY	EACH	1 BOX	SWITZERLAND
□ 440731	RED DAPPEN DISH	2 DISHES PER BOX	EACH	1 BOX	SWITZERLAND
□ 440732	YELLOW DAPPEN DISH	2 DISHES PER BOX	EACH	1 BOX	SWITZERLAND
□ 440733	BLUE DAPPEN DISH	2 DISHES PER BOX	EACH	1 BOX	SWITZERLAND
□ 440734	GREEN DAPPEN DISH	2 DISHES PER BOX	EACH	1 BOX	SWITZERLAND
□ 440735	OPAQUE BOTTLE*	1 BOTTLE PER BOX	EACH	1 BOX	SWITZERLAND
□ 440736	BOTTLE/RED PUMP	1 BOTTLE PER BOX	EACH	1 BOX	SWITZERLAND
□ 440737	BOTTLE/YELLOW PUMP	1 BOTTLE PER BOX	EACH	1 BOX	SWITZERLAND
□ 440738	BOTTLE/BLUE PUMP	1 BOTTLE PER BOX	EACH	1 BOX	SWITZERLAND
□ 440739	BOTTLE/GREEN PUMP	1 BOTTLE PER BOX	EACH	1 BOX	SWITZERLAND
□ 440740	PIPETTES	3 COLOR CODED PIPETTES		1 BOX	SWITZERLAND

All products must be ordered in full case quantities. Please contact Richmond Dental customer service for information about quantity discounts and product category discounts.



RichmondDentalAndMedical.net

855-841-0743

FoamEZ® Mirror Wipes



AVAILABLE IN:
200 count bags

“The FoamEZ® Mirror Wipe makes patient care easier from beginning to end. In addition to being the easiest wipe to remove fog and debris, it retracts tissue gently and reduces positional changes in the mouth.” – WILLIAM L. BARHAM, DDS



- FoamEZ® Mirror Wipe removes fog and debris easily and quickly with one wipe against the mirror
- Eliminates the need to remove mirror from the mouth or turn mirror against (into) the cheek
- Disposable FoamEZ® Mirror Wipe slides easily onto the air/water syringe
- Application of defogger on the FoamEZ® Mirror Wipe prevents fogging during procedures

CODE	PRODUCT	DESCRIPTION	PACKAGING	CASE COUNT	ORIGIN
410100	FOAMEZ® MIRROR WIPE	DISPOSABLE FOAM MIRROR WIPE	200/BAG	10 BAGS	USA

All products must be ordered in full case quantities. Please contact Richmond Dental customer service for information about quantity discounts and product category discounts.

FoamEZ® Splash Hood Liner



AVAILABLE IN:
2 count package



Designed by a dentist, this very practical product from Richmond Dental and Medical can save a practice untold time and money. The FoamEZ® Splash Hood Liner is simple to use. Saving one prosthesis from breakage makes it a very worth while investment.



- FoamEZ® Splash Hood Liner provides protection against chips or cracks in the event the prosthesis releases from the hand during polishing
- FoamEZ® Splash Hood Liner helps keep the Splash Hood clean by covering the inner surface and catching flying debris
- Removal of the FoamEZ® Splash Hood Liner makes for a fast and easy clean up

CODE	PRODUCT	DESCRIPTION	PACKAGING	CASE COUNT	ORIGIN
410000	SPLASH HOOD LINER	FOAM LINER FOR SPLASH HOOD	2 LINERS/BAG	10 BAGS	USA

All products must be ordered in full case quantities. Please contact Richmond Dental customer service for information about quantity discounts and product category discounts.

NEW!

20

**Path-O-Guard®
Earloop
Face Mask
LEVEL III**

Raising Protection to a New Level
A LEVEL III EARLOOP FACE MASK YOU CAN TRUST



AVAILABLE IN:

**Anti-fluid
Anti-fluid, Anti-fog**

COLORS:

**Pink
Blue
Lavender**



- Level III
- Exceeds current ASTM-2100 standard for high barrier surgical face masks
- 3-ply, comfortable and non-irritating
- Latex free
- Bacterial filtration efficiency of $\geq 98\%$ at 1.0 micron
- Available in Anti-Fluid and Anti-Fluid/Anti-Fog
- Available in three colors: blue, pink and lavender



"There was no fogging and they felt extremely soft on the inside."

"I especially liked the anti-fog feature and the confidence of protection with a Level III Mask."

CODE	PRODUCT	DESCRIPTION	PACKAGING	CASE COUNT	ORIGIN
□ 400581	ANTI-FLUID	PINK	50 MASKS/BOX	6 BOXES	USA
□ 400582	ANTI-FLUID	BLUE	50 MASKS/BOX	6 BOXES	USA
□ 400589	ANTI-FLUID	LAVENDER	50 MASKS/BOX	6 BOXES	USA
□ 400583	ANTI-FLUID, ANTI-FOG	PINK	50 MASKS/BOX	6 BOXES	USA
□ 400584	ANTI-FLUID, ANTI-FOG	BLUE	50 MASKS/BOX	6 BOXES	USA
□ 400590	ANTI-FLUID, ANTI-FOG	LAVENDER	50 MASKS/BOX	6 BOXES	USA

All products must be ordered in full case quantities. Please contact Richmond Dental customer service for information about quantity discounts and product category discounts.

NEW!



21

EZ UP[®]
Pellet
Dispensers

Safe and Simple Dispenser for Cotton Pellets

- Innovative design
- Gravity-based expulsion system
- High chemical resistant materials
- Simple to use and easy to disinfect
- Dispenses pellet sizes #1 through #5
- Interchangeable size grid

CODE	PRODUCT	DESCRIPTION	PACKAGING	CASE COUNT	ORIGIN
□ 600095	PELLET DISPENSER	1 DISPENSER, 1 LARGE GRID, 1 SMALL GRID	1	1	SWITZ

All products must be ordered in full case quantities. Please contact Richmond Dental customer service for information about quantity discounts and product category discounts.



www.RichmondDentalAndMedical.net | 855-841-0743

