

# AloePRO™ Synthetic



## Soothe Your Hands and Ease Your Mind.

DASH AloePRO Synthetic exam gloves are your latex-free solution to better feeling gloves.

- Active Aloe™ inner surface soothes dry irritated hands by improving your skin's ability to hydrate itself
- Made from innovative butadiene technology with stress relaxation properties that are superior to nitrile
- Priced lower than most nitrile gloves
- Far superior fit than vinyl gloves
- Six sizes for a precise fit

Available in X-Small, Small, Small-Medium, Medium, Large, X-Large

**DASH**<sup>®</sup>  
**MEDICAL GLOVES**

Where Safety Meets Comfort™



Phone: 1-800-523-2055 Fax: 1-800-523-7795 [www.dashmedical.com](http://www.dashmedical.com)

©2008 DASH Medical Gloves Inc. DASH is a registered trademark of DASH Medical Gloves Inc., registered in the US P&T Office

# AloePRO™ Synthetic Glove Features

**Material:** butadiene  
**Type:** non-sterile  
**Color:** green  
**Glove Interior:** powder free  
**Glove Exterior:** textured surface  
**Shape:** ambidextrous  
**Cuff:** beaded  
**Usage:** single use only  
**Quantities:** 100 gloves per dispenser box, by weight

**Protein Content:** none  
**Quality Certification:** ISO 9001, FDA QSR  
**Tensile Strength**  
   **Before Aging (MPa):** 14  
   **After Aging (MPa):** 14  
**Elasticity/Elongation**  
   **Before Aging:** 500%  
   **After Aging:** 400%

Glove Size	Palm Width (mm)	Length (mm)	Thickness at Palm	Thickness at Finger
X-Small	74 (±4)	250 (±10)	0.10 mm (±0.015)/3.9 mil	0.125 mm (±0.02)/4.9 mil
Small	81 (±4)	250 (±10)	0.10 mm (±0.015)/3.9 mil	0.125 mm (±0.02)/4.9 mil
Small-Medium	89 (±4)	250 (±10)	0.10 mm (±0.015)/3.9 mil	0.125 mm (±0.02)/4.9 mil
Medium	97 (±4)	250 (±10)	0.10 mm (±0.015)/3.9 mil	0.125 mm (±0.02)/4.9 mil
Large	107 (±4)	250 (±10)	0.10 mm (±0.015)/3.9 mil	0.125 mm (±0.02)/4.9 mil
X-Large	115 (±4)	250 (±10)	0.10 mm (±0.015)/3.9 mil	0.125 mm (±0.02)/4.9 mil