


# AloePRO™

## Experience the Healing Power of Aloe.

Now you can experience the healing powers of aloe vera first-hand with DASH AloePRO exam gloves. Coated on the inside with soothing Active Aloe™, AloePRO is designed to help moisturize and soften your hands while you work. And while your hands are being moisturized on the inside, the outside of AloePRO gloves are fully textured to give you the firm grip you need.

Available in X-Small, Small, Small-Medium, Medium, Large, X-Large



CAUTION: This product contains natural rubber latex which may cause allergic reactions.



**DASH**<sup>®</sup>  
MEDICAL GLOVES

Where Safety Meets Comfort™

Phone: 1-800-523-2055 Fax: 1-800-523-7795 [www.dashmedical.com](http://www.dashmedical.com)

©2008 DASH Medical Gloves Inc. DASH is a registered trademark of DASH Medical Gloves Inc., registered in the US P&TM Office

# AloePRO™ Glove Features

**Material:** natural rubber latex  
**Type:** non-sterile  
**Color:** natural  
**Glove Interior:** powder free, aloe vera coated  
**Glove Exterior:** textured surface  
**Shape:** ambidextrous  
**Cuff:** beaded  
**Usage:** single use only  
**Quantities:** 100 gloves per dispenser box, by weight

**Protein Content:** less than 50 µg/g per ASTM D 5712  
**Quality Certification:** ISO 9002, FDA QSR  
**Tensile Strength**  
**Before Aging (MPa):** 18 min  
**After Aging (MPa):** 14 min  
**Elasticity/Elongation**  
**Before Aging:** 650%  
**After Aging:** 500%

Glove Size	Palm Width (mm)	Length (mm)	Thickness at Palm	Thickness at Finger
X-Small	75 (±3)	250 (±10)	0.14 mm (±0.03)/5.5 mil	0.155 mm (±0.03)/6.1 mil
Small	85 (±3)	250 (±10)	0.14 mm (±0.03)/5.5 mil	0.155 mm (±0.03)/6.1 mil
Small-Medium	89 (±3)	250 (±10)	0.14 mm (±0.03)/5.5 mil	0.155 mm (±0.03)/6.1 mil
Medium	96 (±3)	250 (±10)	0.14 mm (±0.03)/5.5 mil	0.155 mm (±0.03)/6.1 mil
Large	107 (±3)	250 (±10)	0.14 mm (±0.03)/5.5 mil	0.155 mm (±0.03)/6.1 mil
X-Large	117 (±3)	250 (±10)	0.14 mm (±0.03)/5.5 mil	0.155 mm (±0.03)/6.1 mil