


Luminance

Aurelia® Luminance™ Fully Textured
Powder Free Latex Examination Gloves



**Maximum Grip
Optimum Comfort**

- Fully textured for maximum grip
- Enhanced polymer coating technology for optimum comfort & increased gripping action
- One of a kind, 100% visually inspected for highest quality assurance



AURELIA
Sign of Comfort™

www.aureliagloves.com | 1-877-AURELIA (1-877-287-3542)

Luminance

Product Description

Aurelia® Luminance™ Powder Free Latex Examination gloves are made from 100% natural Malaysian rubber, ambidextrous, creamy white color, micro-textured and non-sterile. Aurelia gloves conform and are tested to the highest international standards. For more information visit www.aureliagloves.com



Aurelia® Luminance™ Specifications

TENSILE TEST RESULTS (unaged)

	Tensile Strength (MPa)	Elongation at break (%)
Extra Small	30.0	970
Small	31.1	970
Medium	31.0	970
Large	31.2	960
Extra Large	31.3	970
ASTM 3578:05	18 MPa minimum	650% minimum

AGED TENSILE TEST RESULTS (7 Days at 70° C)

	Tensile Strength (MPa)	Elongation at break (%)
Extra Small	28.0	900
Small	28.0	920
Medium	27.5	910
Large	27.3	890
Extra Large	27.9	900
ASTM 3578:05	14 MPa minimum	500% minimum

DIMENSIONS

	Length (inches/mm)	Thickness (Single-wall)	
		Palm (mils/mm)	Finger (mils/mm)
Extra Small Aurelia Luminance ASTM 3578:05	9.4 / 240 8.7 / 220 minimum	5.0 / 0.12 3.2 / 0.08 minimum	6.0 / 0.15 3.2 / 0.08 minimum
Small Aurelia Luminance ASTM 3578:05	9.4 / 240 8.7 / 220 minimum	5.0 / 0.12 3.2 / 0.08 minimum	6.0 / 0.15 3.2 / 0.08 minimum
Medium Aurelia Luminance ASTM 3578:05	9.4 / 240 9.1 / 230 minimum	5.0 / 0.12 3.2 / 0.08 minimum	6.0 / 0.15 3.2 / 0.08 minimum
Large Aurelia Luminance ASTM 3578:05	9.4 / 240 9.1 / 230 minimum	5.0 / 0.12 3.2 / 0.08 minimum	6.0 / 0.15 3.2 / 0.08 minimum
Extra Large Aurelia Luminance ASTM 3578:05	9.4 / 240 9.1 / 230 minimum	5.0 / 0.12 3.2 / 0.08 minimum	6.0 / 0.15 3.2 / 0.08 minimum

ITEM NUMBER

Extra Small	Small	Medium	Large	Extra Large
29885	29886	29887	29888	29889

Ordering Information

To order from your preferred distributor, specify size and item number. Packing: 100 pieces per dispenser box, 10 dispenser boxes per shipping carton.



Authorized Aurelia Distributor:

CAUTION STATEMENT: This product contains natural rubber latex which may cause allergic reactions. Safe use of this glove by or on latex sensitized individuals has not been established.

AURELIA Div. Supermax Inc.

2225 White Oak Circle, Aurora, IL 60502 | Toll Free: 1-877-AURELIA | Tel: 1-630-898-8886 | Fax: 1-630-898-8855
Web: www.aureliagloves.com | Email: info@aueliagloves.com