

MicroPro™

Get a Grip!

DASH uses high quality, low modulus latex that provides a barrier which is so form fitting, you hardly notice you're wearing gloves at all. Unlike some gloves with only a partially textured surface, MicroPro are fully textured on both the fingers and palm giving you the confidence of a firm grip in wet or dry conditions.

Available in X-Small, Small, Medium, Large, X-Large



BOXES ARE
COLOR-CODED
BY SIZE

DASH[®]
MEDICAL GLOVES

Where Safety Meets Comfort™

CAUTION: This product contains natural rubber latex which may cause allergic reactions.



Phone: 1-800-523-2055 Fax: 1-800-523-7795 www.dashmedical.com

©2008 DASH Medical Gloves Inc. DASH is a registered trademark of DASH Medical Gloves Inc., registered in the US P&TM Office

MicroPro™ Glove Features

Material: natural rubber latex
Type: non-sterile
Color: natural
Glove Interior: powder free
Glove Exterior: textured surface
Shape: ambidextrous
Cuff: beaded
Usage: single use only
Quantities: 100 gloves per dispenser box, by weight

Protein Content: less than 50 µg/g per ASTM D 5712
Quality Certification: ISO 9002, FDA QSR
Tensile Strength
 Before Aging (MPa): 18 min
 After Aging (MPa): 16 min
Elasticity/Elongation
 Before Aging: 650%
 After Aging: 500%

Glove Size	Palm Width (mm)	Length (mm)	Thickness at Palm	Thickness at Finger
X-Small	76 (±5)	240 (±5)	0.14 mm (±0.01)/5.5 mil	0.16 mm (±0.01)/6.3 mil
Small	84 (±5)	240 (±5)	0.14 mm (±0.01)/5.5 mil	0.16 mm (±0.01)/6.3 mil
Medium	94 (±5)	240 (±5)	0.14 mm (±0.01)/5.5 mil	0.16 mm (±0.01)/6.3 mil
Large	105 (±5)	240 (±5)	0.14 mm (±0.01)/5.5 mil	0.16 mm (±0.01)/6.3 mil
X-Large	110 (±5)	240 (±5)	0.14 mm (±0.01)/5.5 mil	0.16 mm (±0.01)/6.3 mil