


When was the last time you saw a void on a patient's x-ray?



Not with these! Bottom's up.


AccuDose.® The disposable tubes and plugs that eliminate repetitive pad-to-mouth motion by putting the entire mixture at your fingertips.

With AccuDose, you can quickly and more cost-effectively dispense any type of composite or cement while reducing voids and porosities in any restoration. And because AccuDose tubes and plugs are disposable, you completely eliminate cleanup and cross contamination. No matter what application or viscosity material you need to dispense there is an AccuDose tube and plug to get the job done right.

- Save time by having the entire mixture of material right at your fingertips and eliminate the constant pad-to-mouth motion
- Reduce voids and porosities in any restoration
- Available in a variety of styles for specific applications
- Tubes and plugs are easy to load
- Disposable to eliminate cleanup and cross-contamination
- PhotoBloc™ orange tint, blocks light and extends the life of light-activated materials
- Can be used with all syringes currently on the market
- Trusted and clinically proven for over 35 years



Ordering is easy. Call 800-235-5862 or contact your preferred dealer.
Visit centrixdental.com for more information.



CENTRIX
AccuDose®

Disposable Tubes & Plugs



centrix®

Making Dentistry Easier.™



C-R[®] Syringes, tubes and plugs were invented by Centrix over 35 years ago. Centrix continues to be the world leader in providing delivery systems for all dental materials.

ORDERING INFORMATION:



Low Viscosity

An all purpose tip for medium viscosity, flowable materials such as composites, core materials, glass ionomers and cements. The tip can be clipped to a larger opening if needed.

REF 210021

AccuDose Low Viscosity, 100

REF 210023

AccuDose Low Viscosity, 500



High Viscosity

A tip for high viscosity, heavy materials such as thick composites and core buildup materials. Large opening allows placement of even the heaviest bodied materials into small areas.

REF 280021

AccuDose High Viscosity, 100

REF 280023

AccuDose High Viscosity, 500



19ga & 20ga NeedleTubes[™]

A tip with a 19ga or 20ga bendable blunt needle, for the most precise placement of materials such as cements, base materials, core buildup materials and any other gel type materials. Especially useful for injecting materials into postholes and deep proximal boxes.

REF 290031

AccuDose NeedleTubes, 100
20ga (Photobloc Orange)

REF 290033

AccuDose NeedleTubes, 500
20ga (Photobloc Orange)

REF 290032

AccuDose NeedleTubes, 100
19ga (Melon Green)

REF 290035

AccuDose NeedleTubes, 500
19ga (Melon Green)




RibbonTubes

A tip with a crimped, flat end for dispensing a thin and wide ribbon of material. Especially useful for loading luting cements or composites when placing veneers.

REF 290041

AccuDose RibbonTubes, 100



Placing materials with C-R Tubes & Plugs creating void- and porosity-free restorations



Placing material traditionally with composite placement instruments, spatulas or spiral paste fillers

Ordering is easy.

Call 1-800-235-5862
or contact your preferred dealer.

centrixdental.com



centrix[®]

Making Dentistry Easier.SM



Our Guarantee — If you are not 100% satisfied, simply return the product for a complete refund. It's that easy.

This product is covered by one or more of the following U.S. and foreign patents: 4963093; 5052927; 5122057; 5336088; P19000310-1; 42 00 044.0; 3209773; and 9103646-7.

© 2006 Centrix, Inc. All rights reserved. AccuDose is a registered trademark and NeedleTube is a trademark of Centrix, Inc.