

# Filtek™ Supreme Ultra

## Flowable Restorative



So many uses,  
so simple to use.

# Versatile

## Flowables provide a fast, easy solution for a wide range of cases.

In fact, they have even been called “dental duct tape.” But ours is different. More than a quick fix, you need to be confident your restorative will perform and handle beautifully using your technique. That’s why dentists are choosing Filtek™ Supreme Ultra Flowable Restorative.

Its enhanced properties make it **ideal for multiple indications**—adding versatility to its list of impressive benefits.



### Why dentists use a flowable

#### ... as a liner base

- Excellent handling
- Cavity adaptation
- Stress relief

#### ... as a filling material

- Fast and easy
- Little or no polish needed
- Smooth transition from material to tooth

## Retains its polish and wears better than other leading flowables.\*

The same **true nanotechnology** used by Filtek™ Supreme Ultra Universal Restorative gives our flowable advantages such as:

- Excellent esthetics
- Unsurpassed polish retention\*
- Low wear

### When you need a fast, easy solution, Filtek Supreme Ultra flowable restorative is ideal for a variety of indications.

- Class III and V restorations
- Restoration of minimally invasive cavity preparations (including small, non stress-bearing occlusal restorations)
- Base/liner under direct restorations
- Repair of small defects in esthetic indirect restorations
- Pit and fissure sealant
- Undercut blockout
- Repair of resin and acrylic temporaries

#### As a Class V Restoration



#### As a Minimally Invasive Restoration



#### As a liner under a direct restoration

(Filtek™ Supreme Ultra Universal Restorative is the composite shown.)



Photos courtesy of Dr. Gunnar Reich.

\*Versus competitive materials tested

# Flowable

## Flows when you want it to ... stays where you need it.

How a flowable handles greatly affects the efficiency and outcome of restorative work. It's so important that 3M ESPE purposefully adds a unique handling modifier to its formula to assure it delivers **the flow-on-demand behavior that dentists prefer.**

- Under pressure of extrusion or being manipulated with instruments, the material becomes more liquid ... and quickly flows where it's needed.
- Once placed and the flow ceases, it returns to its original viscosity and holds its shape—there's no slumping or running.

### Exceptional handling confirmed by dentists.

In a field evaluation, **60% of dentists** rated the **flow** of Filtek Supreme Ultra flowable restorative **better than Tetric EvoFlow.\*\***


Dentists also rated Filtek Supreme Ultra flowable restorative **better than Tetric EvoFlow** in these key areas:

- Cavity adaptation
- Holding its shape
- Lack of stickiness

## Filtek Supreme Ultra Flowable Restorative ...

## so many uses, it's virtually indispensable.

We've reduced our shrinkage  
by nearly **20%!**



Bonded disc method

Source: 3M ESPE internal data



## Ordering Information

### Filtek™ Supreme Ultra Flowable Restorative

Refill	Item #	Item #
SHADE	SYRINGE	<b>New</b> CAPSULE
A1	6032A1	6033A1
A2	6032A2	6033A2
A3	6032A3	6033A3
A3.5	6032A3.5	6033A3.5
A4	6032A4	6033A4
B1	6032B1	6033B1
B2	6032B2	6033B2
C2	6032C2	6033C2
D2	6032D2	6033D2
White	6032W	6033W
Extra White	6032XW	6033XW
OA3	6032OA3	6033OA3



Filtek™ Supreme Ultra Flowable Restorative is now available in capsules and offers the same flow-on-demand handling you've enjoyed from the syringes.

**SYRINGE SHADE REFILLS** include 2 – 2g syringes; 20 dispensing tips; Instructions for Use.

**CAPSULE SHADE REFILLS** include 20 – 0.2g capsules; Instructions for Use.

Customer Care Center: 1-800-634-2249 [www.EveryoneHasAShade.com](http://www.EveryoneHasAShade.com)

**3M ESPE**

#### Dental Products

3M Center  
Building 275-2SE-03  
St. Paul, MN 55144-1000  
U.S.A.  
1-800-634-2249

#### 3M Canada

Post Office Box 5757  
London, Ontario N6A 4T1  
Canada  
1-888-363-3685

Please recycle. Printed in U.S.A.  
© 3M 2010. All rights reserved.  
70-2013-0413-9

3M, ESPE and Filtek are trademarks of 3M or 3M ESPE AG. Used under license in Canada. EsthetX, Flow-It, Revolution 2, Tetric EvoFlow, TPH3, Venus Flow and Vertice Flow are not trademarks of 3M or 3M ESPE AG.