

hi-tray



A RIGID DISPOSITION



Clinical
Impression Trays

Zhermack[®]
CLINICAL



hi-tray

The solution for every impression.

The *complete range* of *impression trays* for fixed dentures (*metal* and *plastic*) and for removable dentures (*non-deformable plastic*).



hi-tray metal



hi-tray light clear



hi-tray light plastic



hi-tray light edentulous



hi-tray



RIGID

› The rigidity *prevents distortion* during the setting of the impression.

STERILIZABLE

› *Resistant* to the *high temperatures* in an autoclave.

RETENTIVE

› *High degree of anchorage* of the impression material guaranteed by the conformation of the impression tray.

VERSATILE

› *Complete range* of impression trays in *stainless steel* and *plastic* for all clinical situations.

hi-tray metal

Reliable and durable.

Impression tray for standard applications
In *special BSA stainless steel* with 16.70% chrome.



HI-TRAY METAL

FEATURES

- Possible to cold *sterilise* and in autoclave
- Available in *two versions*:
solid
perforated

ADVANTAGES

- *High resistance to deformation*
- *Greater precision* of the impression
- *“Rim lock” design* for *perfect mechanical retention* of the material to the impression tray
- *Optimum degree of finishing*
- *Long duration*

hi-tray light

clear and plastic

Disposable and versatile.

HI-TRAY LIGHT CLEAR, the impression tray in *clear plastic* for standard applications.

HI-TRAY LIGHT PLASTIC, the impression tray in *grey rigid plastic* for standard applications.

FEATURES

- Made in *transparent* plastic, ideal for chromatic impression materials (*HI-TRAY LIGHT CLEAR*), or *grey* plastic (*HI-TRAY LIGHT PLASTIC*)
- Available in the following *versions*:
small (upper and lower)
medium (upper and lower)
large (upper and lower)
partial (upper, right and left)



HI-TRAY LIGHT CLEAR



HI-TRAY LIGHT PLASTIC

ADVANTAGES

- *Rigid and non-deformable*
- *Disposable*
- *Can be sterilised in autoclave*
- *Adaptable*
- *Perforated design* for *improved mechanical retention*
- *Reinforced grip* with ribbing for *improving the distribution of the pressure* exerted

hi-tray light edentulous

Anatomic and resistant.

Anatomic impression trays in **green rigid plastic** for edentulous arches.



HI-TRAY LIGHT EDENTULOUS

FEATURES

- Made in **green** plastic
- Available in the following **versions**:
 - small** (upper and lower)
 - medium** (upper and lower)
 - large** (upper and lower)

ADVANTAGES

- **Rigid** and **non-deformable**
- **Disposable**
- **Adaptable**
- **Perforated design** for **improved mechanical retention**
- **Reinforced grip** with ribbing for **improving the distribution of the pressure** exerted



hi-tray
Aimed rigidity.



HI-TRAY LIGHT PLASTIC 12 PZ.

		SIZE			
		Small	Medium	Large	Partial
TYPE	Upper	D5GSUP	D5GMUP	D5GLUP	D5GPUP
	Lower	D5GSLO	D5GMLO	D5GLLO	-
	Left	-	-	-	D5GPL
	Right	-	-	-	D5GPR

HI-TRAY LIGHT CLEAR 12 PZ.

		SIZE			
		Small	Medium	Large	Partial
TYPE	Upper	D5CSUP	D5CMUP	D5CLUP	D5CPUP
	Lower	D5CSLO	D5CMLO	D5CLLO	-
	Left	-	-	-	D5CPL
	Right	-	-	-	D5CPR

HI-TRAY METAL

D5K10S	1 kit of 10 trays in upper and lower sizes SOLID
D5K10P	1 kit of 10 trays in upper and lower sizes. PERFORATED
D5K08S	1 kit of 8 trays in upper and lower sizes. SOLID
D5K08P	1 kit of 8 trays in upper and lower sizes. PERFORATED
D5K06S	1 kit of 6 trays in upper and lower sizes. SOLID
D5K06P	1 kit of 6 trays in upper and lower sizes. PERFORATED

HI-TRAY LIGHT EDENTULOUS 12 PZ.

		SIZE		
		Small	Medium	Large
TYPE	Upper	D5ESUP	D5EMUP	D5ELUP
	Lower	D5ESLO	D5EMLO	D5ELLO

HI-TRAY METAL SOLID

		SIZE				
		1	2	3	4	5
TYPE	Superiore	D51SS	D52SS	D53SS	D54SS	D55SS
	Inferiore	D51IS	D52IS	D53IS	D54IS	D55IS

HI-TRAY METAL PERFORATED

		SIZE				
		1	2	3	4	5
TYPE	Superiore	D51SP	D52SP	D53SP	D54SP	D55SP
	Inferiore	D51IP	D52IP	D53IP	D54IP	D55IP

BRANCHES

GERMANY
Zhermack GmbH
Deutschland
Öhlmühle 10 - D-49448 Marl
Tel. +49 - 0 54 43 / 20 33 - 0
Fax +49 - 0 54 43 / 20 33 - 11
info@zhermack.de
www.zhermack.com

POLAND
Zhermapol Sp. z o.o.
ul. Augustówka 14
02 - 981 Warszawa - Polska
Tel. +48 - 22 858 73 11
+48 - 22 858 73 41
Fax +48 - 22 642 07 14
biuro@zhermapol.pl
www.zhermack.pl

USA
Zhermack Inc.
P.O. Box 4195, River Edge
New Jersey 07661-4195
Tel. +1 (877) 819-6206
+1 (732) 389 8540
Fax +1 (732) 389 8543
info@zhermackusa.com
www.zhermackusa.com

COUNTRY CONTACTS

UK AND IRELAND
+44 - (0) 7870 690811
uk@zhermack.com

SPAIN
Tel. 900 99 39 52
info.spain@zhermack.com

FRANCE
Tel. 0800 - 915083
info.france@zhermack.com

RUSSIA
Tel. +7 916 930 1191
info.russia@zhermack.com

MEXICO AND CENTRAL AMERICA
Tel. +52 - 1 - 55 - 91980497
info.mexico@zhermack.com

SOUTH AMERICA
Tel. +57 - 310 - 2601777
info.colombia@zhermack.com

Zhermack SpA

Via Bovazecchino, 100 - 45021 Badia Polesine (RO) ITALY
Tel. +39 - 0425 597611 Fax +39 - 0425 53596 info@zhermack.com www.zhermack.com