

## Gloves – Patterson Dental - Master # 090-2601

For use in the sterilization area, and for operatory disinfection. Puncture and abrasion resistant nitrile compound. Flocked-lined for easier hand insertion and removal. Mold embosses for improved tactile feel and handling of instruments. Heat sterilizable in Chemiclave or steam heat to 270°F up to three times. Color is green. Available in Small, Medium and Large. Three pairs of gloves per package.



Product Numbers	Size	Description
090-2619	Small	Patterson Glove Utility 3/pair
090-2627	Medium	Patterson Glove Utility 3/pair
090-2635	Large	Patterson Glove Utility 3/pair