

Patterson Private Label HP Gun Dispensers

1:1	084-9190
4:1	084-9125
10:1	084-9117

Confirmation Suitability of Moist Heat Sterilization at 273°F (134C°), 3 minutes.

The Patterson HP Gun Dispenser successfully passed an internal test study program simulation of long-term usage in combination with multiple sterilization cycles. The mechanical properties and the functionality of the products are in compliance with the product specification after 1000 sterilization cycles. However occasional slight discoloration and surface effects can be noticed with no impact on functionality.

Confirmation Proof against Standard Disinfectant.

The Patterson HP Gun Dispenser successfully passed an internal test study program simulation of long-term usage in combination with multiple disinfection cycles. The mechanical properties and the functionality of the products are still in compliance with the product specification after 1000 disinfection cycles with standard disinfectants as follows:

Banicide®

Cidex® Plus 28 Day Solution

However occasional slight discoloration and insignificant mechanical strength reduction can be noticed with no impact on functionality.