



Self Seal/Sac pour l'autoclave
Autoclave Bag/Autoscellant

Reorder Number
No de produit
083-8805

Contents/Contenu

- Syringe/Seringue
- Forceps/Davier
- Scalars/Détartreurs
- Elevator/Élévateur
- Rongeur/Pince-gouge
- Scissors/Ciseau
- Handpiece/Pièce à main
- Probes/Sondes
- Mirrors/Miroirs
- College Pliers/Précelles
- Sponges/Compressees
- Gauze Pad/Tampons de gaze

When chemical vapor or steam sterilized, indicator will change color/L'indicateur changera de couleur après la stérilisation à la vapeur chimique ou à l'autoclave

CAUTION: This product contains natural rubber latex which may cause allergic reactions.
ATTENTION: Ce produit contient du latex naturel pouvant causer des réactions allergiques.

Misc/Divers _____

Date/Date _____


Initials/Initiales _____

Distributed by: Distribué par: _____

PATTERSON DENTAL
Corporate Office
1031 Mendota Heights Road
Saint Paul, MN 55120

PATTERSON DENTAL/DENTAIRE
Siège social
Montréal Québec
H3M 3E6

FOLD HERE/PLIER ICI



INDICATOR WILL CHANGE COLOR FROM LIGHT BLUE TO "DIRTY" GRAY

Autoclave Procedure: saturated steam at 250F, 15PSI, 20 minutes. Color change evident when exposed to these conditions.
----- Allow bags to **DRY COMPLETELY PRIOR TO REMOVAL** from sterilizer. -----
Avoid over-handling/mishandling of bags. Be careful when placing heavy or sharp instruments in the bags. Rough handling and misuse can puncture, cut or tear the bags, breaching the sterile field.

REV 9/09

PRESS/PRESSER
FIRMLY/FERMEMENT