

NEW!

DENTSPLY
CAULK



Easily identify setting stages with Jeltrate® Chroma!

Formulated to change colors for each preparation phase



0-30 sec.
Fuschia during mixing

30-45 sec.
Lavendar during tray loading

45 sec-1 min
Light blue during setting. Wait 1 minute for mouth removal.

- 5-day dimensional stability
- 2-minute Fast set
- Easy mixing
- Non-slumping
- Consistent mix provides smooth models

Jeltrate®
Chroma
Dustless Color-Changing
Alginate Impression Material

Ordering Information


605710 **CASE OF 9 CANISTERS**
Fast Set Intro Kit
Reusable Canister contains:
1 Pouch Jeltrate® Chroma Material (454g),
1 Accessory Pack

605700 **CASE OF 12 POUCHES**
Fast Set (454g each)



Jeltrate® Jeltrate® DUSTLESS Jeltrate® Plus DUSTLESS

The high algin content of **Jeltrate® impression materials** provides for a quality impression without the excessive flow you will find in alginate materials with lower algin contents. By eliminating excessive flow, Jeltrate® impression materials allow your patients a level of comfort that cannot be found in other alginate impression materials.



Jeltrate®

other alginate impression material

HIGH ALGIN CONTENT – Why is it important?

- Consistent mix.....Very smooth models
- No remakesSaves time
- High tear strengthEasily removed from undercuts without tearing
- Non-slump formulaPatient comfort

Ordering Information

Jeltrate® Alginate Impression Material

- CASE OF 9 CANISTERS**
 608503 Regular Set – Beige (454g ea)
 608522 Fast Set – Pink (454g ea)
- GROSS PACKAGE**
 608005 12 Twelve-Pouch Packages (21g ea)
 1 Water Measure

NOTE: Pouch Packages not available in Fast Set

- BULK PACKAGE**
 608058 Regular Set (22lbs / 10 kilos)
 608247 Fast Set (22lbs / 10 kilos)

Jeltrate® Dustless Dustless Alginate Impression Material

- CASE OF 8 POUCHES**
 1 Accessory Pack
 1 Empty Canister
- 605100 Regular Set (454g ea)
 605200 Fast Set (454g ea)
- 605220 **CASE OF 20 POUCHES**
 Fast Set (454g ea)

Jeltrate® Plus Antimicrobial Dustless Alginate Impression Material

- CASE OF 9 CANISTERS**
 605502 Regular Set – Beige (454g ea)
 605602 Fast Set – Pink (454g ea)
- GROSS PACKAGE**
 605505 12 Twelve-Pouch Packages (25g ea)
 1 Water Measure

NOTE: Pouch Packages not available in Fast Set

- BULK PACKAGE**
 605503 Regular Set (22lbs / 10 kilos)
 605603 Fast Set (22lbs / 10 kilos)



DENTSPLY
CAULK

38 West Clarke Avenue / Milford, DE 19963
 800.532.2855 / www.dentsply.com / www.caulk.com

All data on file. ©2012 DENTSPLY International. All rights reserved.
 Form #505900 (4/16/12)

