

Patterson	American Ortho	EOP	Sullivan	Schein/Masel	3/M Unitek	Rocky Mtn.	Ortho Org.	Ormco	TP Ortho	Dentaurum
1	1.5	31.5	31.5	31.5	31.5	9	30.5	7	29.5	5
2	2	32	32	32	32	10	31	8	30	6
3	2.5	32.5	32.5	32.5	32.5	11	31.5	9	30.5	7
4	3	33	33	33	33	12	32	10	31	8
5	3.5	33.5	33.5	33.5	33.5	13	32.5	11	31.5	9
6	4	34	34	34	34	14	33	12	32	10
7	4.5	34.5	34.5	34.5	34.5	15	33.5	13	32.5	11
8	5	35	35	35	35	16	34	14	33	12
9	5.5	35.5	35.5	35.5	35.5	17	34.5	15	33.5	13
10	6	36	36	36	36	18	35	16	34	14
11	6.5	36.5	36.5	36.5	36.5	19	35.5	17	34.5	15
12	7	37	37	37	37	20	36	18	35	16
13	7.5	37.5	37.5	37.5	37.5	21	36.5	19	35.5	17
14	8	38	38	38	38	22	37	20	36	18
15	8.5	38.5	38.5	38.5	38.5	23	37.5	21	36.5	19
16	9	39	39	39	39	24	38	22	37	20
17	9.5	39.5	39.5	39.5	39.5	25	38.5	23	37.5	21
18	10	40	40	40	40	26	39	24	38	22
19	10.5	40.5	40.5	40.5	40.5	27	39.5	24	38.5	23
20	11	41	41	41	41	28	40	26	39	24
21	11.5	41.5	41.5	41.5	41.5	29	40.5	27	39.5	25
22	12	42	42	42	42	30	41	28	40	26
23	12.5	42.5				31	41.5	29	40.5	27
24	13	43				32	42	30	41	28

Patterson Brand bands and the corresponding competitor bands are a close approximation. Due to manufacturing tolerances concerning band contour, stretching techniques and material, sizes differ slightly.